Axcel


AXL 5001, AXL 5501, AXL 5002

Éclairage Interventionnel
Surgical Light System
OP-Beleuchtung

Notice technique **FR**

Technical manual **EN**

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QUALITY COMPLIANCE

Certification of MAQUET SAS quality system

LNE/G-MED certifies that the quality system created by MAQUET SAS for the design, manufacturing, marketing, installation and customer servicing of its surgical lights meets the requirements of the following international standards:

- ISO 9001:2000
- ISO 13485:2004

CE Marking/intended use:

Compliance with the requirements of Directive 93/42/EEC dated of June 14th 1993, relating to medical devices has been assessed in accordance with Annex VII of this Directive. This AXCEL™ Surgical lights range is a class I device in accordance with Annex IX of Directive 93/42/EEC.

EC Marking affixed in 2007

CE Marking products recognition:

AXCEL surgical light in single (ceiling mount AXL 5001 or wall mount AXL 5002) or double version (ceiling mount AXL 5501) for minor surgery, emergency, recovery, speciality, examination...

NOTE: with reference to installation rules complying with IEC 60601-2-41 standard, in case of mains failure, AXCEL configurations can not be used in rooms where it is necessary to have emergency 24 V Direct current.

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 To avoid any risk of an electric shock, class I devices must be connected to a power supply system which is earthed.

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 Take care when handling the circuit boards: these boards are supplied in an electrostatic envelope and must be handled with great care. Certain components on the boards may be destroyed by static electricity (ESD, electrostatic discharge). An anti-static wrist strap (OPM 060) must be worn to eliminate this risk.

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CONTENTS

QUALITY REQUIREMENTS	2
SYMBOLS USED	4
1 - AXCEL DESCRIPTION OF PRODUCT RANGE AND CHARACTERISTICS	5
1-1 Description of product range	5
1-2 Technical characteristics	7
1-3 Technical data sheet.....	9
2 - MAINTENANCE	10
2-1 Maintenance	10
2-2 The 4 maintenance keys	10
2-3 Cleaning/Disinfection/ Sterilization	11
2-4 Adjustments	15
2-5 Maintenance kits	16
2-6 Spare parts	16
2-7 Maintenance tools.....	17
2-8 Maintenance procedures.....	17
3 - PARTS LISTS AND DRAWINGS	21
TROUBLESHOOTING.....	55

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**2-1 MAINTENANCE
MAINTENANCE AND INSPECTION**

We have taken particular care with the design and manufacture of the equipment that you have just purchased. In order to keep its initial qualities and operating reliability, regular maintenance and inspections must be performed every year.

These must be performed in compliance with the schedule given:

- *during the guarantee period by an MAQUET SAS technician, or of an approved distributor,*
- *outside the guarantee period by an MAQUET SAS technician, or of an approved distributor or of the hospital's technical department.*

• *Maquet's warranty is only valid and the safety and integrity of the product are only guaranteed if:*

- *All inspection, maintenance and repair operations are performed by Maquet engineers or trained and authorised technical support technicians.*
- *Only original accessories, consumables and spare parts are used.*

MAQUET SAS will not be held responsible in the even of operations which do not comply with these instructions.

2-2 FOUR MAINTENANCE KEYS

IF AN ANOMALY OR DAMAGE IS OBSERVED, REPLACE THE FAULTY PART.

Refer to the corresponding procedures for the detailed assembly or adjustment operations.

A. TIGHTENING

- *Main arm swivel covers*
- *Power supply cover*
- *Cupola cover*
- *Cupola handles*
- *Sterilizable handle*
- *Lamp socket contacts.*

B. ADJUSTMENT

- *Voltage $23,5 V_{EFF} (V_{RMS}^{AC+DC})$*
- *Cupola brake*
- *Balancing the spring arm*

**C. CHECKS
SAFETY POINTS**

- *Tightening of mounting screws on suspension tube and assembly of swivel joints.*
- *S suspensions tube screw are well tightened and stop screw is well located inside the lower tube hole.*
- *Attachment of spring arms (retaining rings).*
- *Cupola stop segment (disassembly and lubrication).*
- *Attachment of covers and handles.*
- ***Replace the spring arm every 6 years.***

OTHER CHECKS

- *Nominal lighting level: 100 000 lux*
- *Measure the earth bonding on all visible conductive metal parts with the earth bonding tester: $R \leq 0,1 \text{ Ohm}$.*
- *Cleanliness of optical core*
- *Condition of lamp socket and lamp*
- *Verticality of suspension tube*
- *Engagement of sterilizable handle.*

ASPECT

- *Cleaning should be performed using water and somap.*
- *Do not use alcohol to clean the surface.*

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2.3 - CLEANING AND MAINTENANCE

CLEANING/DISINFECTION/ STERILIZATION

The user must contact the sanitary specialists at the medical center. The recommended products and procedures must be observed.

In the event of doubt concerning the compatibility of active agents to be used, contact the local MAQUET SAS customer service.

2.3.1 Cleaning and disinfection of the surgical lighting system:

Before cleaning, ensure the power is turned off and the light is cooled down.

General instructions concerning cleaning, disinfection and safety:

- Remove the sterilisable handles.
- Wipe the system with a cloth moistened with a surface cleaner. Follow the manufacturer's dilution and temperature recommendations.
- Rinse the unit with a damp cloth and wipe.
- Wipe evenly with a cloth moistened with disinfectant. Follow the manufacturer's recommendations.
- Rinse with a cloth moistened with water to remove any residue (particularly when cleaners containing aldehydes, quaternary ammonium or surfactants are used).
- Wipe off with a dry cloth.
- Make sure that all liquid cleaning products used have been thoroughly wiped off.

a) Examples of recommended products

- **GETINGE USA Products:** TEC QUAT 256.
- **ANIOS Products:** HEXANIOS G + Rat 0,5% (AMONIUM IV, POLYHEXANIDE) ; ANIOSYME P.L.A (Quartärammonium, Enzyme); SALVANIOS pH10 (Quartärammonium, Guanidinium); ANIOS DDSH (Quartärammonium, Guanidinium).
- **Schülke & Mayr Products:** ANTIFECT PLUS.

b) Examples of prohibited products

Solutions containing glutaraldehyde, phenol, iodine, bleach, alcohol or chloride ions must not be used.

Do not use disinfection by fumigation methods.

2.3.2 Sterilizable handles: cleaning and sterilization

2.3.2.1 Preparation for Cleaning

Use soft cloth to wipe away gross soils from light handle surface, immediately after its use.

Place handles in containment device for transporting and keep soils moist to ease further cleaning.

Avoiding situations that soil the inside surfaces of handles. This can help ease the required cleaning efforts.

Do not scratch the lens (video handles only).

2.3.2.2 Cleaning

Submerge handles while soaking and washing with a detergent¹ solution.

Let handles soak for 15 minutes and then manually wash soils off, using soft cloth and soft brushes.

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2.3.2.3 Inspect Cleanliness of Surfaces

Inspect handle periodically to assure all soils are removed from outer and also inner surfaces. If soils are seen on internal parts, then repeat cleaning or use an ultrasonic cleaning process, as necessary.

Rinsing:

Thoroughly rinse handle using a brisk stream of clean tap water to remove cleaning solutions. Inspect handle for material properties (no cracks, chips, or crazing), cleanliness and assure it functions properly.

2.3.2.4 - Drying

Wipe dry with freshly cleaned soft cloth. Prior to reuse, the light handles are required to be steam sterilized.

2.3.2.5 - Disinfection

Handles could be disinfected with a machine (Getinge washers for example) and rinsed with a temperature maxi. of 93°C.

Stage	Exp. Temp.	Exp. Time
Prewash	18 - 35° C	60 sec
Wash	46 - 50° C	300 sec
Neutralization	41 - 43°C	30 sec
Wash 2	24 - 28°C	30 sec
Rinsing/Disinfection	92 - 93°C	600 sec
Drying		20 min

2.3.2.6 - Sterilization

Cleaned light handles are to be steam sterilized by using the following instructions and sterilization cycle:

Country	Sterilization cycle	Exp. Temp [°C]	Exp. Time [min]	Dry [min]
USA & Canada	Prevacuum ²	132 - 135	10	16
France	ATNC (Prion) (Prevacuum)	134	18	
Autres pays	Prevacuum	Follow local / national regulations requirements		

- Re-inspect each handle to assure it is clean and free of defects before proceeding.
- Wrap handle with sterilization wrapper material (double wrapper or equivalent). Alternatively, handles may be placed in paper/plastic sterilization pouches³ to ease identification after reprocessing and prior to reuse.
- Place on perforated tray with handle opening facing downwards⁴.
- Package with biological and/or chemical indicators for monitoring the sterilization process, per local regulations and facility policies.
- Run sterilization cycle, abiding to sterilizer manufacturer's instructions for load configurations.

To ensure correct sterilization, do not allow any soiling to penetrate inside the handle. With the sterilization parameters indicated, the sterilizable handles are only guaranteed for no more than 50 utilisation cycles. Disposal of the sterilizable handles is subject to the standard rules governing disposal of hazardous products in a hospital environment.

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¹ A non-enzyme-based detergent is recommended (e.g., Getinge USA NeutraWash). Many enzyme based detergents can adversely affect the plastic material used in this handle and should not be used for extended periods of soaking and will need to be thoroughly rinsed off, if used.
² This light handle is a porous instrument.
³ Possible suppliers of sterilization pouches:
Medical Action Industries
SBW Medical
Baxter International
⁴ For air removal and drying purposes

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2.4 - ADJUSTMENTS (Fig. 24 and 25)

*The only adjustments to perform are to be achieved on the spring arm (B).
No adjustment is necessary on main arm (A) nor on lighthouse (C).*

Adjustment of balance weight:

For the adjustment of the balance weight remove screw (1) at the locking sleeve (2) and move the sleeve back.

Spring arm too strong:

Move arm until adjustment nut (3) becomes visible in slot (4).

Turn adjustment nut with the rod (5) counter clock-wise, until spring arm is balanced and the weight is carried properly in every position.

Spring arm too weak:

Move arm until adjustment nut (3) becomes visible in slot (4).

Turn adjustment nut with the rod (5) clock-wise, until spring arm is balanced and the weight is carried properly in every position.

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2.8 MAINTENANCE PROCEDURE

2.8.1 Cleaning the AXCEL cupola (drawing 1)

Before performing any maintenance operation, wait for the cupola to cool down and switch off the light at the main switch.

Clean the cupola without separating the cupola from the light.

Tools required

Flat tip screwdriver, 10mm torque wrench, 9mm socket wrench.

- *Fasten a load to the spring arm-fork swivel joint.*
- *Remove the lamp socket.*
- *Remove the 3 #10 torx screws and the female plug at the optical core.*
- *Remove the cupola seal by hand*
- *Remove the three mounting screws and their spacers at the sterile handle sleeve and remove the underside.*
- *Clean the mirrors using a soft rag and plunge the optical cord in very hot water.*
- *Dry the optical core using a hair dryer or a blower.*
- *Reassemble the unit by performing the disassembly procedure in reverse order*
- *Do not use tools to reinstall the cupola seal.*

2.8.2 Lubrication of cupola-hoop swivel joint (drawings 1 and 7)

Tools required

20 torx wrench, flat tip screwdriver

- *Fasten a load to the spring arm - fork swivel joint.*
- *Remove the #20 torx screw from the swivel joint cover.*
- *Remove the stop segment and check its condition (replace it if necessary).*
- *Turnover the cupola, remove the brake and check its condition (replace it if necessary).*
- *Separate the cupola from the fork.*
- *Clean and lubricate the shaft.*

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- *Apply new grease and reassemble the cupola by performing the disassembly procedure in reverse order.*

2.8.3 Replacement of cupola handle (drawings 1 and 2)

Tools required

Flat tip screwdriver, #20 torx wrench

- *Fasten a load to the spring arm - fork swivel joint.*
- *Remove the cupola seal by hand.*
- *Remove the three mounting screws and their spacers at the sterile handle sleeve and remove the underside.*
- *Remove the 2 torx screws securing the peripheral handle*
- *Replace the handle and reassemble by performing the disassembly procedure in reverse order.*

2.8.4 Replacement of shell-hoop seal (drawing 2)

Tools required

Flat tip screwdriver, #10 torx wrench, #20 torx wrench

- *Remove the optical core (see # 3.6.1).*
- *Separate the cupola from the fork (see # 3.6.2) and place it on a work stand.*
- *Remove the 6 torx #10 screws securing the cupola wing.*
- *Remove the wing from the shell taking care not to misadjust the mirrors*
- *Replace the shell-hoop seal*
- *Reassemble by performing the disassembly procedure in reverse order.*

2.8.5 Replacement of Lexan endpiece (drawing 7)

Tools required

None

- *Carefully remove the Lexan endpiece by pulling it by hand without rotating it.*
- *Slowly uncover the fork and remove the Lexan connector.*
- *Disengage the Lexan plate from the endpiece and reassemble it by performing the disassembly procedure in reverse order.*

2.8.6 Replacing the dimmer (drawing 7)

Tools required

2 mm BTR wrench

- *Remove the Lexan endpiece (see # 3.6.4).*
- *Remove the headless BTR screw.*
- *Disconnect the two connectors.*
- *Remove the dimmer unit.*
- *Reassemble observing the cable positions.*

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Flacher Schraubendreher, #10-Torx-Schlüssel, #20-Torx-Schlüssel

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Tools required

None

- Remove the Lexan endpiece (see # 3.6.4).
- Replace the Lexan plate and reassembly by performing the disassembly procedure in reverse order.

2.8.8 Lubrication of fork-spring arm swivel joint (drawings 7 and 10)

Tools required

Flat tip screwdriver

- Remove the segment cover screw on the spring arm.
- Turn the cover by 180°.
- Remove the stop segment and check its condition (replace it if necessary).
- Remove the fork and cupola assembly with one hand while holding the spring arm with the other.
- Clean and degrease the shaft.
- Apply new grease and reassemble by performing disassembly procedure in reverse order.

2.8.9 Lubrication of main arm - spring arm swivel joint (drawings 8, 9, 10)

Tools required

Flat tip screwdriver, 2mm BTR wrench, retaining ring pliers

- Remove the fork-cupola assembly (see #3.6.7).
- Remove the BTR screw of the main arm and cover.
- Remove the cover by disengaging it vertically so that no stress is applied to the rotating contact.
- Remove the retaining ring, check its condition and replace it if necessary.
- Remove the spring arm.
- Clean, degrease and apply new grease.
- Reassemble by performing disassembly procedure in reverse order.

2.8.10 Replacement of transformer unit (drawings 8 and 9)

Tools required

2.5mm BTR wrench

- Remove the BTR screw securing the access cover and open it.
- Remove the BTR screw securing the transformer bracket.
- Disconnect the network connector located behind the tube.
- Remove the transformer unit by holding it by the bracket and tilting it through the top.
- Disconnect the connectors, secondary side, then disconnect the ground wires and remove the transformer.
- Replace the transformer and reassemble by performing the disassembly procedure in reverse order.

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2.8.11 Replacement of fuse or adjustment of voltage selector (drawings 12 and 13)

Tools required

2.5mm BTR wrench, flat tip screwdriver

- Remove the BTR screw securing the access cover and open it.
- Remove the BTR screw securing the transformer bracket.
- Slightly tilt the power supply by holding it by the bracket.
- Replace the fuse or select the network voltage.
- Reassemble by performing the disassembly procedure in reverse order.

2.8.12 Replacement of rotary casing (drawings 8 and 9)

Tools required

2 and 2.5 mm BTR wrench, retaining ring pliers, small flat tip screwdriver

- Remove the transformer unit (see # 3.6.9).
- Remove the cupolas (see #3.6.7).
- Remove the BTR screw securing the central shaft cover and remove the cover (ditto #3.6.8).
- Remove the retaining ring and its washer.
- Remove the lower arm.
- Remove the retaining ring and its washer from the upper arm.
- Lower the assembly formed by the upper arm and rotary casing along the shaft. Mark the position of the stop.
- Disengage the grommet and remove the arm with the stop and aluminum washer.
- Mark the position of the lugs in their plastic housing.
- Remove the lugs from the housing using a small flat tip screwdriver.
- Remove the two wire guides and the wire grommet.
- Remove the rotary casing.
- Replace the rotary casing.
- Install the grommet and the casing, then install the cable with a loop in the casing.
- Position the aluminum washer and the stop on the shaft.
- Finish the assembly by performing the disassembly procedure in reverse order.

Following any maintenance operation, cleaning of assembly is mandatory. Use a soft degreasing product.

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BTR-Schlüssel 2,5 mm, flacher Schraubendreher

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Caution: Disassembly of certain sub-assemblies can affect operation and safety, for example:

- when servicing the electrical power supply,
- when servicing the suspension arms and balancing system,
- when servicing the optical system of the copulas which are equipped with filters intended to eliminate radiations not visible by the patient. The unit should not be operated, under any circumstances, without these filters.

For this type of service, contact the MAQUET SAS-approved product support service.

TROUBLESHOOTING

ANOMALY	PROBABLE CAUSE	CORRECTIVE ACTION
❑ Lighthouse does not come on.	1 - Bulb failure or missing bulb.	Switch off power supply and replace bulb(s).
	2 - Other cause.	Call MAQUET SAS Service engineering department.
	3 - Mains power failure	Check if another piece of equipment is operating on the same electrical power system.
❑ All lighthouses not lit	- Each lighthouse has its own individual control.	Check that all the power supplies are properly switched on.
❑ Low bulb service life	- Improper bulb or overvoltage condition.	1) Ensure you are only using MAQUET SAS bulbs specified. Do not use bulbs used with surgical lights of previous generation (ANGENIEUX AX4-AX14) which are not compatible for power and size. 2) Verify voltage at bulb terminals. $23,5 V_{EFF} \left(V_{AC+DC}^{RMS} \right)$
❑ Sterilizable handle does not engage correctly on its support.	1 - Average service life of sterilizable handle: 50 sterilization cycles at 134° ± 4°C for 20 minutes.	Check operation of handle lock mechanism (click) and overall handle.
	2 - Maximum service life exceeded / Handle is deformed.	Replace handle
❑ Lighthouse drift	- Suspension tube not levelled. - Unstable ceiling structure	- Check level and ceiling structure. - Call MAQUET SAS Service engineering department.
❑ Lighthouse too loose or too difficult to handle.	- Brake misadjustment. - Lack of lubrication.	- Proceed with Balance Adjustment (See Chap. 3.3). - Call MAQUET SAS Service engineering department.

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SERVICE PROTOCOL	
For use during Installation / Servicing	
<h1 style="margin: 0;">Axcel</h1>	<input type="checkbox"/> Contract Order n°..... <input type="checkbox"/> Chargeable <input type="checkbox"/> Other:.....
<input type="checkbox"/> Installation Delivery date: Installation date: Commissioning date:	
<input type="checkbox"/> Servicing Preventive / Curative Valid until:	
Hospital:	
Building:	Floor level:
Service:	Operating room:
Product(s) identification	
Cupola n°1	Ref. number: S/N:
Cupola n°2	Ref. number: S/N:
Cupola n°3	Ref. number: S/N:
Overall assessment	
<input type="checkbox"/> Configuration fully operational	
<input type="checkbox"/> Configuration shall not be used until all deficiencies are corrected	
<input type="checkbox"/> Free of direct risk but deficiencies detected. May be corrected on short term.	
<input type="checkbox"/> Configuration no longer safe. Taking out of service is recommended.	
Remarks	
.....	
Signature	
Name:	Name:
Date:	Function:
Installer/Technician authorized by MAQUET	Customer

Measuring equipment used

IEC 62353 compatible control units

Type: **Voltmeter**

S/N:

Calibration valid until:

Type: **Luxmeter**

S/N:

Calibration valid until:

Type: **Megohmmeter**

S/N:

Calibration valid until:

Type:

S/N:

Calibration valid until:

Electrical safety test

The protective earth resistance is to be measured between point **0** and the point **3**

Pass

Fail

Protective earth resistance $\leq 300 \text{ m}\Omega$


Measured: $\text{m}\Omega$

The continuity test is to be performed between point **3** and the point **4**

Pass



Fail

Continuity test

0


MEASURING POINT DEPENDING OF THE INSTALLATION


Ensure proper calibration of your measuring equipment before EACH reading.

Please remove the identification sticker from the packaging and place it in one of the below window frame.


Cupola 1 with S/N

Cupola 2 with S/N

Cupola 3 with S/N

Note: if applicable, please staple the delivery note with this document.

Maintenance

- Replacement of : S/N :
- Replacement of : S/N :
- Replacement of : S/N :
- MAQUET batteries replaced

Remarks

Illumination verification

→ Minimum value : 40 000 lux, maximum value allowed : 160 000 lux

- Using the OPM077, measure the bulb voltage.
U = 23,5 ± 1V (AC + DC) ou U = 20,8 ± 1V (DC)

U =V

Cupola S/N < 20 000 → 80 000 lux

Cupola S/N > 20 000 → 100 000 lux

- N/A □ Pass □ Fail
- □ □
- □ □
- □ □

Cupola S/N :

Mesuread : lux

Cupola S/N :

Mesuread : lux

Cupola S/N :

Mesuread : lux

Switching over from ambient / main illumination

Setting

- Balance of the spring arm
- Vertical stop of the spring arm

 **A spring arm wears out in time and should be replaced every 6 years**

Replaced?
Yes No

- Test the equipment is easy to handle

Checking

- Hardness of the suspension by shaking the configuration
- Verticality of the suspension tube
- 6 fixing screws. If replaced, apply soft locktight glue.
- Fixing of the plastic covers on the spring arm
- Fixing and proper movement of the shutters
- Safety segment is positioned correctly (dismantle and lubricate if necessary)

 **A safety segment wears out in time and should be replaced every 2 years. Replaced ?**

Yes No

- Safety sleeve present and positioned with its fixing screw
- Assembly of the joints for each arm. Tightening of the screw and positioning of the retaining ring
- N/A Sterilizable handle holder
- Sterilizable handle engages correctly and stays in place
- Fixing of the outer handles
- Remove the fork bumper rubber cap and check the tightening of the circlip
- Limit stop on the yoke for the cupola rotation is functional
- N/A Keypad operates correctly
- System switches to battery mode and the ON/OFF button LED on the keypad turns orange and switches back to mains
- N/A Opening and closing of the lamp cover
- Replace the bulb cover
- Replace the bulb
- Light patch, adjusting the mirrors if necessary
- WPS is connected to the ground
- All boards are properly fixed inside the WPS
- N/A The correct setting on the regulator boards respective to the cupola
- Software version in START UP menu. If inferior to V5.40 → obsolete, inform the customer

Tightening

- Collar around the power supply connector at the ceiling tube
- Wires on the power supply connector at the ceiling tube
- Ceiling covers + proper positioning
- Wires on the connectors' terminal heads inside the WPS
- Wires on the battery pack
- All visible screws
- Presence of the stopper on the concerned suspension

N/A

Appearance

- All seals hold correctly and are not worn
- General appearance of the underside (no scratches, no cracks)
- No corrosion anywhere
- No paint chip anywhere
- Appearance of the keypad
- Attachment of the WPS on the wall
- Appearance of the WPS cover
- Holding of the cover and the hinges of the WPS
- Cleaning the complete configuration
- Check the label bearing the CE marking

N/A

We do not recommend the use of alcoholized solution

SERVICE PROTOCOL	
For use during Installation / Servicing	
<h1 style="margin: 0;">Axcel</h1>	<input type="checkbox"/> Contract Order n°..... <input type="checkbox"/> Chargeable <input type="checkbox"/> Other:.....
<input type="checkbox"/> Installation Delivery date: Installation date: Commissioning date:	
<input type="checkbox"/> Servicing Preventive / Curative Valid until:	
Hospital:	
Building:	Floor level:
Service:	Operating room:
Product(s) identification	
Cupola n°1	Ref. number: S/N:
Cupola n°2	Ref. number: S/N:
Cupola n°3	Ref. number: S/N:
Overall assessment	
<input type="checkbox"/> Configuration fully operational	
<input type="checkbox"/> Free of direct risk but deficiencies detected. May be corrected on short term.	
<input type="checkbox"/> Configuration shall not be used until all deficiencies are corrected <input type="checkbox"/> Configuration no longer safe. Taking out of service is recommended.	
Remarks	
.....	
Signature	
Name:	Name:
Date:	Function:
Installer/Technician authorized by MAQUET	Customer

Measuring equipment used

IEC 62353 compatible control units

Type: Voltmeter	S/N:	Calibration valid until:
Type: Luxmeter	S/N:	Calibration valid until:
Type: Megohmmeter	S/N:	Calibration valid until:
Type:	S/N:	Calibration valid until:


Electrical safety test

The protective earth resistance is to be measured between point **0** and the point **3**

Pass **Fail** Protective earth resistance $\leq 300 \text{ m}\Omega$ Measured: $\text{m}\Omega$


The continuity test is to be performed between point **3** and the point **4**

Pass **Fail** Continuity test

0


MEASURING POINT DEPENDING OF THE INSTALLATION


Ensure proper calibration of your measuring equipment before EACH reading.

Please remove the identification sticker from the packaging and place it in one of the below window frame.


Cupola 1 with S/N

Cupola 2 with S/N

Cupola 3 with S/N

Note: if applicable, please staple the delivery note with this document.

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Cupola S/N :

Mesuread : lux

Switching over from ambient / main illumination

Setting

- Balance of the spring arm
- Vertical stop of the spring arm

A spring arm wears out in time and should be replaced every 6 years

Replaced?

Yes No

- Test the equipment is easy to handle

Checking

- Hardness of the suspension by shaking the configuration
- Verticality of the suspension tube
- 6 fixing screws. If replaced, apply soft locktight glue.
- Fixing of the plastic covers on the spring arm
- Fixing and proper movement of the shutters
- Safety segment is positioned correctly
(dismantle and lubricate if necessary)

Yes No

A safety segment wears out in time and should be replaced every 2 years. Replaced ?

- Safety sleeve present and positioned with its fixing screw
- Assembly of the joints for each arm. Tightening of the screw and positioning of the retaining ring

N/A

- Sterilizable handle holder
- Sterilizable handle engages correctly and stays in place
- Fixing of the outer handles
- Remove the fork bumper rubber cap and check the tightening of the circlip
- Limit stop on the yoke for the cupola rotation is functional

N/A

- Keypad operates correctly
- System switches to battery mode and the ON/OFF button LED on the keypad turns orange and switches back to mains

N/A

- Opening and closing of the lamp cover
- Replace the bulb cover
- Replace the bulb
- Light patch, adjusting the mirrors if necessary
- WPS is connected to the ground
- All boards are properly fixed inside the WPS

N/A

- The correct setting on the regulator boards respective to the cupola
- Software version in START UP menu. If inferior to V5.40 → obsolete, inform the customer

Tightening

- Collar around the power supply connector at the ceiling tube
- Wires on the power supply connector at the ceiling tube
- Ceiling covers + proper positioning
- Wires on the connectors' terminal heads inside the WPS
- Wires on the battery pack
- All visible screws
- Presence of the stopper on the concerned suspension

Appearance

- All seals hold correctly and are not worn
- General appearance of the underside (no scratches, no cracks)
- No corrosion anywhere
- No paint chip anywhere
- Appearance of the keypad
- Attachment of the WPS on the wall
- Appearance of the WPS cover
- Holding of the cover and the hinges of the WPS
- Cleaning the complete configuration
- Check the label bearing the CE marking

We do not recommend the use of alcoholized solution

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