
Surgical light

Lámpara de quirófano

Technical Manual **EN**

Manual técnico **ES**

CONTENTS

SYMBOLS USED IN THE MANUAL	3
I. G8 EVOLUTION	4
QUALITY COMPLIANCE G8e.....	4
1 GENERAL CHARACTERISTICS G8e.....	5
2 CHARACTERISTICS OF THE SUSPENSION G8e	6
3 CLEANING/DISINFECTION/STERILIZATION G8e.....	8
4 MAINTENANCE G8e	10
4.1 Preventive maintenance.....	10
4.2 First level maintenance	10
4.3 Maintenance tools (OPM).....	12
4.4 Recommended spare parts and consumables.....	12
5 MAINTENANCE PROCEDURES G8e.....	13
5.1 Start-up test	13
5.2 Replacing a reflector	15
5.3 Removing the protective underside	17
5.4 Removing the cross member cover/accessing the electronics	17
5.5 Reinstalling the cross member cover	17
5.6 Removing the handle, centre reflector and optical core	17
6 G8e DIAGRAMS	18
6.1 G8e wiring	18
6.2 G8e central shaft	19
6.3 G8e cupola	20
6.4 G8e filter board.....	21
6.5 G8e PSU and microcontroller board	21
6.6 G8e control keypad	22
6.7 G8e motor-operated bulb holder	23
6.8 G8e fan module.....	23
7 DIAGRAMS AND PARTS LISTS	24
8 PROBLEMS AND TROUBLESHOOTING.....	36
APPENDIX 1	
G8 filter board.....	38
APPENDIX 2	
G8e PSU and microcontroller board	41
APPENDIX 3	
G8e control keypad control card	49
APPENDIX 4	
G8e control keypad card	50
APPENDIX 5	
G8e bulb holder starter board	51
APPENDIX 6	
G8e fan module.....	51

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I. G8 EVOLUTION

QUALITY STANDARDS COMPLIANCE

Certification of MAQUET SA quality system

AFAQ and G-MED certify that the system implemented by MAQUET for the design and production of surgical lights complies with the requirements of the following standards:

- ISO 9001 Version 2000
(AFAQ certificate N° QUAL/1989/32c)
- NF EN ISO 13485 Version 2004
(LNE/G-MED Certificate No. 1893/13485/1).

Standards

The G8 Evolution surgical light is designed to comply with the following standards :

- IEC 60601-1
- IEC 60601-2-41
- EN 60601-1-2

CE Marking/Intended use

Compliance with the requirements of Directive 93/42/EEC dated of June 14th 1993, relating to medical devices has been assessed in accordance with Annex VII of this Directive. The G8 Evolution surgical light is a class I device in accordance with Annex IX of Directive 93/42/EEC.

The range includes single, double or triple ceiling-mounted configurations, with or without a pre-installed video system (HANAUVISION). Possible options include a double fork, flat screen support and separate camera support.

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1 GENERAL CHARACTERISTICS OF G8 EVOLUTION SURGICAL LIGHTS

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* Factory set value

Note:

- The tolerated values are guaranteed on purchase of the product.
- The non-tolerated values are measured on a sample coming out of production.

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2 CHARACTERISTICS OF THE SUSPENSION G8e

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





3 CLEANING / DISINFECTION / STERILISATION G8e

Users must contact their hospital's sanitary specialists. The recommended products and procedures must be applied. Should there be any doubt concerning the compatibility of active agents to be used, contact the local MAQUET customer service.

3.1 Cleaning and disinfecting the surgical light

-  Check that the power is switched off and the light has cooled down before starting cleaning.

3.1.a General instructions concerning cleaning, disinfection and safety:

- Remove the sterilisable handles.
- Wipe the system with a cloth moistened with a surface cleaner. Follow the manufacturer's dilution and temperature recommendations.
- Rinse with a cloth moistened with water. Wipe dry.
- Wipe evenly with a cloth moistened with disinfectant. Follow the manufacturer's recommendations.
- Rinse with a cloth moistened with water to remove any residue (particularly when cleaners containing aldehydes, quaternary ammonium or surfactants are used).
- Wipe off with a dry cloth.
- Make sure that all liquid cleaning products used have been thoroughly wiped off.

4.1.2 Examples of recommended products

GETINGE USA products: TEC-QUAT 256.

ANIOS products: SURFA'SAFE; HEXANIOS G + R at 0.5%; ANIOSYME P.L.A.; SALVANIOS pH10; ANIOS DDSH.

Schülke & Mayr products: Antifect Plus

4.1.3 Examples of unsuitable products

-  Solutions containing glutaraldehyde, phenol, iodine, bleach, alcohol or chloride ions must not be used.

 Do not use disinfection by fumigation methods.

3.2 Cleaning and sterilising the handles

3.2.1 Before cleaning

- Use a soft cloth immediately after use to wipe away soiling from the handle surface.
- Store handles in a place that keeps them moist to make further cleaning easier.
- Take care to store them in such a way that the inside does not get soiled.

3.2.2 Cleaning

- Soak the handles in a detergent solution.¹
- Soak for 15 minutes to allow the solution to act, then clean by hand with a soft brush and a lint-free cloth.
- During cleaning, check regularly that the handles are fully clean and that no soiling remains on the inside or outside.
- If any soiling remains, repeat cleaning or use an ultrasonic cleaning process.
- Rinsing: Rinse thoroughly in clean water to completely eliminate the detergent solution.
- Drying: Wipe with a clean lint-free cloth.

¹ Never soak the handles in enzyme-based detergents as they may damage the material used; rinse thoroughly if these detergents are used.

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3.2.3 Disinfection

Handles may be disinfected by machine (Clean MAQUET) and rinsed at a maximum temperature of 93°C.

Typical recommended cycles:

Stage	Temperature	Time
Pre-washing	18 - 35° C	60 sec
Washing	46 - 50° C	300 sec
Neutralisation	41 - 43°C	30 sec
Washing 2	24 - 28°C	30 sec
Rinsing/Disinfection	92 - 93°C	600 sec
Drying		20 min

3.2.4 Sterilisation

After cleaning the handles must be steam sterilised as set out below:

Country	Sterilisation cycle	Temperature [°C]	Time [min]	Drying [min]
USA & Canada	Prevacuum ²	132 - 135	10	16
France	ATNC (Prion) (Prevacuum)	134	18	
Other countries	Prevacuum	Comply with national regulations		

- Check that each handle is clean before continuing the process.
- Wrap the handles with sterilisation wrapper material (double wrapper or equivalent).

Handles may also be placed in paper or plastic sterilisation bags³ for easier identification and reuse.

- Place the handles on steriliser trays⁴ with the opening downwards.
- Package with biological and/or chemical indicators for monitoring the sterilisation process. Follow applicable regulations.
- Run the sterilisation cycle according to the steriliser manufacturer's instructions.

 To ensure correct sterilisation do not allow any soiling to penetrate inside the handle.

 Handles are not guaranteed beyond 350 sterilisation cycles with the above sterilisation parameters.

- Dispose of the handles in the same way as other hazardous products in a hospital environment.

² This handle is made of a porous material.

³ Possible sterilisation bag suppliers:
Medical Action Industries
SBW Medical
Baxter International

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4.2.c Annual inspection (must be performed by an authorised technician)

- Check the bulb holder. Replace the bulbs.
- Check that the limit stops are in place on the flat screen support arms and cupolas.
- Clean the internal optical items (mirrors and lenses).
- Check that the light operates for one hour if the main power supply is cut (operation on backup supply).
- Replace the fan subassembly.

Safety items

Check the following points:

- Check that there are no cracks on the upper pivot of the DF Acrobat 3000 spring arm ref. **568101988** (if manufactured between 2001 and 2006).
- If there is the slightest sign of a crack, replace the spring arm.
- If no cracks are visible, install kit ref. 368104900.
- Attachment screws on suspension tube correctly tightened, seals in position.
- Arm(s) correctly assembled: safety of assembly.
- Limit stop segments in place on cupola attachment points. Disassemble cupolas and lubricate sleeve.
- Attachment of all covers and caps on cupolas and arms.
- All visible screws correctly tightened.
- Camera or dummy unit firmly attached: set screws tightly fastened inside cross member.

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Other checks

- Nominal light output level: see technical data.
- Earth continuity: max. 0.1 Ohm
- Verticality of the suspension tube; arms do not drift.
- Balancing system adjusted correctly: brakes tight and spring arms adjusted.
- Sterilisable handle locking mechanism.
- Adjustment of rotation limit stops.

Appearance:

Check the following points:

- General condition of arm covers.
- Condition of lenses.
- Condition of seals.
- Suspension system clean.

Dismantling certain elements may affect operation and safety.

Examples include:

- Servicing the electrical power supply,
- Servicing the suspension arm and balance system,
- Servicing the optical system of cupolas equipped with filters designed to eliminate radiation not visible to the patient. Surgical lights must never be used without these filters.

Contact the authorised MAQUET after-sales service department for this type of inspection.

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4.4 Recommended spare parts and consumables G8e

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Recommended stock of spare parts (for an installed base of 5 or 20 G8e cupolas)

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6 G8e MAINTENANCE PROCEDURES

6.1 Start-up test

- a) Each time the unit is turned on, the discharge bulb is fired 4 times then backup takes over (4 clicks are heard).
- b) If the starter is in working order, the bulb flashes with each click though it may not light up. The bulb may need to be changed if this happens.
- c) If the starter is worn out, even a new bulb will not necessarily flash with each click. Change the bulb holder if this happens.
- d) If you do not hear any clicks, the starter (bulb holder) is either defective or the control electronics are damaged.

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Fig. 6

5.2 Replacing a reflector

Contact an authorised MAQUET service centre to have this part replaced.

Figure 1

- a) Remove the segments from between the reflectors. Remove as many segments as necessary until there is sufficient space to remove and replace the reflector.

Figure 2

- b) Remove the outer seal (larger diameter) from around the reflector to be replaced.

Figure 3

- c) If the reflector to be replaced is adjacent to the camera fastener, remove the fastener's two M2x4 screws.

Figure 4

- d) Push the reflector away from the cupola. Using a flat, blunt, rigid tool, disengage it from its slot under the inner seal.

 Take care not to damage the inner seal while removing it.

Figure 5

- e) Lightly coat the inner seal with liquid soap to make installing the new reflector easier.
- f) Clear the empty space and position the new reflector against the outer stop. Slip it underneath the inner seal little by little with a flat, blunt, rigid tool.

 Take care not to damage the inner seal while removing it.

- g) Put the two M2x4 screws back in the camera fastener (if removed).

Figure 6

- h) Push the reflectors toward the inside of the cupola and put the outer seal back into place by hand.
- i) Apply a bead of silicone glue between the reflectors and insert the segments.

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Fig. 7

Fig. 8

5.3 Removing the protective underside

Figure 1

- 1) Remove the four M4x8 screws.

Figure 2

-  When putting the underside back in place, take care to put the two bottom and top seals back in place correctly so as not to damage the underside.

5.4 Removing the cross member cover/ accessing the electronics

Figure 3

- 2) Raise the hinged cover and hold it in place in order to avoid scratching the cross member cover.
- 3) Remove the eight M4x8 screws and the external tooth lock washers.

5.5 Reinstalling the cross member cover

Figure 4

- 4) Place the cross member cover on the cupola. Position the white seal between the protective underside and the cross member cover.
- 5) Lift the outer edges of the two black seals above the cross member cover. Position the cross member cover against the cross member and sandwich the visible ends of the black seals between.

Figure 5

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5.6 Removing the handle, centre reflector and optical core

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Removing the optical core:

Figure 8

- 11) Remove the two M8x0.75 nuts and the two M8 external tooth lock washers (via the top of the cupola).

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
6.1 G8e wiring

6.1 Cableado G8e


- 1 Central shaft (example of a three-arm shaft)
- 2 Brush holder
- 3 Main arm wiring
- 4 Spring arm wiring
- 5 Double fork wiring
- 6 Control keypad wiring

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* To ensure the backup bulb's stated performance levels.

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
- ① marrón — Hanauvision
- ② malva — Hanauvision
- ③ PE verde/amarillo
- ④ N negro — 27 V AC/DC*
- ⑤ L blanco — 27 V AC/DC*

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6.3 G8e cupola

6.3 Cúpula G8e


- 1** Fans and temperature control board (activation at 50°C)
368111555
- 2** PSU board and microcontroller board
368109555
- 3** Motor-operated bulb holder
368112998
- 4** Filter board
368108555
- 5** Hanauvision camera connector
- V 1-3 = 24 V power supply
- V 2 = earth
- V 5-6 = Autofocus presence
- V 11-12 = Data presence

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**6.4 G8e filter board
368108555**

**6.4 Tarjeta filtro G8e
368108555**

- 1 Four M4x8 screws for fastening board to cross member
- 2 Fan power supply connector (24 V)
- 3 Cupola power supply connector (24 V)
- 4 Cover open/closed detection connector
- 5 PSU and microcontroller power supply
- 6 Backup bulb power supply (24 V)
- 7 Motor power supply
- 8 Main bulb status signal to microcontroller board
- F1 (6.3A) protects the PSU board and the main bulb
- F2 (6.3A) protects the backup lamp
- F4 (1A) protects the camera


See the wiring diagrams in Appendix 1

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
- 1 Six M4x8 screws for fastening board to cross member
 - 2 Main lamp power supply
 - 3 Unused connector
 - 4 Main bulb status signal to filter board
 - 5 Control keypad connector
 - 6 PSU and microcontroller power supply (24 V)
- See the wiring diagrams in Appendix 2

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**7.6 G8e control keypad
368107998**

**7.6 Teclado de control G8e
368107998**


- 1 Control keypad connection
 - 2 Communication with microcontroller board
- See the wiring diagrams in Appendices 3 and 4

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**6.7 G8e motor-operated bulb holder
368112998**

**6.7 Soporte lámpara motorizado G8e
368112998**

- 1 Failover-to-backup mechanism
 - 2 Main bulb
 - 3 Backup bulb
 - 4 Main bulb starter
 - 5 Main bulb power supply (80 V ± 5)
 - 6 Motor (12 V) and backup bulb power supply (23.5 V ± 1)
- See the wiring diagram in Appendix 5

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

7 TROUBLESHOOTING

Anomaly	Likely cause	Corrective action
Cupola drifts	Suspension tube not vertical	Check tube verticality and ceiling structure
Cupola too stiff or does not stay in position	Balancing incorrectly adjusted	Adjust the balancing
	Ceiling structure unstable	Call MAQUET technical department
	Brakes not tight enough	Adjust the brakes
Light flickers (normal during cold start-up)	Poor contact, problems with electronics or ageing of main bulb	Call MAQUET technical department
Insufficient light output level	Main bulb faulty, electronic unit faulty, low light output selected	Check bulbs and replace if necessary, check bulb position, call technical department
Light beam not round	Incorrect type of bulb(s)	Check bulbs and replace if necessary
	Incorrect bulb or bulb holder position	Check bulb position, call technical department
Bulb service life too short	Non-compliant bulbs used	Check that the bulbs recommended by MAQUET are used
Lifecycle indicator fully lit even after replacing bulb	Lifecycle indicator has not been reset	Reset the indicator
Greenish light	Main bulb too old, main bulb faulty	Check bulbs and replace if necessary, check bulb position
	Light output level too high	Select a lower light output level
Cupola does not turn off	Off button not held down long enough	Press the Off button and hold for at least 2 seconds
	Cupola has just been turned on	Wait until the light reaches its service temperature (60 seconds) before turning it off
Cupola does not turn on	Power supply cut	Check power supply and fuses
	No bulb, bulb incorrectly installed	Check bulbs and replace if necessary, check that bulbs are correctly installed
	Electronic unit faulty	Call technical department
Cupola cover crazed or cracked	Cleaning solution too strong, unit cleaned or disinfected when hot	Use a different cleaning solution
Service life of sterilisable handles too short	Sterilisation procedure too harsh	Check sterilisation procedure
Sterilisable handles damaged or cracked	End of sterilisable handle service life	Replace sterilisable handles

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



Tarjeta filtro G8e Flyback






APPENDIX 1 D

G8e filter board
Halogen lamp control

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
APPENDIX 2 A

G8e PSU and microprocessor boards
Wiring schematic

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
APPENDIX 2 B

G8e PSU and microprocessor boards
Boost & Flyback

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
APPENDIX 2 C

G8e PSU and microprocessor boards
High-power section UPS

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
APPENDIX 2 D

G8e PSU and microprocessor boards
Flyback

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
APPENDIX 2 E

G8e PSU and microprocessor boards
Control diagram

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
APPENDIX 2 F

G8e PSU and microprocessor boards
Microcontroller

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
APPENDIX 2 G

G8e PSU and microprocessor boards
Control section UPS

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
APPENDIX 2 H

G8e PSU and microprocessor boards
Flyback and boost control

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
APPENDIX 2 I

G8e PSU and microprocessor boards
Control/protection

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



ANEXO 6

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II. G8

1 GENERAL CHARACTERISTICS OF G8 SURGICAL LIGHTS

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3 MAINTENANCE G8

3.1 Mechanical and functional checks

To preserve your surgical light's original performance and reliability, annual maintenance and inspections should be performed as follows:

- by a MAQUET technician or MAQUET-approved distributor during the guarantee period,
- by a MAQUET technician or MAQUET-approved distributor or by the hospital's technical maintenance department after the guarantee period expires,

3.1.a Daily inspection

- Check that the bulbs operate correctly.
- Check that the sterilisable handle clicks and locks in place correctly.
- Check that the arms move normally and do not sag.

3.1.b Monthly inspection

- Check that the light operates correctly from the backup power supply if a power cut were to occur.
- Check that the limit stops are in place.
- Test the fallback to the spare bulb.

3.1.c Annual inspection (must be performed by an authorised technician)

- Check the bulb holder. Replace the bulbs.
- Check that the limit stops are in place on the flat screen support arms and cupolas.
- Clean the internal optical items (mirrors and lenses).
- Check that the light operates for one hour if the main power supply is cut (operation on backup supply).
- Replace the fan subassembly every 2 years.

Safety items

Check the following points:

- Attachment screws on suspension tube correctly tightened, seals in position.
- Arm(s) correctly assembled: safety of assembly.
- Limit stop segments in place on cupola attachment points. Disassemble cupolas and lubricate sleeve.
- Attachment of all covers and caps on cupolas and arms.
- All visible screws correctly tightened.
- Camera or dummy unit firmly attached: set screws tightly fastened inside cross member.
- Check the two welding beads on the mounting lug for the spring arm balancing spring (see HHE location). If the onset of cracking is detected, the spring arm should be replaced. If the serial number corresponds to an affected part (see HHE table), kit ARD368104900 should be installed.

Other checks

- Nominal light output level: see technical data.
- Earth continuity: max. 0.1 Ohm
- Verticality of the suspension tube; arms do not drift.
- Balancing system adjusted correctly: brakes tight and spring arms adjusted.
- Sterilisable handle locking mechanism.

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- Adjustment of rotation limit stops.

Appearance:

Check the following points:

- General condition of arm covers.
- Condition of lenses.
- Condition of seals.

-  Suspension system clean.
-  Dismantling certain elements may affect operation and safety. Examples include:
 - Servicing the electrical power supply,
 - Servicing the suspension arm and balance system,
 - Servicing the optical system of cupolas equipped with filters designed to eliminate radiation not visible to the patient. Surgical lights must never be used without these filters.

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3.2 Cupola electrical checks

The dependability of this surgical light depends primarily on the quality of the connections inside the cupola.

Light on:

- Check the shape of the illuminated field: The discharge lamp should be operating and located at the centre of the optical core.
- Check the operation of the +/- keys for setting the four illumination levels. The indicator on the control panel should flash for approx. 60 seconds then be lit continuously.
- Check the level indicating the bulb service life.

Light off:

Remove the bulb cover:

- Check that the eyelets (1) are correctly attached to the bulb cover and tighten them using a wrench.
- Check that the bulb cover (2) fits perfectly in the main cover.

Remove the bulb holder:

- Check the position of the discharge bulb and reset if necessary (bad contacts may occur).
- Use a small screwdriver to open all the pins on connector (3).
- Manually check the operation of the locking system (magnet above the bulbs).
- Check the position of the bulbs on the bulb holder (discharge bulb and halogen bulb).
- Replace the bulb holder every two years.

Remove the main cover:

- Use two 7 mm wrenches to check that the female contact nuts (4) are correctly tightened.
- Open all the white connectors (5) and use a small screwdriver or sharp tool to open the eyelets.
- Visually check the condition of the bulb holder cables (check that nothing has burnt).
- Check the position of the MOS (6) which should be aligned

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vertically.

- Check the position and condition of the fuse (7).
- Manually check that the fan (8) turns freely. Replace it every two years.
- Check the position of the insulation sheet and cut it if it covers the fan.
- Check the position of the red reset button (9).



Fit the bulb cover:

- Turn the light on and check that the fan (8) operates quietly and regularly.
- Turn the light off and replace the main cover, bulb holder and bulb cover.

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4.3 Recommended spare parts and consumables G8

Recommended stock of spare parts (for an installed base of 5 or 20 G8 cupolas)

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5.1 Replacing a reflector

See pages 16 and 17

5.2 Removing the protective underside

Figure 1

- 1) Remove the four M4x8 screws.

Figure 2

When putting the underside back in place, take care to put the two bottom and top seals back in place correctly so as not to damage the underside.

5.3 Removing the cross member cover/ accessing the electronics

Figure 3

- 2) Open the lamp access hatch.
- 3) Remove the eight M4x8 securing screws from the cross member.

5.4 Reinstalling the cross member cover

Figure 4

- 4) Place the cross member cover on the cupola. Position the white seal between the protective underside and the cross member cover.

Figure 5

- 5) Lift the outer edges of the two black seals above the cross member cover. Position the cross member cover against the cross member and sandwich the visible ends of the black seals between.

Figure 6

- 6) Press on the cross member cover then lightly tighten the eight M4x8 screws, starting with the four near the underside.
- 7) Place the two black seals completely underneath the cross member.
- 8) Finish tightening the eight M4x8 screws then adjust the position of the white seal between the cross member cover and the protective underside (See Fig. 2).

5.5 Removing the handle, centre reflector and optical core

Removing the handle and the centre reflector:

Figure 7

- 9) Cut the autofocus button wires (two grey wires, pins 5 and 6 on the Sub-D15 connector inside the cross member).

 When putting the handle back in place, connect the two autofocus button wires to the Sub-D15 connector with a connecting block.

Figure 8

- 10) Release the handle by removing the four M5x10 screws.

Removing the optical core:

Figure 9

- 11) Remove the two M8x0.75 nuts and the two external tooth lock washers (via the top of the cupola).
- 12) Put the lamp access hatch back on the cross member (See Fig. 3).

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


6.3 G8 cupola

6.3 Cúpula del G8


- 1 Failover to backup card
- 2 Lifecycle LED reset
- 3 Bulb holder
- 4 Electronic control card 'A'
- 5 Safety switch
- 6 Fan
- 7 Electronic control card 'B'
- 8 Hanaurovision connector
- V 1-3 = 12 V power supply
- V 2 = earth
- V 5-6 = Autofocus presence
- V 11-12 = Data presence

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6.4 G8 electronic circuit board 'A'
568100188

6.4 Tarjeta electrónica "A" del G8
568100188


- 1 Main bulb power supply connector
- 2 Electronic circuit board 'B' connector
- 3 Camera connector, 24 V AC/DC
(centre pin not used)
- 4 Control keypad input voltage connector
(centre pin not used)
- 5 Extra connector, 24 V AC/DC
(centre pin not used)
- 6 Control keypad light intensity connector
- 7 Control keypad connector
- 8 Fan connector, 24 V DC
(3 centre pins not used)
- 9 Lifecycle LED reset connector

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6.5 G8 electronic circuit board 'B'
56078999

6.5. Tarjeta electrónica "B" del G8
56078999


- 1 Electronic circuit board 'A' connector
- 2 Input voltage connector, positions 1 and 3, 24 V AC/DC
- 3 Backup bulb connector, positions 2 and 4 (centre pin not used)
- 4 250 V fuse
- 5 Input voltage connector for backup bulb failover test 24 V AC/DC (centre pin not used)
- 6 Extra connector 24 V AC/DC (centre pin not used)
- 7 Backup bulb failover test connector (2 centre pins not used)
- 8 Light intensity input voltage connector 24 V AC/DC (centre pin not used)

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**6.6 Failover to backup card
56078959**

**6.6 Tarjeta de paso a emergencia
56078959**


- 1 Backup bulb connector
- 2 Input voltage connector, 24 V AC/DC
- 3 Backup bulb failover test connector

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**6.7 G8 control keypad
56078642**

**6.7 Teclado de control G8
56078642**


- 1 Input voltage connector, 24 V AC/DC
- 2 Control keypad connector
- 3 Electronic circuit board 'A' connector

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**6.8 G8 bulb holder
568100245**

**6.8 Soporte lámpara G8
568100245**

- 1 Main bulb
- 2 Backup bulb
- 3 Connector:
V 1 - 2 = Main bulb power supply
V 5 - 6 = Backup bulb power supply

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
To find your closest representative:

Please go to
www.maquet.com

Su contacto local:

Le rogamos nos visite en
www.maquet.com