

USER MANUAL - EN

MAQUET
GETINGE GROUP

HLX 3000

SURGICAL WORKPLACES



Copyright

All rights reserved. May not be duplicated, adapted or translated, in whole or in part, without prior written permission except where allowed by copyright law.

© Copyright Maquet S.A.

Subject to technical changes.

The illustrations and technical specifications provided in this manual may, on account of future product developments, differ slightly from the actual product supplied.

01.04.2010 | Ed2G



CONTENTS

| | | |
|--|--|----|
| Quality standards compliance | | 4 |
| Warnings | | 5 |
| Symbols used in this manual | | 6 |
| Symbols used on the product | | 6 |
| 1 Introduction | | 7 |
| 2 Description | | 9 |
| 2.1 Example of double configuration | | 9 |
| 2.2 HLX 3000 product range | | 10 |
| 3 Use | | 11 |
| 3.1 Energix WPS power supply unit | | 11 |
| 3.2 Optional LCD | | 12 |
| 3.3 Optional control via a remote PC | | 12 |
| 3.4 Optional backup power supply | | 13 |
| 4 Positioning | | 14 |
| 4.1 Moving the lighthead | | 14 |
| 4.2 Pre-positioning examples by speciality | | 15 |
| 4.3 Degrees of rotation | | 17 |
| 4.4 Fitting the sterilisable handle | | 18 |
| 5 Cleaning/Disinfection/Sterilisation | | 19 |
| 5.1 Cleaning and disinfecting the surgical light | | 19 |
| 5.2 Cleaning and sterilising the handles | | 20 |
| 6 Maintenance | | 22 |
| 6.1 First level maintenance | | 22 |
| 6.2 Annual maintenance | | 23 |
| 7 Replacing the bulbs | | 24 |
| 8 Accessories | | 25 |
| 9 General characteristics | | 26 |
| 10 EMC declaration | | 27 |
| 11 Troubleshooting | | 31 |

QUALITY STANDARDS COMPLIANCE

CERTIFICATION OF MAQUET SA'S QUALITY SYSTEM

LNE/G-MED certifies that the quality system developed by MAQUET SA for design, implementation, sales, installation and after-sales service of surgical lights complies with the requirements of the following international standards:

- ISO 9001:2008
- ISO 13485:2004

REFERENCE STANDARDS

The HLX 3000 surgical light is designed to comply with the following standards:

- EN 1041: 2008 Information supplied by the manufacturer with medical devices
- EN ISO 14971:2000/A1:2003 Medical devices - Application of risk management to medical devices
(ISO 14971:2000)
- EN 60601-1:1990 Medical electrical equipment - Part 1: General requirements for safety
Amendment A1:1993 to EN 60601-1:1990
Amendment A2:1995 to EN 60601-1:1990
- EN 60601-1-1:2001 Medical electrical equipment - Part 1-1: General requirements for safety
- Collateral standard: safety requirements for medical electrical systems
- EN 60601-1-2:2007 Medical electrical equipment - Part 1-2: General requirements for safety
- Collateral standard: Electromagnetic compatibility - Requirements and tests
- EN 60601-1-4:1996 Medical electrical equipment - Part 1-4: General requirements for safety
- Collateral standard: Programmable electrical medical systems
Amendment A1:1999 to EN 60601-1-4:1996
- EN 60601-1-6:2004 Medical electrical equipment - Part 1-6: General requirements for safety
- Collateral standard: Usability
- EN 60601-2-41:2000 Medical electrical equipment - Part 2-41: Particular requirements for the safety of surgical luminaires and luminaires for diagnosis
- EN 980:2008 Graphical symbols for use in the labelling of medical devices
- FCC part 15 Radiofrequency device

This product has been verified for compliance with the following additional standards:

- CAN/CSA-C22.2 No. 601.1-M90 (R2005)
(includes national differences for Canada),
- EN 60601-1:1990 + A1:1993 + A2:1995
- UL 60601-1, first edition, 2006-04-26
(includes national differences for the USA)

This device does not meet FCC/UL requirements unless it is equipped with a MAQUET SA electrical power supply.

CE MARKING

Compliance with the requirements of European Directive 93/42/EEC dated 14 June 1993 relating to medical devices was assessed as described in Annex VII of the Directive. HLX 3000 surgical lights belong to Class I as described in Annex IX of Directive 93/42/EEC.

WARNINGS



WARNING

Changes or modifications not expressly approved by MAQUET SA could void the user's authority to operate the equipment.



WARNING

Light is a form of energy that can dry out tissue, particularly if light beams from more than one lighthouse are superimposed. Users must be vigilant and set appropriate illumination levels for each operation and patient, in particular for long operations.



WARNING

Light is a form of energy that, on account of certain wavelengths emitted, may not be suitable for certain pathologies.



WARNING

The surgical light is a high-intensity light source. Do not look directly into it.



WARNING

Do not use the device in the presence of flammable anaesthetic gases.



WARNING





Do not use in an MRI environment.










WARNING

To avoid all risks of asepsis, never touch the control keypad on the fork during a surgical procedure unless it is protected by a single-use cover.

SYMBOLS USED IN THIS MANUAL

| Symbols | Meaning |
|---|--|
|  | Mandatory May affect patient or user safety |
|  | Recommendation Risk of damage to device or accessories |
|  | CE label The device complies with the requirements of European directive 93/42/EEC relating to medical devices. |
|  | Medical Equipment Classified with respect to electric shock, fire and mechanical hazards only in accordance with UL 60601-1, IEC 60601-2-41 and CSA C22.2 No. 601-M90. |

SYMBOLS USED ON THE PRODUCT

| Symbols | Meaning |
|---|--|
|  | Caution Read the documentation for the unit thoroughly. |
|  | Alternating current |
|  | Direct current |
| REF. SN. | Product technical description and serial numbers. |
| $22.8V_{EFF} (V_{RMS}^{AC+DC})$ | Rectified true RMS voltage on bulb terminals |
|  | Hot surface |
|  | Medical Equipment Classified with respect to electric shock, fire and mechanical hazards only in accordance with UL 60601-1, IEC 60601-2-41 and CSA C22.2 No. 601-M90. |
|  | CE label The device complies with the requirements of European directive 93/42/EEC relating to medical devices. |
|  | Do not dispose of this device as unsorted municipal waste. Take it to a collection facility for value enhancement, recycling or re-use. |

1 INTRODUCTION

1.1 INTENDED USE

- The HLX 3000 range is designed for use in medical applications to illuminate patients' bodies during operations, diagnosis or treatment.
- The HLX 3000 range consists of single-arm, dual-arm or triple-arm ceiling-mounted lights with single or double forks. Some models are prewired for video (HANAUVISION). Possible options: comfort bracket, flat screen holder, separate camera holder, dimmer.

1.2 SPECIAL FEATURES

- Excellent shadow dilution
- Wide range of movement
- Compatible with laminar flows
- Upgraded multimedia
- Excellent adaptation of the operating field to all types of surgery.

1.3 INAPPROPRIATE USE

- The unit is intended only for use as mentioned in the user manual and is suitable only for such use. Any other use may result in danger to users and/or damage to the device or its environment.
- For safety reasons, do not alter or adapt the device in any way without MAQUET SA's approval.

1.4 INTENDED USERS

- The device may be operated only by trained medical staff.
- The device must be cleaned by specialised personnel.

1.5 INTENDED EQUIPMENT

The models named in this document may be equipped with other manufacturers' end appliances (e.g. flat screens). For information on operating them, please refer to each manufacturer's user manuals.

1.6 ENVIRONMENTAL REQUIREMENTS

TRANSPORT AND STORAGE

- Ambient temperature: -10°C to 60°C
- Relative humidity: 20-75%
- Atmospheric pressure: 500 to 1060 hPa

OPERATION

- Ambient temperature: 10°C to 40°C
- Relative humidity: 20-75%
- Atmospheric pressure: 700 to 1060 hPa
- The HLX 3000 system must be installed and commissioned in accordance with the EMC information provided in chapter 10.
- Portable and mobile RF communication devices may affect the correct operation of this equipment.

1.7 MAINTENANCE

MAQUET SA's warranty and the safety and integrity of the operation of the product are guaranteed only if:

- All inspection, maintenance and repair operations are performed by a MAQUET engineer or a trained and authorised technical support technician.
- Only original accessories, consumables and spare parts are used.
- Maintenance operations and checks are performed and documented at least once a year as specified by the preventive maintenance programme.

1.8 USER MANUAL

- This user manual is an integral component of the equipment.
- It must be kept near the equipment for future reference.
- It must be read in full before attempting to use the equipment.

1.9 FCC PART 15

(USA only) This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at personal expense.

1.10 MANUFACTURED BY

MAQUET SA
 Parc de Limère
 Avenue de la Pomme de Pin
 CS 10008 ARDON
 45074 ORLÉANS CEDEX 2
 FRANCE
 Telephone: +33 (0) 2 38 25 88 88
 Fax: +33 (0) 2 38 25 88 00
www.maquet.com

2 DESCRIPTION

2.1 EXAMPLE OF DOUBLE CONFIGURATION





Figure 1


- 1 Main arm
- 2 Spring arm
- 3 Comfort bracket
- 4 Fresnel lens
- 5 Sterilisable handle
- 6 Side positioning handle
- 7 Underside
- 8 Shell
- 9 Cover
- 10 Lamp holder
- 11 Optional dimmer
- 12 Power on/off
- 13 Keypad used to increase lighting level
- 14 Keypad used to decrease lighting level
- 15 LED-indicator of the lighting level
- 16 Indicator of normal operation

2.2 HLX 3000 PRODUCT RANGE


- The HLX 3000 surgical lighting system is available in three different versions:



HLX 3003




HLX 3005









HLX 3004 TV
pre-wired for a HANAUVISION camera

3 USE

3.1 ENERGIX WPS POWER SUPPLY UNIT




The light is turned on and off via the ENERGIX WPS power supply units.


| | | |
|---|--------------------------------|---|
|  | <p>On/Off button</p> | <ul style="list-style-type: none"> ■ Turns lighthouse on and off with a single press. ■ Lighthouse lights up gradually to the last intensity value stored in its memory. ■ Factory set to 100% illumination. |
|  | <p>LED off</p> | <ul style="list-style-type: none"> ■ Power supply off |
|  | <p>Green LED</p> | <ul style="list-style-type: none"> ■ Power supply on |
|  | <p>Red LED</p> | <ul style="list-style-type: none"> ■ Backup power supply on (WPS XX1 only) |
|  | <p>Illumination adjustment</p> | <ul style="list-style-type: none"> ■ To change levels, press once or press and hold. ■ 8 illumination levels (2 per LED) |
|  | <p>Flashing yellow LED</p> | <ul style="list-style-type: none"> ■ Bulb failure or ■ Ambient light |

 **WARNING**

Do not store objects on the wall power supply units.



3.2 OPTIONAL LCD



The LCD is used to:

- Check the service life of consumables (bulbs, bulb holder, batteries),
- Perform routine tests (backup tests),
- Troubleshoot malfunctions.

3.3 OPTIONAL CONTROL VIA A REMOTE PC

| | | |
|--|--|--|
| | <ul style="list-style-type: none"> ■ Light turned on via the ENERGIX unit | <ul style="list-style-type: none"> ■ Link not operational ■ Settings are adjusted via the ENERGIX unit |
| | <ul style="list-style-type: none"> ■ Light turned on via the PC | <ul style="list-style-type: none"> ■ Operational link ■ Settings are adjusted via the PC <div style="background-color: #0056b3; color: white; padding: 5px; margin-top: 10px;"> </div> <p>The link will be cut if adjustments are made via the ENERGIX unit.</p> |

3.4 OPTIONAL BACKUP POWER SUPPLY



WARNING

In the event of a power cut, only lightheads connected to a backup power supply will remain on.

- The power supplies may be optionally equipped for use with your operating room's 24 V backup power supply.
- If the mains supply is restored after being cut, the light turns on in surgical light mode.

| | | |
|---|---|--|
| | <p>Red LED</p> | <p>Backup power supply on (WPS XX1 only)</p> |
| <p>Battery backup test (daily check)</p> | | |
| | <ul style="list-style-type: none"> ■ Turn on the lighthead. ■ Press and hold the On/Off button for 3 seconds. | |
| | <ul style="list-style-type: none"> ■ The backup batteries take over. ■ The LED turns from green to red. ■ The lighthead automatically switches back to the mains supply 3 seconds later. | |
| <p>Battery capacity test (monthly check)</p> | | |
| | <ul style="list-style-type: none"> ■ Press and hold the «+», «-» buttons and the On/Off button for 2 seconds. | |
| | <ul style="list-style-type: none"> ■ The backup batteries take over. ■ Allow the batteries to discharge (one hour per lighthead). | |
| <p>BATTERY GOOD</p> | <ul style="list-style-type: none"> ■ Battery charge sufficient. | |
| <p>BATTERY BAD</p> | <ul style="list-style-type: none"> ■ Replace the batteries. | |
| | <ul style="list-style-type: none"> ■ Press the On/Off button to interrupt the test or to switch back to the mains supply. | |

4 POSITIONING

4.1 MOVING THE LIGHTHEAD

- The light should be positioned prior to any procedures to avoid having to move it more than necessary later on. Correctly positioning the light for each operation will limit the chances of it coming into contact with other objects (IV pole, pendant, etc.).




Figure 1

- Recommended distance between the underside and the operating field: 1 m.




Figure 2

- Use the removable handle at the centre or the handle on the side to position the lighthouse.




Figure 3

- The rotation of the handle allows to focus the light according to the distance to the operating field.



RECOMMENDATION





- Do not turn the lighthouse towards the ceiling. It will automatically switch off.
- Do not use the light suspension system to carry or raise objects.
- Do not hang from the light.






WARNING

To avoid all risk of asepsis, only the medical team may touch the sterilisable handle during procedures.


4.2 PRE-POSITIONING EXAMPLES BY SPECIALITY

| Pre-positioning examples | Surgical specialities |
|---|---|
|  <p data-bbox="327 705 403 730">Figure 1</p> | <p data-bbox="911 526 1007 595">General surgery</p> |
|  <p data-bbox="327 1084 403 1108">Figure 2</p> | <p data-bbox="911 875 1129 1014">Urology, organ transplants, gynaecology, delivery</p> |
|  <p data-bbox="327 1471 403 1496">Figure 3</p> | <p data-bbox="911 1274 1362 1379">General surgery, abdominal surgery, digestive surgery, thoracic surgery</p> |
|  <p data-bbox="327 1848 403 1872">Figure 4</p> | <p data-bbox="911 1695 1038 1727">Proctology</p> |

| Pre-positioning examples | Surgical specialities |
|---|--|
|  <p data-bbox="225 611 304 638">Figure 5</p> | <p data-bbox="810 450 979 483">Neurosurgery</p> |
|  <p data-bbox="225 992 304 1019">Figure 6</p> | <p data-bbox="810 779 1166 925">Plastic and reconstructive surgery, face transplants, oral and maxillofacial surgery</p> |
|  <p data-bbox="225 1382 304 1408">Figure 7</p> | <p data-bbox="810 1182 1002 1290">Otolaryngology, ophthalmology, dermatology</p> |


Note: The positions illustrated are given for guidance only. Each operator will choose the correct position of the light for his/her work habits.

4.3 DEGREES OF ROTATION



Double fork version

- Unlimited rotation on all axes, except on the lighthouse (stop at 420 °)



Single fork version

4.4 FITTING THE STERILISABLE HANDLE



WARNING

Check that the sterilisable/disposable handle is compatible with the product.

Before reusing the handle after it has been sterilised, always:

- Check for cracks.
- Check that the handle clicks into place correctly in the surgical light.




Figure 1

Fitting the sterilisable handle

- Insert the handle into the mount until it clicks into place.

Removing the sterilisable handle

- Press down on the push button while removing the handle.

5 CLEANING / DISINFECTION / STERILISATION

Users must contact their hospital's sanitary specialists. The recommended products and procedures must be applied. Should there be any doubt concerning the compatibility of active agents to be used, contact the local MAQUET customer service department.

5.1 CLEANING AND DISINFECTING THE SURGICAL LIGHT



RECOMMENDATION

Check that the power is switched off and the light has cooled down before starting cleaning.

GENERAL INSTRUCTIONS CONCERNING CLEANING, DISINFECTION AND SAFETY

- Remove the sterilisable handles.
- Wipe the equipment with a cloth moistened with a surface cleaner. Follow the manufacturer's dilution, application time, and temperature recommendations.
- Rinse the unit with a cloth and clean water. Wipe dry.
- Wipe evenly with a cloth moistened with disinfectant. Follow the manufacturer's recommendations.
- Use a cloth to rinse with clean water in order to remove residues (in particular products containing aldehydes, quaternary ammonium or surfactants).
- Wipe with a dry cloth.
- Make sure no liquid residue is left on the device after cleaning.

RECOMMENDED PRODUCTS

Getinge USA product: TEC-QUAT 256.

Anios products: SURFA'SAFE; 0.5% HEXANIOS G + R; ANIOSYME P.L.A.; SALVANIOS pH 10; ANIOS DDSH.

Schülke & Mayr products: ANTIFECT PLUS.

PROHIBITED PRODUCTS



WARNING

- Solutions containing glutaraldehyde, phenol, iodine, bleach, alcohol or chloride ions must not be used.
- Fumigation methods are unsuitable for disinfecting the unit and must not be used.

5.2 CLEANING AND STERILISING THE HANDLES

BEFORE CLEANING

- Use a soft cloth immediately after use to wipe away soiling from the handle surface.
- Store handles in a place that keeps them moist to make further cleaning easier.
- Take care to store them in such a way that the inside does not get soiled.

CLEANING

- Immerse the handles in a detergent solution.¹
- Soak for 15 minutes to allow the solution to act, then clean by hand with a soft brush and a lint-free cloth.
- During cleaning, check regularly that the handles are fully clean and that no soiling remains on the inside or outside.
- If any soiling remains, repeat cleaning or use an ultrasonic cleaning process.
- Rinsing: Rinse thoroughly in clean water to completely eliminate the detergent solution.
- Drying: Wipe with a clean lint-free cloth.

DISINFECTION

Handles may be disinfected by machine (e.g. Getinge) and rinsed at a maximum temperature of 93°C.

Typical recommended cycles:

| Step | Temperature | Time |
|----------------|-------------|--------|
| Pre-wash | 18-35°C | 60 sec |
| Wash | 46-50°C | 5 min |
| Neutralisation | 41-43°C | 30 sec |
| Wash 2 | 24 - 28°C | 30 sec |
| Rinse | 92-93°C | 10 min |
| Dry | | 20 min |

¹ A non-enzyme-based detergent is recommended. Enzymatic detergents may damage the handles. Never soak the handles in these detergents for prolonged periods. Rinse thoroughly.

STERILISATION

After cleaning, the handles must be steam sterilised as set out below:

| Countries | Sterilisation cycle | Temperature [°C] | Time [min] | Drying [min] |
|-----------------|-----------------------------|----------------------------------|------------|--------------|
| USA & Canada | Prevacuum ¹ | 132 - 135 | 10 | 16 |
| France | ATNC (Prion) (Prevacuum) | 134 | 18 | |
| Other countries | Prevacuum ¹ | Comply with national regulations | | |

- Check that each handle is clean before continuing the process.
- Wrap the handles with sterilisation wrapper material (double wrapper or equivalent). They may also be placed in paper or plastic sterilisation bags², for easier identification and reuse.
- Place the handles on steriliser trays with the opening downwards.³
- Package with biological and/or chemical indicators for monitoring the sterilisation process, in accordance with applicable regulations.
- Run the sterilisation cycle according to the steriliser manufacturer's instructions.



RECOMMENDATION

- To ensure correct sterilisation do not allow any soiling to penetrate inside the handle.
- Handles are guaranteed for no more than 350 sterilisation cycles when the above sterilisation parameters are used.
- Dispose of sterilisable handles in the same way as other hazardous products in a hospital environment.

¹ This handle is made of a porous material.

² Possible sterilisation bag suppliers :

Medical Action Industries









SBW Medical

Baxter International

³ For air removal and faster drying.

6 MAINTENANCE

6.1 FIRST LEVEL MAINTENANCE

| Daily checks (user) | |
|---|---|
| Solutions are listed in Section 11 | |
|  | <ul style="list-style-type: none"> ■ Check the lighthouse heads for chipped paint, impact marks and any other damage. |
|  | <ul style="list-style-type: none"> ■ Check the stability/drift of the main arms and the spring arms. |
|  | <ul style="list-style-type: none"> ■ Check that the spring arm remains in position. ■ Three checkpoints: bottom, middle, top. |
|  | <ul style="list-style-type: none"> ■ Check whether the sterilisable handle clicks and locks in place correctly; replace it if not. |
|  | <ul style="list-style-type: none"> ■ Check that the bulbs operate correctly. |
|  | <ul style="list-style-type: none"> ■ Check whether the backup power supply turns on and the light operates correctly if a power cut were to occur (page 13). |
| Monthly checks | |
| (HLX 3000 power supply units with battery backup only) | |
|  | <ul style="list-style-type: none"> ■ Check the lighthouse heads for chipped paint, impact marks and any other damage. |
|  | <ul style="list-style-type: none"> ■ Check whether the backup power supply turns on and the light operates correctly if a power cut were to occur (red LED). ■ Check the capacity of the batteries (see page 13). |

6.2 ANNUAL MAINTENANCE

To preserve your surgical light's original performance and reliability and ensure its safety, annual maintenance and inspections should be performed by:

- a MAQUET technician or
- an authorised MAQUET distributor, or
- a hospital technician trained in servicing MAQUET SA surgical lights.

Note: To find the training solution that best matches your needs, please contact your nearest MAQUET office.



RECOMMENDATION

- We strongly recommend taking out a MAQUET maintenance contract for all maintenance work on the light.
-



WARNING

The operation and safety of the device may be affected by the removal of certain components during servicing operations. For example:

- When servicing the electrical power supply.
- When servicing the suspension arm and balance system.
- Servicing the optical system of cupolas equipped with filters designed to eliminate radiation not visible to the patient. Surgical lights must never be used without these filters.

Contact the authorised MAQUET after-sales service department for this type of inspection.

7 REPLACING THE BULB




WARNING

- Do not change the light bulb while operating.
- The unit must operate with the cover closed. When performing maintenance, beware of hot surfaces which are indicated by the following icon:




RECOMMENDATION


- Replace bulbs every 1000 hours.
- Switch off the power supply and allow the lighthouse to cool for 15 minutes.
- Only use MAQUET XIR 001 bulbs.
- Handle the bulb with care using a clean, dry cloth.
- Never touch light bulb with bare hands. Grease on bulbs can shorten their life or break them.




- Press the button and lift the upper cover.




- Open the two clasps of the lamp holder.




- Remove the lamp holder.
- Remove the bulb.




- Insert the new bulb fully in the bulb holder.




- Check that the bulb is correctly seated.



- Reinstall the bulb holder in the lighthouse.






- Close the clasps.



- Close the cover.

8 ACCESSORIES

| Accessories | Code | Part number | |
|---|---|--------------|---------|
|  | Set of 3 sterilisable handles | 5 679 17 999 | HLX 005 |
|  | Pack of 1 spare bulb IRC 40 W - 22,8 V | 5 679 01 947 | XIR 001 |
|  | Pack of 10 spare bulbs IRC 40 W - 22,8 V | 5 679 01 948 | XIR 010 |

9 GENERAL CHARACTERISTICS

(IN ACCORDANCE WITH STANDARD IEC 60601-2-41 AND IEC 60601-1)

| Specifications | Unit | HLX 3003 | HLX 3005 | HLX 3004TV | |
|---|--|----------------------|------------|------------|----|
| Nominal illumination* (Ec) | lx ± 20% | 80,000 | 140,000 | 100,000 | |
| Diameter d10 | cm (inch) | 16.5 (6.5) | 16.5 (6.5) | 16.5 (6.5) | |
| Diameter d50 | cm (inch) | 9 (3.5) | 9 (3.5) | 9 (3.5) | |
| Illumination depth | cm (inch) | 93 (36) | 148 (58.2) | 71 (27.9) | |
| Colour temperature** (Ra) | K ± 10% | 4,300 | 4,300 | 4,300 | |
| Colour rendering index | NA | 93 | 93 | 93 | |
| Dilution des ombres | With one mask | % | 40 | 65 | 87 |
| | With two masks | % | 45 | 55 | 51 |
| | At base of tube | % | 99 | 93 | 97 |
| | With one mask, at base of tube | % | 39 | 59 | 81 |
| | With two masks, at base of tube | % | 44 | 50 | 45 |
| Radiant energy | mW.m ⁻² . lx ⁻¹ | 4.5 | 4.5 | 4.5 | |
| Irradiance (Ee) | W/m ² | 360 | 495 | 450 | |
| Electrical classification | | Class I | | | |
| Degree of protection against harmful ingress of water | - | Ordinaire | | | |
| Degree of protection against electric shock | | No applied parts | | | |
| Methods of sterilization or disinfection | - | See chapter 5 | | | |
| Mode of operation | - | Continuous operation | | | |

Note:

* Measured values: HLX 3003:73,100 lx; HLX 3005: 116,000 lx; HLX 3004TV: 87,400 lx

** Measured values: HLX 3003: 4,143 K; HLX 3005: 4,148 K; HLX 3004TV: 4,128 K

10

EMC DECLARATION

(IN ACCORDANCE WITH STANDARD EN 60601-1-2: 2001)

Table 201 - Guidance and manufacturer's declaration — electromagnetic emissions

The HLX 3000 system is intended for use in the electromagnetic environment specified below. The customer or the user of the HLX 3000 system should ensure that it is used in such an environment.

| Immunity test | Conformity | Electromagnetic environment – guidance |
|---|-------------------|--|
| Radio-frequency emissions CISPR 11 | Group 1 | The HLX 3000 system uses radio-frequency energy for its internal functions only. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| Radio-frequency emissions CISPR 11 | Class A | The HLX 3000 system is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network which supplies buildings used for domestic purposes. |
| Harmonic current emissions }IEC 61000-3-2 | N/A | |
| Voltage fluctuations and flicker IEC 61000-3-3 | N/A | |

Table 202 - Guidance and manufacturer's declaration — electromagnetic immunity


The HLX 3000 system is intended for use in the electromagnetic environment specified below. The customer or the user of the HLX 3000 system should ensure that it is used in such an environment.

| Immunity test | IEC 60601 Test level | Compliance level | Electromagnetic environment – guidance |
|---|---|---|--|
| Electrostatic discharge (DES) IEC 61000-4-2 | ± 6 kV contact ± 8 kV air | ± 6 kV contact ± 8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient/burst IEC 61000-4-4 | ± 2 kV for power supply lines ± 1 kV for input/output lines | ± 2 kV for power supply lines ± 1 kV for input/output lines | Mains power quality should be that of a typical commercial or hospital environment. |
| Surges IEC 61000-4-5 | ± 1 kV differential mode ± 2 kV common mode | ± 1 kV differential mode ± 2 kV common mode | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | < 5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (95% dip in U_T) for 5 cycles | < 5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (95% dip in U_T) for 5 cycles | Mains power quality should be that of a typical commercial or hospital environment. If the user of the HLX 3000 system requires continued operation during power mains interruptions, it is recommended that the HLX 3000 system be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 3 A/m | 3 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |

Note: V_T is the AC mains voltage prior to application of the test level.

Table 204 - Guidance and manufacturer's declaration — electromagnetic immunity

The HLX 3000 system is intended for use in the electromagnetic environment specified below. The customer or the user of the HLX 3000 system should ensure that it is used in such an environment.

| Immunity test | IEC 60601 Test level | Compliance level | Electromagnetic environment – guidance |
|---|---|---------------------|--|
| Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3 | 3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz | 3 Vrms 3 V/m | <p>Portable and mobile RF communications equipment should be used no closer to any part of the HLX 3000 system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = [3.5 / I_{T1}] \sqrt{P} = 1,17 \sqrt{P}$ $d = [3.5 / E_{E1}] \sqrt{P} \quad 80 \text{ MHz} - 800 \text{ MHz} = 1,17 \sqrt{P}$ $d = [7 / E_{E1}] \sqrt{P} \quad 800 \text{ MHz} - 2,5 \text{ GHz} = 2,34 \sqrt{P}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the HLX 3000 system is used exceeds the applicable RF compliance level above, the HLX 3000 system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the HLX 3000 system.

^b Over the frequency range of 150 kHz to 80 MHz, the field strength should be less than 3 V/m.

Table 206 - Recommended separation distances between portable and mobile RF communications equipment and the HLX 3000 system

The HLX 3000 system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the HLX 3000 system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the HLX 3000 system as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output of transmitter W | Separation distance according to frequency of transmitter m | | |
|--|--|--|---|
| | 150 kHz to 80 MHz $d = [1.17] \sqrt{P}$ | 80 MHz to 800 MHz $d = [1.17] \sqrt{P}$ | 800 MHz to 2.5 GHz $d = [2.34] \sqrt{P}$ |
| 0.01 | 0.12 | 0.12 | 0.24 |
| 0.1 | 0.37 | 0.37 | 0.74 |
| 1 | 1.17 | 1.17 | 2.34 |
| 10 | 3.69 | 3.69 | 7.38 |
| 100 | 11.70 | 11.70 | 23.40 |

For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

11 TROUBLESHOOTING

| Problem | Likely cause | Corrective action |
|---|--|--|
| <ul style="list-style-type: none"> ■ The lighthouse does not turn on | <ul style="list-style-type: none"> ■ Faulty or missing bulb | <ul style="list-style-type: none"> ■ Cut off the power supply and replace the bulb(s) |
| | <ul style="list-style-type: none"> ■ Power cut | <ul style="list-style-type: none"> ■ Check if other devices on the same supply network are operating |
| | <ul style="list-style-type: none"> ■ Other reason | <ul style="list-style-type: none"> ■ Call MAQUET technical department |
| <ul style="list-style-type: none"> ■ None of the lighthouses light up | <ul style="list-style-type: none"> ■ Each lighthouse has its own control | <ul style="list-style-type: none"> ■ Check the LED on the keypad on each lighthouse (green LED) |
| <ul style="list-style-type: none"> ■ Bulb service life too short | <ul style="list-style-type: none"> ■ Improper bulb or overvoltage condition | <ul style="list-style-type: none"> ■ Make sure that you are using the required MAQUET bulbs exclusively ■ Check voltage at the bulb pins: $22.8V_{EFF} (V_{AC+DC}^{RMS})$ |
| <ul style="list-style-type: none"> ■ The sterilisable handle does not click into place correctly | <ul style="list-style-type: none"> ■ Sterilisation parameters (temperature, time) exceeded | <ul style="list-style-type: none"> ■ Check whether the locking mechanism operates correctly (audible click) and check the entire handle |
| | <ul style="list-style-type: none"> ■ Its maximum service life has expired or the handle is twisted or bent. | <ul style="list-style-type: none"> ■ Replace the handle |
| <ul style="list-style-type: none"> ■ The lighthouse drifts | <ul style="list-style-type: none"> ■ Suspension tube not vertical | <ul style="list-style-type: none"> ■ Check tube verticality and ceiling structure |
| | <ul style="list-style-type: none"> ■ Ceiling structure unstable | <ul style="list-style-type: none"> ■ Call MAQUET technical department |
| | <ul style="list-style-type: none"> ■ Brake incorrectly adjusted | <ul style="list-style-type: none"> ■ Adjust the brake |
| <ul style="list-style-type: none"> ■ Lighthouse too flexible or too rigid to manoeuvre | <ul style="list-style-type: none"> ■ Brake incorrectly adjusted | <ul style="list-style-type: none"> ■ Adjust the brake screw |
| | <ul style="list-style-type: none"> ■ Insufficient lubrication | <ul style="list-style-type: none"> ■ Call MAQUET technical department |

MAQUET

GETINGE GROUP

MAQUET SA
Parc de Limère
Avenue de la Pomme de Pin
CS 10008 ARDON
45074 ORLÉANS CEDEX 2, France
Telephone: +33 (0) 2 38 25 88 88
Fax: +33 (0) 2 38 25 88 00

To find your closest representative:
Please go to
www.maquet.com

The GETINGE Group is a leading global provider of equipment and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. Equipment, service and technologies are supplied under the brands ArjoHuntleigh for patient handling and hygiene, disinfection, DVT prevention, medical beds, therapeutic surfaces and diagnostics, GENTINGE for infection control and prevention within Healthcare and Life Science and MAQUET for Surgical Workplaces, Cardiovascular and Critical Care.

MAQUET® is a registered trademark of MAQUET S.A GmbH & Co. KG. The MAQUET specifications provided herein are for guidance only. MAQUET reserves the right to change them without notice. • User manual • Reference: 0132102 EN Ed 2G •