

Hemodynamic Normal Values


		Central Venous Oxygenation - Oxygenation Balance (Oxygen load of the venous blood after passing through the organs)	ScvO ₂ *	70-80 %	
		O ₂ Consumption (Consumption of O ₂ by organs)	VO ₂ I	125-175 ml/min/m ²	
Oxygen Delivery	O ₂ Delivery (Delivery of O ₂ via blood to organs)		DO ₂ I	400-650 ml/min/m ²	
	Hemoglobin (Oxygen transporter in blood)		Hb**	8.7-11.2 mmol/l (Male) 7.5-9.9 mmol/l (Female)	
	Arterial / capillary oxygen saturation (Oxygen load of arterial blood)		SaO ₂ /SpO ₂	96-100 %	
Blood Flow	Flow	Cardiac Index (Trend, Cal, td, PC)	CI	3.0-5.0 l/min/m ²	
	Chronotropy	Heart Rate/Pulse Rate	HR/PR	60-100 1/min	
	Stroke Volume	Preload	Stroke Volume Index (Output per heart beat)	SVI	40-60 ml/m ²
			Global Enddiastolic Volume Index (Volume of blood in the heart)	GEDI	680-800 ml/m ²
			Intrathoracic Blood Volume Index (Volume of blood in heart & lungs)	ITBI	850-1000 ml/m ²
			Stroke Volume Variation (Dynamic fluid responsiveness)	SVV***	<10 %
			Pulse Pressure Variation (Dynamic fluid responsiveness)	PPV***	<10 %
	Afterload	Systemic Vascular Resistance Index (Resistance of vascular system)	SVRI	1700-2400 dyn*s*cm ⁵ *m ²	
			MAP	70-105 mmHg	
	Contractility	Global Ejection Fraction (Ratio of stroke volume & preload)	GEF	25-35%	
Left Ventricular Contractility (Increase of arterial pressure over time)			dPmx	Trend info - mmHg/s	
Cardiac Function Index (Ratio of CI and preload)			CFI	4.5-6.5 1/min	
Cardiac Power Index (Global cardiac performance)			CPI	0.5-0.7 W/m ²	
Lung	Extravascular Lung Water Index (Lung edema)		ELWI	3.0-7.0 ml/kg	
	Pulmonary Vascular Permeability Index (Permeability of lung tissue)		PVPI	1.0-3.0	

Absolute values (non-indexed values) are only usable in trend screens and have no normal range. *A high-normal / high ScvO₂ can be a sign of insufficient O₂ utilization **14-18 g/dl (Male); 12-16 g/dl (Female) ***SVV and PPV are only applicable in fully ventilated patients with a tidal volume ≥ 8 ml/kg PBW (predicted body weight) and without cardiac arrhythmias

Hemodynamic Decision Model



NOTE: PULSION Medical Systems is a medical device manufacturer and does not practice medicine. PULSION does not recommend these values for use on a specific patient. This decision model is not obligatory.



V+ = volume loading
V- = volume withdrawal
Cat = catecholamine / cardiovascular agents

Please reevaluate your clinical decisions and the set target parameters.

Targeted Values

• GEDI (ml/m ²) (if ELWI >10 → 700-800)	> 700
• GEF (%)	> 25
• CFI (1/min)	> 5
• ELWI (ml/kg) (slow response)	≤ 10

• Volume Responsiveness?
(Passive Leg Raising / Endexpiratory Occlusion Test / Volume Challenge / SVV / PPV?)

• Contractility Problem?
(GEF / CFI / Echo?)

This information is intended for an international audience outside the US. This information is aimed exclusively at healthcare professionals or other professional audiences and are for informational purposes only, is not exhaustive and therefore should not be relied upon as a replacement of the Instructions for Use, service manual or medical advice. Getinge shall bear no responsibility or liability for any action or omission of any party based upon this material, and reliance is solely at the user's risk. Any therapy, solution or product mentioned might not be available or allowed in your country. Information may not be copied or used, in whole or in part, without written permission by Getinge.

Pulsion Medical Systems SE · Hans-Riedl-Str. 17 · 85622 Feldkirchen · Germany · +49 89 45 99 14-0 · zentrale.pulsion@getinge.com · DMS-0006646 · 03/2024

www.getinge.com