

XEN

Chemistries
for the ultimate
peace of mind

Chamber Shine



Cart Washer



Maintenance



Manual



Sterilizers



Washer/
Disinfectant

Instructions for Use

XEN Chamber Shine is a foam-based cleaner that removes scale and mineral deposits from the internal stainless steel surface of washer and sterilizer chambers. Chamber Shine is provided in a foaming manual-trigger spray and is for manually cleaning the inside of chambers. It is not for use when instruments are present and decontamination is taking place.



Order no.	Description
6036005001	Chamber shine 1 x 16 fl. oz.



Directions for Use in Automated Washing Equipment

1. Ensure correct PPE is worn and all skin is covered.
2. Turn the trigger spray to the ON position.
3. Actuate the spray and cover the chamber walls with foam, working from the back of the chamber to the front.
4. Allow the product to soak for up to 30 minutes. If necessary, an acid resistant brush can be used to improve the cleaning. Once the soak time has elapsed, rinse the chamber surfaces thoroughly with water or by wiping with a wet cloth, this time working from the front of the chamber to the back. This cleaning can be repeated if necessary.
5. When finished, run a wash cycle without an instrument load to ensure complete rinsing.



Directions for Use in Sterilizers/Autoclave Chambers

1. Avoid the product entering the drain area.
2. Ensure correct PPE is worn and all skin is covered.
3. Turn the trigger spray to the ON position.
4. Actuate the spray and cover the chamber walls with foam, working from the back of the chamber to the front.
5. Allow the product to soak for up to 30 minutes. If necessary, an acid resistant brush can be used to improve the cleaning.
6. Once the soak time has elapsed, rinse the chamber surfaces thoroughly by wiping with a wet cloth, this time working from the front of the chamber to the back. This cleaning can be repeated if necessary.

Caution

- Do not apply product to hot surfaces or acid sensitive materials.
- Do not mix with other chemicals.
- Wear correct PPE and follow applicable guidelines. Additional information on SDS.
- Check material compatibility before use.

Dispose of contents and container in accordance with local, regional, national, and international regulations.

Contains phosphoric acid, nitric acid, and xanthan gum. For more ingredient information, visit www.getinge.com

Storage Information

Store in a cool, dry place protected from frost and away from alkalis and strong oxidizing agents. Store upright in original containers. Optimum storage temperature 41°F–77°F/5°C–25°C.

	Aluminum	Anodized Aluminum	Titanium	Copper	Brass	Gold	Nickel	Cobalt Chrome	Matte chrome plated	Stainless steel	Glass	Ceramics	Rubber	Plastic	Teflon™/PTFE	Polycarbonate
Chamber Shine										•						

- = The product is suitable for the material
- = Pay special attention to the guidelines from the instruments/device manufacturer

Teflon™ is a trademark of Chemours.

Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries. Copyright 2024 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Manufacturer: Quadralene Ltd, Unit 3, St. Modwen Park, Andresseway Way, Derby, DE21 6YH, United Kingdom.

Sales Office, US: 1 Geoffrey Way, Wayne, NJ 07470, USA

Sales Office, Canada: 90 Matheson Blvd West, Suite 300, Mississauga, Ontario L5R 3R3, Canada

www.getinge.com



MCV00110912 REVA