

# XEN

Chemistries  
for the ultimate  
peace of mind

## Chamber Shine



Cart Washer



Maintenance



Manual



Sterilizers



Washer/  
Disinfectant

## Product Specification

### Description

A foam-based, manual-trigger spray cleaner. Designed to remove scale, mineral deposits, and rust from the internal stainless steel surfaces of washer and sterilizer chambers. Can also be used on stainless steel loading equipment or any other stainless steel surface. It is not for use on instruments or when decontamination is taking place.

### Application

For cleaning internal stainless steel surfaces sterile processing equipment. Refer to equipment and instrument/device manufacturers' instructions before use.

### Key Features

- Recommended for cleaning all stainless-steel washer and sterilizer chambers
- Ready to Use—Quick and easy to use. No diluting or mixing necessary
- Restores stainless steel finish—Removes corrosion and contaminants, prolonging equipment lifespan
- Removes mineral deposits, limescale, and rust—Improves appearance and performance, potentially reducing equipment downtime
- Foaming manual-trigger spray—Remains on walls longer than liquid sprays for better penetration of scale and deposits
- Environmentally friendly—Biodegradable formulation with recyclable packaging



### Ingredients

Water, phosphoric acid, nitric acid, rheology modifier, surfactants, and corrosion inhibitor.

### Technical Data

Specific Gravity / Density (68°F/20°C)	1.024 (typical)
pH as supplied (typical)	1-2
Expiration/Use By Date	24 months from date of manufacture in original sealed containers
Product Dimensions	16 fl. oz. = 4.3"L x 1.8"W x 8.5"H (110 x 46 x 216 mm)

## Material Compatibility

Always follow the cleaning instructions from the manufacturer of the items and instruments and check material for compatibility before use. The following chart serves only as a guide.

	Aluminum	Anodized Aluminum	Titanium	Copper	Brass	Gold	Nickel	Cobalt Chrome	Matte chrome-plated	Stainless steel	Glass	Ceramics	Rubber	Plastic	Teflon™/ PTFE	Polycarbonate
Chamber Shine										•						

- = The product is suitable for the material
- = Pay special attention to the guidelines from the instrument/device manufacturer

Teflon™ is a trademark of Chemours.

## Storage Information

Store in a cool, dry place protected from frost and away from alkalis and strong oxidizing agents. Store upright in original containers. Optimum storage temperature 41°F–77°F/5°C–25°C

## Cautions

- Do not apply to hot surfaces or acid-sensitive materials
- Rinse thoroughly
- Intended for use on stainless steel only
- Autoclaves: Avoid product entering drain
- Always wear correct PPE
- Refer to Safety Data Sheet (SDS) for further safety and handling

## Environmental Awareness

Biodegradable formulation and recyclable packaging.

Dispose of contents and container in accordance with local, regional, national, and international regulations.

## Ordering Information

Size	Order Number
6 x 16 fl. oz.	6036005001

## Service

For technical support in regards to appropriate detergent use, please contact your local Getinge representative.

Or click or scan here:



Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries. Copyright 2023 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Manufacturer: Quadralene Ltd, Unit 3, St. Modwen Park, Andresseway Way, Derby, DE21 6YH, United Kingdom  
Sales Office, US: 1 Geoffrey Way, Wayne, NJ 07470, USA  
Sales Office, Canada: 90 Matheson Blvd West, Suite 300, Mississauga, Ontario L5R 3R3, Canada

[www.getinge.com](http://www.getinge.com)

**GETINGE**

MCV00110901 REVA