

XEN

Chemistries
for the ultimate
peace of mind

Enzymatic+



Orthopedics



Robotic
Instruments



Flexible
Scopes



Rigid
Scopes



Cart Washer



Manual



Ultrasonic



Washer/
Disinfecter

Product Specification

Description

With its low-foam, pH-neutral, ultra-concentrated formulation, XEN Enzymatic+ delivers superb cleaning for a wide range of medical instruments using manual methods, ultrasonics, automated washer-disinfectors, or cart washers. This concentrated, transparent solution has improved foam control, proven biofilm removal properties, and leaves no visible residue.

Application

For cleaning of surgical-grade instruments/reusable items, endoscopes, and other devices. Refer to equipment and instrument/device manufacturers' instructions and AAMI ST79 and AORN cleaning recommendations before use.

Key Features

- Multiuse—Single formulation works for manual, ultrasonic, and automated cleaning
- pH-neutral—Recommended by AAMI and AORN to prevent pitting and discoloration of medical instruments
- Low foaming—Improved foam control optimizes performance of and prevents damage to ultrasonics and automated washers.
- Quadruple-enzyme formulation—Protease, lipase, amylase, and cellulase enzymes attack a range of soils
- Safe for staff—Solution is transparent across temperature range for manual cleaning, increasing sharps' visibility
- Controls water hardness—Neutralizes minerals, leaving no visible residue inside chamber during testing
- Removes biofilm—Testing showed ability to break down biofilm, which can hide bacteria from sterilization

Ingredients

Protease, lipase, amylase & cellulase enzymes, surfactants, sequestering agents and corrosion inhibitors.



Technical Data

Specific Gravity/Density (68°F/20°C)	1.038 (typical)
pH as supplied (typical)	8.6
pH at 0.26 oz./gal. (2 ml/L) water (typical)	8.4
Viscosity (concentrate, 68°F/20°C)	5 cSt (typical)
Expiration/Use By Date	24 months from date of manufacture in original sealed containers
Dosing (Manual and Sonic)	0.26–0.51 oz./gal. (2.0–4.0 ml/L)
Temperature Range (Sonic)	77°F–113°F (25°C–45°C)
Dosing (Automated)	0.03–0.26 oz./gal. (0.2–2.0 ml/L)
Temperature Range (Automated)	112°F–140°F (50°C–60°C)
Product Dimensions	5 L = 7.6" L x 5.7" W x 10" H (193 x 146 x 255 mm)
	5 gal. = 11.6" L x 9.4" W x 15.4" H (295 x 240 x 390 mm)
	15 gal. = 15.7" Dia. x 24" H (399 x 609 mm)

GETINGE

Material Compatibility

Always follow the cleaning instructions from the manufacturer of the items and instruments and check material for compatibility before use. The following chart serves only as a guide.

	Aluminum	Anodized Aluminum	Titanium	Copper	Brass	Gold	Nickel	Cobalt Chrome	Matte-chrome plated	Stainless steel	Glass	Ceramics	Rubber	Plastic	Teflon™/ PTFE	Polycarbonate
Enzymatic+	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•

• = The product is suitable for the material

◦ = Pay special attention to the guidelines from the instrument/device manufacturer

Teflon™ is a trademark of Chemours.

Storage Information

Store in a cool, dry place protected from frost and away from acids and strong oxidizing agents. Store upright in original containers. Optimum storage temperature 41°F–77°F/5°C–25°C

Cautions

- Check material for compatibility before use
- Do not mix with other chemicals
- Always wear correct PPE
- Refer to Safety Data Sheet (SDS) for further safety and handling

Environmental Awareness

Biodegradable formulation and recyclable packaging.

Dispose of contents and container in accordance with local, regional, national, and international regulations.

Ordering Information

Size	Order Number
1 x 5 L	6036002901
1 x 5 gal.	6036003001
1 x 15 gal.	6036003101

Service

For technical support in regards to appropriate detergent use, please contact your local Getinge representative.

Or click or scan here:



Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries. Copyright 2023 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Manufacturer: Quadralene Ltd, Unit 3, St. Modwen Park, Andresseway Way, Derby, DE21 6YH, United Kingdom

Sales Office, US: 1 Geoffrey Way, Wayne, NJ 07470, USA

Sales Office, Canada: 90 Matheson Blvd West, Suite 300, Mississauga, Ontario L5R 3R3, Canada

www.getinge.com

GETINGE 

MCV00110894 REVA