# Prismalix™

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Dear customer,

*Thank you for choosing the **MAQUET SA PRISMALIX surgical light**.*

***PRISMALIX** is an entirely new surgical lighting concept implementing innovative sophisticated optical and mechanical systems.*

***MAQUET SA** has designed the **prismalix** range to ensure the most comfortable working conditions for the surgical team and optimum safety conditions for the patient.*

***PRISMALIX** is the result of **MAQUET SA's** on-going contact with our customers and the preview experience of our engineering and design team. As a result, **PRISMALIX** ensures:*

- ***Superb shadow dilution.***
- ***Easy movement handling.***
- ***Original treatment of the light volume.***

*To get the most from the exceptional features of the **PRISMALIX surgical light**, please carefully read this manual.*

*Don't hesitate to contact **MAQUET SA** with any questions you may have concerning **PRISMALIX** or our other products.*

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## Quality compliance

### Certification of MAQUET SA quality system

LNE/G-MED certifies that the quality system created by Maquet SA for the design, manufacturing, marketing, installation and customer servicing of its surgical lights meets the requirements of the following international standards:

- ISO 9001:2000
- ISO 13485:2004

**PRISMALIX** is designed to fulfill the following applicable standards:

- **EN ISO 14971:2000**  
*Medical devices - Application of risk management to medical devices (ISO 14971:2000)*
- **EN ISO 14971:2000/A1:2003**
- **EN 60601-1:1990**  
*Medical electrical equipment - Part 1: General requirements for safety Amendment A1:1993 to EN 60601-1:1990 Amendment A2:1995 to EN 60601-1:1990 Amendment A13:1996 to EN 60601-1:1990*
- **EN 60601-1-2:2001**  
*Medical electrical equipment - Part 1-2: General requirements for safety - Collateral standard: Electromagnetic compatibility - Requirements and tests*
- **EN 60601-1-4:1996**  
*Medical electrical equipment - Part 1-4: General requirements for safety - Collateral standard: Programmable electrical medical systems Amendment A1:1999 to EN 60601-1-4:1996*
- **EN 60601-1-6:2004**  
*Medical electrical equipment - Part 1-6: General requirements for safety - Collateral standard: Usability*
- **EN 60601-2-41:2000**  
*Medical electrical equipment - Part 2-41: Particular requirements for the safety of surgical luminaires and luminaires for diagnosis*

This product was investigated according to the following additional standards: CAN/CSA-C22.2 No. 601.1-M90 (R2005) (includes National Differences for Canada), EN 60601-1:1990 + A1:1993 + A2:1995 + A13:1996, UL 60601-1, 1st Edition, 2006-04-26 (includes National Differences for USA).

### CE Marking

Compliance with the requirements of Directive 93/42/EEC dated of June 14th 1993, relating to medical devices has been assessed in accordance with Annex VII of this Directive. This **PRISMALIX** Surgical lights range is a class I device in accordance with Annex IX of Directive 93/42/EEC.

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In order to guarantee the quality of our products, it is necessary to contact MAQUET SA in case of use or environment change.

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## CONTENTS

QUALITY COMPLIANCE .....	4
SYMBOLS USED .....	6
WARNINGS.....	7
INSTRUCTIONS FOR SAFE USE .....	8
THE PRISMALIX CONCEPT.....	11
DESCRIPTION AND CHARACTERISTICS.....	15
I) UTILIZATION	
A) STARTUP.....	19
B) FUNCTION AND USE OF STERILIZABLE HANDLE .....	21
C) PREPOSITIONING THE SURGICAL LIGHT.....	23
D) MANEUVERING THE SURGICAL LIGHT .....	25
1) S and S.A. suspensions.....	25
2) A.C.S. suspensions .....	25
3) Handling the A.C.S. suspensions.....	27
II) CLEANING AND MAINTENANCE	
A) CLEANING AND DISINFECTION .....	28
B) MAINTENANCE.....	31
1) Preventive maintenance Four keys to proper maintenance .....	31
2) First level maintenance .....	31
GENERAL CHARACTERISTICS .....	32
SPARE PARTS, ACCESSORIES.....	33
EMC DECLARATION .....	35
TROUBLESHOOTING .....	46
ANNEXES INSPECTION SHEET.....	53

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\* rms voltage: see installation manual of ENERGIX units

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## WARNINGS

 Light is an energy which can potentially dry tissue. The user must adapt the lighting level to the needs of the operation to be performed, in particular when several projectors are used in combination.

 Because of the light intensity of the surgical light in normal mode, it is dis-advised trying to look it facely without protection.

 In the event of a mains power failure, only the cupolas with power supplies connected to a backup power supply system will remain operational.

 When changing a failed bulb :

- Switch off the power supply and allow the projector to cool for 5 minutes.
- Only recommended bulbs should be used.
- Handle the bulb with care using a clean, dry cloth.
- Never touch light bulb with bare hands.

 After each sterilization and before each new use of the sterilizable handle :

- Check that there are no cracks and that the lock operates correctly.
- Fit the handle on the light and lock it by rotating (click).

 Only the surgical team is authorized to handle the sterilizable handle to ensure sanitation.

 Do not use the unit in the presence of inflammable anesthetic gases.

 The unit must be operated with the cover closed.

 When performing maintenance, take necessary precautions with respect to hot surfaces which are indicated by appropriate markings.

 Do not use the wall mounted power supply units as a storage surface.

 Disassembly of certain sub-assemblies can affect operation and safety, for example:

- when servicing the electrical power supply,
- when servicing the suspension arms and balancing system,
- when servicing the optical system of the cupolas which are equipped with filters intended to eliminate radiations not visible by the patient. The unit should not be operated, under any circumstances, without these filters.

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## INSTRUCTIONS FOR SAFE USE

### Dear User

Please note that certain operations may only be carried out by staff with the appropriate qualifications:

- **The ceiling unit may only be operated by trained medical staff.**
- *This unit uses state of the art technology and is completely safe. Nevertheless dangers can arise, especially when it is operated by unqualified staff, when it is not well operated or not used as intended.*
- **Cleaning of the appliance must be carried out by qualified sanitary specialists.**
- *Please read these operating instructions carefully before using the unit. You will then enjoy all the advantages it can offer and protect yourself and others from harm.*
- *Please follow these instructions to train staff in the operation, care and safe use of the unit.*
- *The unit is only intended for use as mentioned in the operating instructions and is only suitable for such use. Any other use can be dangerous to life and limb and/or cause damage to the product or to other property of the operator.*
- *For safety reasons, no alterations or adaptations may be made on your own initiative without the agreement of MAQUET.*
- *For your safety, contact your supplier should any problems arise that are not fully addressed by these instructions.*
- *The contents of the user manual may be changed at any time without notice.*

### Additional equipment:

- *The models named in this document may be equipped with other manufacturers' end appliances (e.g. monitors). For information on operating them, please refer to each manufacturer's user manuals.*
- *The ceiling service unit is designed for continuous operation.*

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#### Transport and storage:

- *Ambient temperature: -25°C to 70°C*
- *Relative humidity: 10 to 75%*
- *Atmospheric pressure: 500 - 1060 hPa*
- *Do not store outside.*
- *Do not subject to severe vibrations.*

#### Operating conditions:

- *Ambient temperature: 10°C to 40°C*
- *Relative humidity: 30 to 75%*
- *Atmospheric pressure: 700 - 1060 hPa*
- *The appliance is not suitable for use in explosive atmospheres.*
- *PRISMALIX must be installed and commissioned in accordance with the EMC information provided on page 35.*
- *Portable RF communication devices may affect the operation of this equipment.*

#### Maintenance

- *Maquet's warranty is only valid and the safety and integrity of the product are only guaranteed if:*
  - *All inspection, maintenance and repair operations are performed by Maquet engineers or trained and authorised technical support technicians.*
  - *Only original accessories, consumables and spare parts are used.*

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**The user must adapt energy level regarding kind of patient and operation :** Patient's photosensitivity could be developed by physical conditions by tissue type where operation is occurred and by medicine action.

*In case of several cupolas combination, he must take care of protecting patient avoiding a long exposure to an energy concentration which are not suitable with the kind of patient and surgery according to international standard IEC 60601-2-41.*

### 1) EXCEPTIONAL VISUAL COMFORT FOR THE SURGICAL TEAM

*The optical core of the PRISMALIX has been computer-optimized to control light down to the smallest particle and direct it to the surgical field.*

*Energy distribution in the volume is fully controlled through multiplication of the number of prisms and their individual angles. All of the light beams coming from the prisms (2184 for 4000 cupola, 3055 for 6000 cupola and 4355 for 8000 cupola) are superposed to build a deep (80 cm) uniform light volume designed to enhance working comfort.*

### 2) A NEW WAY OF TREATING THE LIGHT VOLUME

*The increasingly sophisticated surgical techniques have brought MAQUET SA, original inventor of prismatic surgical lights, to make new breakthroughs in the "prism logic" with its new patented AVL technology (Fig. 1).*

*By simply rotating the central sterile handle, the AVL system (adjustable volume of light) aligns the diameter of the light column with the surgical field to maintain a high level of consistency of light in the entire volume, from the skin to the deepest cavities.*

### 3) UNEQUALLED SURFACE AND INDEPTH SHADOW DILUTION

*The cupolas have very large illuminating surfaces despite their remarkably compact easy-to-handle design.*

*Each of the light sources formed by the total reflection prisms (Fig. 2) acts independently of the others to compensate for obstruction of any given part of the overall light beam. Near the handle, the nearly-axial beam incidence facilitates illumination of cavities and their walls.*

### 4) ACCURATE COLOR SHADE RESTORATION OF TISSUES OPERATED ON

*The color temperatures specially selected by MAQUET SA for PRISMALIX facilitate extended viewing of red tissue. This feature is made possible through the use of special thermal filters combined with optical treatments reflecting the IR emissions of the bulb. The cupolas also have the highest color rendering quality available.*

### 5) PERFECT CONTROL OF THE LIGHT AND HEAT ENERGY IN A CLOSED CUPOLA

*A combination of thermal bridges, optical filters and internal vents ensure optimum surgical team comfort by reducing the nearby heat and radiation.*

*This closed concept, implementing smooth and rounded shapes, thoroughly observes strict, current cleaning and decontamination requirements.*

*High output of the PRISMALIX optical system contributes to efficient distribution of light in the field without hot spots or loss of energy (TABLE 1, Fig. 1).*

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**6) A SHARED VIEW OF THE SURGICAL FIELD FOR THE ENTIRE SURGICAL TEAM**

*Video recording of a surgical procedure is often reserved for teaching surgical practices. MAQUET SA has made this feature accessible to everyone for other applications, as standard equipment:*

- *Optimization of surgical operation planning, synchronization of surgical team interventions (anesthesiologist and surgeon, instrument person and surgeon) using the fixed focus PRISMAVISION CFF monitoring camera.*
- *Future applications such as videoconference and telemedicine using the variable focus PRISMAVISION ZOOM camera.*

**7) A "FEATHER-WEIGHT" CUPOLA**

Proper positioning the PRISMALIX cupolas is facilitated by:

- *Compact, light-weight design.*
- *Suspension geometry ensuring complete rotation of the satellite cupola around the main cupola.*
- *A large number of swivel joints for optimum handling and directional control.*

**8) A REVOLUTIONARY SUSPENSION**

PRISMALIX lights equipped with Automatic Compensated Suspension (ACS) (Fig. 4) provide:

- *Angular rotation of 210° facilitating easy handling. The cupola can go beyond the vertical position without encountering limit stops.*
- *Two action radii for main cupolas providing greater flexibility in swivel movements.*
- *Immediate, absolute stability after positioning.*

**9) A NEW ELECTRICAL POWER SUPPLY CONCEPT**

*The ENERGIX WPS power supply units are specifically designed by MAQUET SA to supply its complete range of lights installed throughout the world. These power supply units feature:*

- *Stable lighting, insensitive to power supply system fluctuations.*
- *Easy and fast to install by their ergonomic, modular, compact design.*
- *Optimized service life of bulbs through use of " slow start " function and electronic current and voltage regulation features.*
- *Battery charger integrated in power supply unit.*
- *Specific MAQUET SA power unit available on option.*
- *User interface implements touch-sensitive keys and LCD display.*
- *Compatible with most demanding disinfection standards.*
- *Electronic management of consumables: bulbs, bulb sockets and batteries.*
- *Programming of safety checks and maintenance inspections.*
- *Integrated maintenance software.*
- *Remote control by R5232 serial link.*

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PRX 4001 S.A.X, PRX 4401 S.A.X,  
PRX 4441 S.A.X.

**S SERIES**

*The S series surgical light (Standard Suspension) comprises (Fig. 5.1) the following :*

- *Spacer tube (1).*
- *Main arm (2) with balancing arm.*
- *Cupola arm (3).*
- *Type 4000 cupola (4).*

*The S series comprises the following models:*

*PRX 4001 S, 4401 S (ceiling-mounted), PRX 4002 S (wall-mounted), PRX 4003 S (mobile).*

**S.A. SERIES**

*The S.A. series surgical light (Enhanced Suspension) comprises (Fig. 5.2):*

- *A single, double or triple suspension comprising one, two or three main arms (1), slanted or straight.*
- *One or two spring-compensated balancing arms (2).*
- *A combination of type 6000 (3) and type 4000 (4) cupolas.*

*The S.A. series comprises the following models:*

*PRX 4001 S.A, PRX 4401 S.A,  
PRX 4441 S.A, PRX 6001 S.A,  
PRX 6401 S.A, PRX 6101 S.A,  
PRX 6441 S.A, PRX 6411 S.A,  
PRX 6601 S.A, PRX 3001 S.A,  
PRX 4301 S.A, PRX 4331 S.A,  
PRX 4431 S.A, PRX 6301 S.A,  
PRX 6431 S.A, PRX 2001 S.A,  
PRX 4201 S.A, PRX 4421 S.A,  
PRX 6201 S.A, PRX 6421 S.A,  
PRX 4101 S.A, PRX 4411 S.A.*

**A.C.S SERIES**

*The A.C.S. series (Automatic Compensated Suspension) comprises (Fig. 5.3):*

- *A single, double or triple suspension comprising one, two or three main arms (1).*
- *One or two compensation units (4).*
- *A balancing arm (2) when the light carries a 4000 satellite (3).*
- *One or two cupola arms (6) used to secure the 6000 or 8000 cupolas (5)*

*The A.C.S. series comprises the following models:*

*PRX 6001 A.C.S, PRX 6401 A.C.S,  
PRX 6101 A.C.S, PRX 6441 A.C.S,  
PRX 6411 A.C.S, PRX 8401 A.C.S,  
PRX 8441 A.C.S, PRX 8601 A.C.S,  
PRX 8431 A.C.S, PRX 6201 A.C.S,  
PRX 6421 A.C.S, PRX 8201 A.C.S,  
PRX 8421 A.C.S, PRX 6301 A.C.S,  
PRX 6431 A.C.S, PRX 8411 A.C.S,  
PRX 8441 A.C.S, PRX 8301 A.C.S,*

*The S.A.X. series comprises the following models:*

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**CUPOLAS**

*Outside handles (2) facilitate optimum prepositioning before the surgical team arrives. By using multiple handles (3 handles on the 6000 cupola and 4 handles on the 8000 cupola), all areas of the light are fully accessible. The type 4000 cupola, is easily prepositioned, raised or lowered by its handle.*

*External handles (1) used with the 6000 and 8000 cupola arms provide easy up-down movement of the cupola.*

*PRISMALIX cupolas are equipped with new generation PRISMALIX undersides (4) providing more powerful lighting, excellent depth of field and complete shadow dilution.*

*The sterilizable handle (3) is used by the surgeon to direct the light with respect to his view of the surgical field. By twisting handle also ensures optimum volume control of the light with respect to the size of the surgical field (see chapter I, para. B "FUNCTIONS AND USE OF STERILIZABLE HANDLE").*

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*For a large number of surgical operations, the copulas are positioned parallel to the operating table.*

**In special cases,** the surgical team may want to use the copula in a **vertical position.**

*In this case, with the 6000 and 8000 copulas, as a safety measure and to avoid any unwanted concentration of heat on the illuminating underside of the copula, we suggest that you chose one of the following solutions:*

- *If the lighting level can be reduced to 2/3 of the maximum intensity, the copula can keep the vertical position.*

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- *If maximum lighting intensity must be maintained, slightly incline the copula with respect to the vertical, keeping the graduator at the maximum output.*

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
**Reminder:** The 6000 and 8000 copulas are equipped with a device that automatically switches off the light when the copula goes past the vertical position by an angle of 10° (copula illuminating ceiling).

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
- 1 - Power on
- 2 - Key used to increase lighting level
- 3 - Key used to decrease lighting level
- 4 - Lighting level indicator
- 5 - Failed bulb indicator / transition to Ambient light
- 6 - Green indicator light indicates normal operation. The indicator light goes red when operating from batteries (mains power failure).

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## I - USE

-  - Do not use the unit in the presence of inflammable anesthetic gases.
- To benefit from the manufacturer's warranty, only recommended power supplies must be used.
- Because of the light intensity of the surgical light in normal mode or in ambient mode, it is disadvised trying to look it facely without protection.

### A) Start Up

#### Startup and lighting intensity control

The light is switched on and controlled using the ENERGIX WPS power supply units (see corresponding manuals).

When button (1) is pressed, the green indicator light (6) comes on to indicate that the power supply is in normal operation.

The concerned cupola comes on gradually to stabilize at between 30 and 100% of its maximum lighting capacity in accordance with the last value stored (lighting value on last shutdown).

The lighting level indicator (4) has 4 LED, each of which comes on at half then completely, to provide 8 possible lighting levels for the cupola.

The lighting level is adjusted using the (+) key (2) and (-) (3) key.

The lighting is adjusted by successively pressing, or holding the button pressed, to obtain the desired lighting level.

#### Operation of backup battery option (6) (See corresponding documents)

As maybe required, the power supplies can be equipped to receive the 24 V backup system for the operating theater.

In the event of a mains power failure and transition to the backup power supply, indicator light (6) goes red and the cupola operates at its maximum level.

In the event of a mains power failure, only the cupolas with power supplies connected to a backup power supply system will remain operational.

Following a cutout and return of the mains power, the light will start up directly in surgical light mode.

#### Operation of lamp failure indicator

The yellow indicator (5) flashes when the power supply is on and the bulb has failed and when operating with ambient light.

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### 1) INSTALLING AND REMOVING HANDLE

#### INSTALLING HANDLE (Fig. 9.1)

Engage the sterilizable handle on its support to "click" it in to its locked position.

Turn the handle on its shaft until it clicks again.

#### REMOVING THE HANDLE (Fig. 9.2)

To remove the sterilizable handle, push and hold the button while pulling the handle downward to remove it from its support.

### 2) FUNCTION

#### CUPOLA WITH CAMERA (Fig. 9.3)

Cupolas equipped with camera do not come with a light volume control system.

By rotating the handle, the image on the monitor screen is rotated so that it is always properly aligned with respect to the actual position of the operating team and/or observers.

The camera starts automatically when PRISMALIX is switched on. The image is only generated provided the PRISMALIX has been connected to a monitor by means of a cable.

The camera is easy to use as it is directly aligned with the surgical light when directed to the surgical field. The camera lens is designed to cover a 25 cm diameter field at a distance of one meter. The size of the image depends on the distance between the camera and the field observed: the closer the surgical field, the larger the field image. The image is sharp for all distances between the handle and the surgical field between 70 cm and 170 cm (1.7 m).

#### CUPOLA WITHOUT CAMERA (Fig. 9.4) Page 20

The surgeon can vary the light volume by simply rotating the sterilizable handle. In this way, the surgeon can decrease or increase the diameter of the light volume. A limit stop indicates the end of travel for the volume adjustment.

For example:  
**INITIAL SURGICAL PROCEDURE, INCISION, SUPERFICIAL WOUND**  
-> ENLARGED LIGHT SPOT  
**IN-DEPTH WORK**  
-> CONCENTRATED LIGHT SPOT

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### C) PREPOSITIONING THE SURGICAL LIGHT

*The surgical light must be prepositioned prior to any procedure to minimize subsequent handling. By prepositioning the surgical light appropriately for each procedure, potential interactions with possible obstacles are limited (Serum holder, distribution arm, etc.).*

*The cupolas can be handled using the peripheral handles (4). Macro handling is facilitated by the lights ability to rotate 360° around the main axis (Fig. 10.1).*

**For example**, for a thoracic operation, the main cupola (2) will be positioned above the back rest plate of the table (1) with the satellite cupola (3) positioned at an angle such that the light beam provides the surgical team with additional light, avoiding heads and bodies.

*The illustration showing the position of the two cupolas with respect to each other is given as an initial prepositioning **example only**, being understood that each surgical procedure can require different positioning.*

*The prepositioning operation should be performed, preferably, before the patient is brought into the operating theater. The handles on the cupola and the cupola arms (A.C.S Series) can be used for this purpose (Fig. 10.2).*

 **Once the sterile handle is installed at the center of the cupola, it should only be touched by the surgical team members.**

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**D) HANDLING THE SURGICAL  
LIGHTS**

**1) S. AND S.A. SUSPENSIONS**

*The movements enabled by the different suspensions are described in Figures 11.1 and 11.2.*

*The cupola arm has an inclination of 45° upward and 50° downward.*

*The angle can be adjusted to avoid hitting the ceiling.*

**2) A.C.S. SUSPENSIONS**

*The movements enabled by these suspensions are described in Figure 11.3.*

*The 6000 and 8000 cupola arms have an inclination of 15° upward and 180° downward.*

*Note that the cupola position, with the balancing arm entirely folded against itself, can be offset with respect to the anchoring axis (See pages 26/27).*

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


*All PRISMALIX models can be positioned with a high degree of precision. The ACS system provides additional advantages in this respect. In addition to their high flexibility the uniqueness of the ACS series lights is due to the exceptional swivel capability provided by the **compensation unit**.*

*With this system the light can go from an extended position A to a folded position B with no limit stop (Fig. 12.1).*

*The movement is performed using handles on the cupolas and the cupola arms.*

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*To take into account the position of the operating table and the patient with respect to the lighting axis, it is sometimes useful to **offset** the cupola axis in reference to its anchoring axis.*

*Figures 12.2 and 12.3 illustrate these positions, clearly showing how the ACS unit allows the cupola to turn along action radius 1 when in the normal position, and along action radius 2 in the deployed position.*

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## II - CLEANING AND MAINTENANCE

### A) CLEANING/DISINFECTION STERILIZATION

The user must contact the sanitary specialists at the medical center. The recommended products and procedures must be observed.

In the event of doubt concerning the compatibility of active agents to be used, contact the local MAQUET SA customer service.

**2.1 - Cleaning and disinfection of the surgical lighting system :**  
Before cleaning, ensure the power is turned off and the light is cooled down.

**General instructions concerning cleaning, disinfection and safety :**

- Remove the sterilisable handles.
- Wipe the system with a cloth moistened with a surface cleaner. Follow the manufacturer's dilution and temperature recommendations.
- Rinse with a cloth moistened with water. Wipe dry.
- Wipe evenly with a cloth moistened with disinfectant. Follow the manufacturer's recommendations.
- Rinse with a cloth moistened with water to remove any residue (particularly when cleaners containing aldehydes, quaternary ammonium or surfactants are used).
- Wipe off with a dry cloth.

#### a) Examples of recommended products

- **GETINGE USA Products:** TEC QUAT 256
- **ANIOS Products:** SURFA'SAFE; HEXANIOS G + Rat 0,5%; ANIOSYME P.L.A ; SALVANIOS pH10 ; ANIOS DDSH
- **Schülke & Mayr Products:** Antifect Plus

#### b) Examples of prohibited products

 Solutions containing glutaraldehyde, phenol, iodine, bleach, alcohol or chloride ions must not be used.

 Do not use disinfection by fumigation methods.

### 2.2 - Sterilizable handles: cleaning and sterilization

#### 2.2.1 - Preparation for Cleaning

Use soft cloth to wipe away gross soils from light handle surface, immediately after its use.

Place handles in containment device for transporting and keep soils moist to ease further cleaning.

Avoiding situations that soil the inside surfaces of handles. This can help ease the required cleaning efforts.

Do not scratch the lens (video handles only).

#### 2.2.2 - Cleaning

Submerge handles while soaking and washing with a detergent<sup>1</sup> solution.

Let handles soak for 15 minutes and then manually wash soils off, using soft cloth and soft brushes.

#### 2.2.3 - Inspect Cleanliness of Surfaces

Inspect handle periodically to assure all soils are removed from outer and also inner surfaces.

If soils are seen on internal parts, then repeat cleaning or use an ultrasonic cleaning process, as necessary.

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**Rinsing:**

Thoroughly rinse handle using a brisk stream of clean tap water to remove cleaning solutions. Inspect handle for material properties (no cracks, chips, or crazing), cleanliness and assure it functions properly.

**2.2.4 - Drying**

Wipe dry with freshly cleaned soft cloth. Prior to reuse, the light handles are required to be steam sterilized.

**2.2.5 - Disinfection**

Handles could be disinfected with a machine (Getinge washers for example) and rinsed with a temperature maxi. of 93°C.

Stage	Exp. Temp.	Exp. Time
Prewash	18 - 35° C	60 sec
Wash	46 - 50° C	300 sec
Neutralization	41 - 43°C	30 sec
Wash 2	24 - 28°C	30 sec
Rinsing/Disinfection	92 - 93°C	600 sec
Drying		20 min

**2.2.6 - Sterilization**

Cleaned light handles are to be steam sterilized by using the following instructions and sterilization cycle:

Country	Sterilization cycle	Exp. Temp [°C]	Exp. Time [min]	Dry [min]
USA & Canada	Prevacuum <sup>2</sup>	132 - 135	10	16
France	ATNC (Prion) (Prevacuum)	134	18	
Autres pays	Prevacuum	Follow local / national regulations requirements		

- Re-inspect each handle to assure it is clean and free of defects before proceeding.
- Wrap handle with sterilization wrapper material (double wrapper or equivalent). Alternatively, handles may be placed in paper/plastic sterilization pouches<sup>3</sup> to ease identification after reprocessing and prior to reuse.
- Place on perforated tray with handle opening facing downwards<sup>4</sup>.
- Package with biological and/or chemical indicators for monitoring the sterilization process, per local regulations and facility policies.
- Run sterilization cycle, abiding to sterilizer manufacturer's instructions for load configurations.

**To ensure correct sterilization, do not allow any soiling to penetrate inside the handle.**

**With the sterilization parameters indicated, the sterilizable handles are only guaranteed for no more than 50 utilisation cycles.**

**Disposal of the sterilizable handles is subject to the standard rules governing disposal of hazardous products in a hospital environment.**

<sup>1</sup> A non-enzyme-based detergent is recommended (e.g., Getinge USA NeutraWash). Many enzyme based detergents can adversely affect the plastic material used in this handle and should not be used for extended periods of soaking and will need to be thoroughly rinsed off, if used.  
<sup>2</sup> This light handle is a porous instrument.  
<sup>3</sup> Possible suppliers of sterilization pouches:  
Medical Action Industries  
SBW Medical  
Baxter International  
<sup>4</sup> For air removal and drying purposes

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## B) MAINTENANCE

### 1) PREVENTIVE MAINTENANCE:

SEE MAINTENANCE KEYS IN CUSTOMER TECHNICAL  
MANUEL (CTM).

👉 **Maintenance and inspection (to be performed by a certified technician).**

To preserve your surgical lights' original performance and reliability, annual maintenance and inspections should be performed as follows:

- during the warranty period, by an MAQUET SA technician or an MAQUET SA-certified distributor,
- outside the warranty period, by an MAQUET SA technician or MAQUET SA-certified distributor, or the hospital's engineering department.

MAQUET SA will not accept liability resulting from handling operations which do not comply with the present instructions.

### 2) FIRST LEVEL MAINTENANCE

#### a) Daily inspection

- Check that bulbs operate correctly (Fig. 13.2 and 13.3).
- Check that sterilizable handle clicks and locks in correctly (see chapter I.B).

#### b) Monthly inspection

- Check operation of light on backup power supply system if applicable used in the event of primary mains supply cutout.

#### c) Replacement of factory parts

**PRELIMINARY REMARK:** The bulb holders must be replaced on a yearly basis. MAQUET SA guarantees the optical performance of its PRISMALIX surgical lights provided they are equipped with bulbs ECL0001 solely, to the exclusion of any other type of bulb.

**Bulb replacement procedure:**  
(Fig. 13.1)

👉 **Do not touch electrical connections and the patient at the same time.**

- 1- Identify the faulty bulb. Start by checking the type of unit involved by reading the nameplate on the cupola.
- 2- Switch off the projector's electrical power supply.
- 3- Turn the quarter-turn knob and open the upper cover.
- 4- **CAUTION: Before starting any maintenance operation, allow the projector bulb holder to cool for 5 minutes.**
- 5- Remove the lamp holder taking care not to knock the hot globe of the bulb against hard surface.
- 6- **CAUTION: Never handle the bulb and the bulb holder with your bare hands. Always use a dry, non-greasy cloth to remove a faulty bulb. Retract the spring clamp to remove the bulb.**
- 7- Remove the bulb from its cover. Take it by fins taking care to not touch the onion-shaped dome. Place the bulb on the leafs checking the right positioning with the foolproof plug and then clamp the bulb using the spring clamps.
- 8- Reinstall the bulb holder, checking its right insert.
- 9- Close the cover and lock the assembly by the quarter-turn knob.

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**NOTES:**

(1) All measures taken at 1 m below lowest point of illuminating underside with concentrated spot.

(2) Central illumination (Ec) is the minimum illumination guaranteed on purchase of the product.

(3) Periphery of light spot where illumination reaches 10% of central illumination.

(4) Circle where illumination reaches 50% of central illumination.

(5) All of these values are defined with respect to the central illumination (Ec) measurements.

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





**Table 201 - Guidance and Manufacturer's Declaration - Electromagnetic Emissions**

PRISMALIX surgical lights are intended for use in the electromagnetic environment specified below. PRISMALIX customers or users should ensure that they are used in such an environment.

<b>Immunity test</b>	<b>Compliance</b>	<b>Electromagnetic environment — directives</b>
RF emissions CISPR 11	Group 1	PRISMALIX lights use RF energy only for their internal functions. Their RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	PRISMALIX lights are suitable for use in all establishments other than domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations and flicker emissions IEC 61000-3-3	Not applicable	

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**Table 202 - Guidance and manufacturer's declaration - electromagnetic immunity**

PRISMALIX surgical lights are intended for use in the electromagnetic environment specified below. PRISMALIX customers or users should ensure that they are used in such an environment..

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment — guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	< 5% $U_T$ (dip > 95% of $U_T$ ) for 0.5 cycles  40% $U_T$ (dip = 60% of $U_T$ ) for 5 cycles  70% $U_T$ (dip = 30% of $U_T$ ) for 25 cycles  < 5% $U_T$ (dip = 95% of $U_T$ ) for 5 cycles	< 5% $U_T$ (dip > 95% of $U_T$ ) for 0.5 cycles  40% $U_T$ (dip = 60% of $U_T$ ) for 5 cycles  70% $U_T$ (dip = 30% of $U_T$ ) for 25 cycles  < 5% $U_T$ (dip = 95% of $U_T$ ) for 5 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the PRISMALIX requires continuous operation during mains interruptions, it is recommended that it be powered from an uninterruptible power supply.
Power frequency (50/60 Hertz) magnetic fields IEC 61000-4-8	3 A	3 A	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note:  $U_T$  is the voltage of the AC supply network before applying the test level.

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**Table 204 - Guidance and manufacturer's declaration - electromagnetic immunity**

PRISMALIX surgical lights are intended for use in the electromagnetic environment specified below. PRISMALIX customers or users should ensure that they are used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment — guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 V rms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 V rms 3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the PRISMALIX, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = \left[ \frac{3.5}{V1} \right] \sqrt{P} = 1.17 \sqrt{P}$ $d = \left[ \frac{3.5}{E1} \right] \sqrt{P} \text{ from 80 MHz to 800 MHz} = 1.17 \sqrt{P}$ $d = \left[ \frac{7}{E1} \right] \sqrt{P} \text{ from 800 MHz to 2.5 GHz} = 2.34 \sqrt{P}$ <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, and <math>d</math> is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic survey on site<sup>a</sup>, should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the PRISMALIX is used exceeds the applicable RF compliance level above, the PRISMALIX should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the PRISMALIX.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m..

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**Table 206 - Recommended separation distance between portable and mobile RF communication devices and PRISMALIX lights.**

PRISMALIX is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the PRISMALIX can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the PRISMALIX as recommended below, according to the maximum output power of the communications equipment.

Maximum rated output power of transmitter <b>W</b>	Separation distance depending on transmitter frequency <b>m</b>		
	150 kHz to 80 MHz $d = [ 1.17 ]\sqrt{P}$	80 MHz to 800 MHz $d = [ 1.17 ]\sqrt{P}$	800 MHz to 2.5 GHz $d = [ 2.34 ]\sqrt{P}$
	0,01	0.12	0.12
0.1	0.37	0.37	0.74
1	1.17	1.17	2.34
10	3.69	3.69	7.38
100	11.70	11.70	23.40

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters ( $m$ ) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts ( $W$ ) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

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Caution: Disassembly of certain sub-assemblies can affect operation and safety, for example:

- when intervening on the electrical power supply,
- when intervening on the suspension arms and balancing system,
- when intervening on the optical system of the cupolas which are equipped with filters intended to eliminate radiations not visible by the patient. The unit should not be operated, under any circumstances, without these filters.

For this type of intervention, contact the MAQUET SA-approved product support service.

## TROUBLESHOOTING

ANOMALY	PROBABLE CAUSE	CORRECTIVE ACTION
<input type="checkbox"/> <b>Cupola does not come on though green power supply indicator light is lit.</b>	1 - Bulb failure or missing bulb.	Switch off power supply and replace bulb(s).
	2 - Other cause	Call engineering department.
<input type="checkbox"/> <b>Light volume variation control does not operate.</b>	1 - Power supply indicator light is red: mains system has failed and unit is operating on hospital 24V backup system.	1- Wait for return to normal power system or startup of hospital power generator.
	2 - Power cut.	2- Check if another equipment is connected.
	3 - Other cause.	3- Call engineering department.
<input type="checkbox"/> <b>All cupolas not lit</b>	- Each cupola has its own individual control.	Check that all the power supplies are properly switched on (green indicator light).
<input type="checkbox"/> <b>Power supply does not come on.</b>	1 - Mains power failure	Check if another piece of equipment is operating on the same electrical power system.
	2 - Other cause	Call engineering department.
<input type="checkbox"/> <b>6000 or 8000 cupola goes off when handled.</b>	1 - Protection device trips when cupola exceeds vertical position.	Return cupola to normal utilization position (see cupola description, page 15).
	2 - Other cause	Call engineering department.
<input type="checkbox"/> <b>Unit goes off completely during an operation</b>	- Power supply failure or abnormal cupola position.	Check for presence of the main power system on the control unit: - green indicator light should be lit. If not lit, call engineering department. - if lit, check position of cupola which should not go past the vertical position.
<input type="checkbox"/> <b>Low bulb service life</b>	- Bulbs do not comply or overvoltage condition is present.	1) Ensure you are using bulbs recommended by MAQUET SA. Do not use bulbs used with surgical lights of previous generation (PRISMATIC) which are not compatible for power and size. 2) Call engineering department to verify the voltage of batteries (24 Volts) of emergency supply . 3) Call engineering department to verify AC voltage at bulb terminals.  $23 V_{EFF} \left( V_{AC+DC}^{RMS} \right)$
<input type="checkbox"/> <b>Sterilizable handle does not engage correctly on its support.</b>	1 - Average service life of sterilizable handle: 100 sterilization cycles at 134° ± 4°C for 20 minutes.	Check operation of handle lock mechanism (click) and overall handle.
	2 - Average service life exceeded / Handle is deformed.	Call engineering department.

## TROUBLESHOOTING

ANOMALY	PROBABLE CAUSE	CORRECTIVE ACTION
❑ <i>Cupola drift</i>		<i>Call engineering department.</i>
❑ <i>Cupola too loose or too difficult to handle.</i>		<i>Call engineering department.</i>
❑ <i>Light volume variation control does not operate.</i>	<i>1 - Only certain systems are equipped with light volume variation system.</i>	<i>Check if cupola is equipped or not with a camera.  If cupola is equipped with a camera, cupola does not have a light volume variation function (see chapter I, para B2).  If cupola is not equipped with a camera, call engineering department.</i>
	<i>2 - Other cause</i>	<i>Call engineering department.</i>
❑ <i>Difference between small and large light spot is difficult to distinguish.</i>	<i>1 - User is not positioned at correct distance with respect to light spot.</i>	<i>Check operation of light volume variation control by taking position at 1 meter from illuminated surface. - When handle is turned clockwise to its limit stop position, light spot is enlarged (Wide). - When handle is turned counter-clockwise to its limit stop position, light spot is reduced (Narrow).</i>
	<i>2 - Other cause</i>	<i>If mechanism does not operate, call engineering department.</i>
❑ <i>Lighting level drops abruptly though power supply has not been interrupted (6000 and 8000 cupola)</i>	<i>- Bulb has failed during operation. - 6000 and 8000 cupolas have several bulbs. Certain bulbs do not operate.</i>	<i>Switch off power supply and replace bulb(s) (at end of surgical procedure).</i>
❑ <i>Voltmeter indicates a maximum value of around 20 V</i>	<i>Voltmeter not suited to shape of voltage measured</i>	<i>Refer to technical manual or contact our product support service</i>

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*Our network is at your disposal if you have any queries.*

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