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

Post surgical complications increase hospital length of stay by +75%²

Post surgical complications increase readmission rates by +45%²

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
Examples from high to low risk of procedure: aortic aneurysm, oesophageal resection, rectal resection, knee arthroscopy

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Type	elective	emergency
Duration	<2h	>4h
Blood loss	minimal	>25%

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ASA	I	II	III	IV	V
Heart					
CHD	∅		moderate	severe	
EF	normal		reduced	<30%	
NYHA		I II		III IV	
Kidney					
GFR	normal		reduced	<30%	
Lung					
FEV1	normal		reduced	<30%	

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



Labels in the diagram:
 - Patient-specific calibration factor (determined with thermodilution) points to 'cal'.
 - Area under the pressure curve points to the integral symbol.
 - Shape of pressure curve points to the fraction's denominator.
 - Heart rate points to 'HR'.
 - Compliance points to 'C(p)'.

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$$CO = \frac{(T_b - T_i) \times V_i \times K}{\int \Delta T_b \times dt}$$

Labels in the diagram:
 - Blood temperature points to T_b .
 - Injectate temperature points to T_i .
 - Injectate volume points to V_i .
 - Correction constant* points to 'K'.
 - Area under the thermodilution curve points to the denominator integral.

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Literature

1. Salzwedel et al., Perioperative goal-directed hemodynamic therapy based on radial arterial pulse pressure variation and continuous cardiac index trending reduces postoperative complications after major abdominal surgery: a multicenter, prospective, randomized study. *Crit Care* 2013; 17(5): R191.
2. Michard et al., Potential return on investment for implementation of perioperative goal-directed fluid therapy in major surgery: a nationwide database study. *Perioper Med* 2015, 4(11): 1-8.
3. Adapted from: Bellamy, Wet, dry or something else? *Br J Anaesth* 2016, 97(6): 755-757.
4. Adapted from: Cannesson, Arterial Pressure Variation and Goal-Directed Fluid Therapy. *J Cardiothorac Vasc Anesth*. 2010, 24(3):487-497.
5. Nygren et al. Guidelines for Perioperative Care in Elective Rectal/Pelvic Surgery: Enhanced Recovery After Surgery (ERAS®) Society Recommendations. *World J Surg* 2013, 37: 285-305.
6. Nelson et al., Guidelines for pre- and intra-operative care in gynecologic/ oncology surgery: Enhanced Recovery After Surgery (ERAS®) Society recommendations — Part I. *Gynecol Oncol*. 2016, 140(2): 313-322.
7. Wesseling KH et al. A simple device for the continuous measurement of cardiac output. *Adv Cardiovasc Phys* 1983; 5: 16-52.
8. Baudendistel LJ et al. Evaluation of extravascular lung water by single thermal indicator. *Crit Care Med* 1986; 14(1):52-56.
9. Michard F. et al., Global end-diastolic volume as an indicator of cardiac preload in patients with septic shock. *Chest* 2003, 124(5): 1900-1908.
10. Adler C. et al., Fluid therapy and acute kidney injury in cardiogenic shock after cardiac arrest. *Resuscitation* 2013, 84(2): 194-199.
11. Brown LM, Calfee CS, Howard JP, Craig TR, Matthay MA, McAuley DF. Comparison of thermodilution measured extravascular lung water with chest radiographic assessment of pulmonary oedema in patients with acute lung injury. *Ann Intensive Care* 2013; 3(1): 25.
12. Saugel B, Ringmaier S, Holzapfel K, Schuster T, Phillip V, Schmid RM, Huber W. Physical examination, central venous pressure, and chest radiography for the prediction of transpulmonary thermodilution-derived hemodynamic parameters in critically ill patients: A prospective trial. *J Crit Care* 2011; 26(4): 402-10.
13. Lemson J, van Die LE, Hemelaar AE, van der Hoeven JG. Extravascular lung water index measurement in critically ill children does not correlate with a chest X-ray score of pulmonary edema. *Crit Care* 2010; 14(3): R105.
14. Sakka S. G. et al., Assessment of cardiac preload and extravascular lung water by single transpulmonary thermodilution. *Intensive Care Med* 2000, 26(2): 180-187.
15. Pearse R. M. et al., Changes in central venous saturation after major surgery, and association with outcome. *Crit Care* 2005, 9(6): R694-699.

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