# SEVOFLURANE CONSUMPTION IN MANUAL VS AUTOMATIC MODES

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# Sevoflurane consumption in manual vs automatic modes

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(1) Department of Anesthesiology, Reanimation and Intensive Care, AZ Sint Jan Brugge-Oostende, Brugge, Belgium. (2) Department of Anesthesia, Ghent University, Ghent, Belgium. (3) Department of Anesthesia and Critical Care Medicine, Maria Middelaers Hospital, Ghent, Belgium. (4) Department of Cardiovascular Sciences, KULeuven, Leuven, Belgium. (5) Department of Anesthesiology, UZLeuven, Leuven, Belgium. (6) Department of Anesthesiology, KULeuven - University of Leuven, Leuven, Belgium

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» all anesthesia machines should be upgraded as soon as possible with automatic delivery technology to minimize atmospheric pollution with volatile anesthetics. «

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» minimal adjustments in daily practice can dramatically reduce the environmental impact of anesthesia «

» AGC improved stability and arguably additionally improves safety, and should therefore be advocated also from a clinical perspective. «

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» On an ecological level, it should be emphasized that the patient receives the same level of anaesthesia, with even increased safety, but with a lower cost for both society and the biosphere. «

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1. Charlesworth M, Swinton F. Anaesthetic gases, climate change, and sustainable practice. *Lancet Planet Health*. 2017;1(6):e216–7. [https://doi.org/10.1016/S2542-5196\(17\)30040-2](https://doi.org/10.1016/S2542-5196(17)30040-2).
2. Masson-Delmotte V, Zhai P, Pirani A, Connors SL, Péan C, Berger S, Caud N, Chen Y, Goldfarb L, Gomis MI, Huang M, Leitzell K, Lonnoy E, Matthews JBR, Maycock TK, Waterfield T, Yelekçi O, Yu R, Zhou B (eds.). IPCC, 2021: summary for policymakers. In: climate change 2021: the physical science basis. contribution of working group I to the sixth assessment report of the intergovernmental panel on climate change. Cambridge University Press. In Press.
3. Ryan SM, Nielsen CJ. Global warming potential of inhaled anesthetics: application to clinical use. *Anesth Analg*. 2010;111:92–8. <https://doi.org/10.1213/ANE.0b013e3181e058d7>.
4. Van Norman GA, Jackson S. The anesthesiologist and global climate change: an ethical obligation to act. *Curr Opin Anaesthesiol*. 2020;33(4):577–83. <https://doi.org/10.1097/ACO.0000000000000887>.
5. Laster MJ, Fang Z, Eger EI 2nd. Specific gravities of desflurane, enflurane, halothane, isoflurane, and sevoflurane. *Anesth Analg*. 1994;78(6):1152–3. <https://doi.org/10.1213/00000539-199406000-00022>.
6. Brattwall M, Warrén-Stomberg M, Hesselvik F, Jakobsson J. Brief review: theory and practice of minimal fresh gas flow anesthesia. *Can J Anaesth*. 2012;59(8):785–97. <https://doi.org/10.1007/s12630-012-9736-2>.
7. Ghijssels I, De Cooman S, Carette R, Peyton P, De Wolf A, Hendrickx J. Performance of an active inspired hypoxic guard. *J Clin Monit Comput*. 2016;30:63–8. <https://doi.org/10.1007/s10877-015-9684-9>.
8. Struys M, Kalmar A, De Baerdemaeker L, Mortier E, Rolly G, Manigel J, et al. Time course of inhaled anaesthetic drug delivery using a new multifunctional closed-circuit anaesthesia ventilator. In vitro comparison with a classical anaesthesia machine. *Br J Anaesth*. 2005;94:306–17. <https://doi.org/10.1093/bja/aei051>.
9. Carette R, De Wolf A, Hendrickx J. Automated gas control with the Maquet FLOW-i. *J Clin Monit Comput*. 2016;30:341–6. <https://doi.org/10.1007/s10877-015-9723-6>.
10. Moran P, Barr D, Holmes C. Saving sevoflurane: automated gas control can reduce consumption of anesthetic vapor by one third in pediatric anesthesia. *Pediatr Anesth*. 2019;29:310–4. <https://doi.org/10.1111/pan.13587>.
11. Suttner S, Boldt J. Low-flow anaesthesia. Does it have potential pharmacoeconomic consequences? *Pharmacoeconomics*. 2000;17(6):585–90. <https://doi.org/10.2165/00019053-200017060-00004>.
12. Lucangelo U, Garufi G, Marras E, Ferluga M, Turchet F, Bernabè F, Comuzzi L, Berlot G, Zin WA. End-tidal versus manually-controlled low-flow anaesthesia. *J Clin Monit Comput*. 2014;28(2):117–21. <https://doi.org/10.1007/s10877-013-9516-8>.
13. Feldman JM, Hendrickx J, Kennedy RR. Carbon dioxide absorption during inhalation anesthesia: a modern practice. *Anesth Analg*. 2021;132(4):993–1002. <https://doi.org/10.1213/ANE.0000000000005137>.
14. Faul F, Erdfelder E, Buchner A, Lang AG. Statistical power analyses using G\*Power 3.1: tests for correlation and regression analyses. *Behav Res Methods*. 2009;41(4):1149–60. <https://doi.org/10.3758/BRM.41.4.1149>.
15. Ripple WJ, Wolf C, Newsome TM, Barnard P, Moomaw WR. World scientists' warning of a climate emergency. *Bioscience*. 2020;70(1):8–12. <https://doi.org/10.1093/biosci/biz088>.
16. Kennedy R, Hendrickx J, Feldman J. There are no dragons: lowflow anaesthesia with sevoflurane is safe. *Anaesth Intensive Care*. 2019;47(3):223–5. <https://doi.org/10.1177/0310057X19843304>.
17. Kharasch E, Frink E, Zager R, Bowdle T, Artru A, Nogami W. Assessment of low-flow sevoflurane and isoflurane effects on renal function using sensitive markers of tubular toxicity. *Anesthesiology*. 1997;86:1238–53. <https://doi.org/10.1097/00000542-199706000-00004>.
18. Gradus RHJM, Nillesen PHL, Dijkgraaf E, Van Koppen RJ. A cost-effectiveness analysis for incineration or recycling of dutch household plastic waste. *Ecol Econ*. 2017;135:22–8. <https://doi.org/10.1016/j.ecolecon.2016.12.021>.
19. Climate change 2007: Working Group III: mitigation of climate change. [https://archi.ve.ipcc.ch/publi/catio/ns\\_and\\_data/ar4/wg3/en/ch7s7-4-5-2.html](https://archi.ve.ipcc.ch/publi/catio/ns_and_data/ar4/wg3/en/ch7s7-4-5-2.html). Accessed 23 Dec 2021
20. Bressler RD. The mortality cost of carbon. *Nat Commun*. 2021;12(1):4467. <https://doi.org/10.1038/s41467-021-24487-w>.
21. De Medts R, Carette R, De Wolf AM, Hendrickx JFA. Desflurane usage during anesthesia with and without N<sub>2</sub>O using FLOW-i automatic gas control with three different wash-in speeds. *J Clin Monit Comput*. 2018;32(4):763–9. <https://doi.org/10.1007/s10877-017-0031-1>.
22. De Cooman S, De Mey N, Dewulf BB, Carette R, Deloof T, Sosnowski M, De Wolf AM, Hendrickx JF. Desflurane consumption during automated closed-circuit delivery is higher than when a conventional anesthesia machine is used with a simple vaporizer-O<sub>2</sub>-N<sub>2</sub>O fresh gas flow sequence. *BMC Anesthesiol*. 2008;17(8):4. <https://doi.org/10.1186/1471-2253-8-4>.
23. Young CC, Harris EM, Vacchiano C, Bodnar S, Bukowy B, Elliott RRD, Migliarese J, Ragains C, Trethewey B, Woodward A, Gama de Abreu M, Girard M, Futier E, Mulier JP, Pelosi P, Sprung J. Lung-protective ventilation for the surgical patient: international expert panel-based consensus recommendations. *Br J Anaesth*. 2019;123(6):898–913. <https://doi.org/10.1016/j.bja.2019.08.017>.
24. Sherman J, Le C, Lamers V, Eckelman M. Life cycle greenhouse gas emissions of anesthetic drugs. *Anesth Analg*. 2012;114(5):1086–90. <https://doi.org/10.1213/ANE.0b013e31824f6940>.

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[www.getinge.com/jp](http://www.getinge.com/jp)

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