
Surgical light

Operacinis apšvietimas

User's manual **EN**

Naudojimo instrukcija **LT**

CONTENTS

| | |
|--|-----------|
| INTRODUCTION | 3 |
| QUALITY STANDARDS COMPLIANCE | 4 |
| SYMBOLS USED ON PRODUCT | 5 |
| SYMBOLS USED IN MANUAL | 5 |
| WARNINGS | 6 |
| 1 INSTRUCTIONS FOR SAFE USE | 8 |
| 2 DESCRIPTION | 10 |
| 2.1 Discharge lamps | 10 |
| 2.2 Electronic ballast | 10 |
| 2.3 Voltage source | 10 |
| 2.4 Specific light output level variation features | 11 |
| 2.5 Turning the light on and off (starter mechanism) | 11 |
| 2.6 Conditions of use | 11 |
| 3 USE | 13 |
| 3.1 Before starting | 13 |
| 3.2 Turning on the light | 13 |
| 3.3 Positioning the light | 15 |
| 3.4 Pre-positioning examples | 16 |
| 3.5 Setting the light output level | 19 |
| 3.6 Indicators on control keypad | 19 |
| 3.7 Turning off the cupola | 19 |
| 4 CLEANING / DISINFECTION / STERILISATION | 20 |
| 4.1 Cleaning and disinfecting the surgical light | 20 |
| 4.2 Cleaning and sterilising the handles | 20 |
| 5 MAINTENANCE | 22 |
| 5.1 Preventive maintenance | 22 |
| 5.2 First level maintenance | 22 |
| 5.3 Recycling | 23 |
| 6 CHANGING BULBS | 25 |
| 6.1 Removing bulbs | 25 |
| 6.2 Inserting bulbs | 27 |
| 6.3 Testing failover to the spare bulb | 27 |
| 7 GENERAL CHARACTERISTICS | 28 |
| 8 ACCESSORIES | 29 |
| 9 ECM DECLARATION | 31 |
| 10 TROUBLESHOOTING | 37 |

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ĮVADAS

Dear customer,

Thank you for choosing the G8 Evolution surgical light.

MAQUET has designed this product line to ensure the most comfortable working conditions for surgical teams and optimal safety conditions for patients.

The G8 Evolution follows on from its illustrious predecessor, offering even more comfort and efficiency.

- new discharge bulb with ceramic burner, for improved colorimetry and service life
- new motorised bulb holder, enabling reversible switching between the main and backup bulbs
- new electronics, enabling motorised bulb holder control
- new optical filters, reducing the energy radiation particularly in backup mode
- new mirrors, identical for improved spot definition
- new hinged aluminium bulb access hatch, for an improved finish and ergonomics when replacing bulbs
- new product colour for compatibility with the new accessories.

Please read this manual thoroughly in order to fully understand the features of your G8 Evolution light and maintain it properly.

Do not hesitate to contact MAQUET if you have any questions about the G8 Evolution surgical light or our other products.

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QUALITY STANDARDS COMPLIANCE

Certification of MAQUET SA quality system

LNE/G-MED certifies that the quality system developed by MAQUET SA for design, implementation, sales, installation and after-sales service of surgical lights complies with the requirements of the following international standards:

- ISO 9001:2000
- NF EN ISO 13485:2004

Reference standards

The HLX 3000 surgical light was designed to comply with the following standards:

- EN ISO 14971:2000 Medical devices - Application of risk management to medical devices (ISO 14971:2000)
- EN ISO 14971:2000/A1:2003
- EN 60601-1:1990 Medical electrical equipment - Part 1: General requirements for safety
Amendment A1:1993 to EN 60601-1:1990
Amendment A2:1995 to EN 60601-1:1990
Amendment A13:1996 to EN 60601-1:1990
- EN 60601-1-2:2001 Medical electrical equipment - Part 1-2: General requirements for safety - Collateral standard: Electromagnetic compatibility - Requirements and tests
- EN 60601-1-4:1996 Medical electrical equipment - Part 1-4: General requirements for safety - Collateral standard: Programmable electrical medical systems
Amendment A1:1999 to EN 60601-1-4:1996
- EN 60601-1-6:2004 Medical electrical equipment - Part 1-6: General requirements for safety - Collateral standard: Usability
- EN 60601-2-41:2000 Medical electrical equipment - Part 2-41: Particular requirements for the safety of surgical luminaires and luminaires for diagnosis

CE Marking/Intended use

Compliance with the requirements of Directive 93/42/EEC relating to medical devices has been assessed and approved by LNE/G-MED. The G8 Evolution surgical light is a class I device in accordance with Annex IX of Directive 93/42/EEC.

The range includes single, double or triple ceiling-mounted configurations, with or without a pre-installed video system (HANAUVISION). Possible options include a double fork, flat screen support and separate camera support.

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WARNINGS

 Light energy can potentially dry tissue. Users must adapt lighting levels to the needs of the operation to be performed, in particular when several cupolas are used in combination.

 Do not look directly at the light source due to its high intensity.

 In the event of a power failure, only cupolas connected to a backup power supply system will remain operational.

 When changing a blown bulb:

- Switch off the power supply and leave the cupola to cool with the cover open for 45 minutes.
- Only use the recommended model of MAQUET bulb.
- Handle bulbs carefully using a clean, dry cloth.
- Never touch bulbs with bare hands.

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 Before each operation check that the sterilisable handle (light colour on black handle) is in place. During operations this handle must only be used by the surgical team in order to guarantee aseptic conditions.

 After each sterilisation and before using the sterilisable handle again:

- Check that there are no cracks.
- Check that the handle operates correctly on the light.

 Do not use the unit in the presence of inflammable anaesthetic gases.

 The unit must be operated with the cover closed. When conducting maintenance, take appropriate precautions to avoid touching surfaces labelled as hot.

 Dismantling certain elements may affect operation and safety. Examples include:

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 Do not use the surgical light suspension arm to suspend or lift objects. Do not hang on the light.

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Do not use in an MRI room.

Certain light wavelengths may be incompatible with certain pathologies.

Fumigation methods are unsuitable for disinfecting the unit and must not be used.

In the event of a mains power failure, only the cuspolas with power supplies connected to a backup power supply system will remain operational.

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1 INSTRUCTIONS FOR SAFE USE

Please note that certain operations may only be carried out by staff with the appropriate qualifications:

- **The ceiling unit may only be operated by trained medical staff.**
- This unit uses state of the art technology and is completely safe. Nevertheless dangers can arise, especially when it is operated by unqualified staff, when it is not well operated or not used as intended.
- **Cleaning of the appliance must be carried out by qualified sanitary specialists.**
- Please read these operating instructions carefully before using the unit. You will then enjoy all the advantages it can offer and protect yourself and others from harm.
- Please follow these instructions to train staff in the operation, care and safe use of the unit.
- The unit is only intended for use as mentioned in the operating instructions and is only suitable for such use. Any other use can be dangerous to life and limb and/or cause damage to the product or to other property of the operator.
- For safety reasons, no alterations or adaptations may be made on your own initiative without the agreement of MAQUET.
- For your safety, contact your supplier should any problems arise that are not fully addressed by these instructions.
- The contents of the user manual may be changed at any time without notice.

Additional equipment

- The models named in this document may be equipped with other manufacturers' end appliances (e.g. monitors). For information on operating them, please refer to each manufacturer's user manuals.

Transport and storage

- Ambient temperature: -25°C to 70°C
- Relative humidity: 10 to 75%
- Atmospheric pressure: 500 - 1060 hPa
- Do not store outside.
- Do not subject to severe vibrations.

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Operating conditions

- Ambient temperature: 10°C to 40°C
- Relative humidity: 30 to 75%
- Atmospheric pressure: 700 - 1060 hPa
- The appliance is not suitable for use in explosive atmospheres.
- G8 Evolution must be installed and commissioned in accordance with the EMC information provided on page 30.
- Portable RF communication devices may affect the operation of this equipment.
- The unit is designed for continuous operation.

Maintenance

- Maquet's warranty is only valid and the safety and integrity of the product are only guaranteed if:
 - All inspection, maintenance and repair operations are performed by Maquet engineers or trained and authorised technical support technicians.
 - Only original accessories, consumables and spare parts are used.

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2 DESCRIPTION

2.1 Discharge lamps

The G8 Evolution surgical light is equipped with a discharge lamp for the main bulb and a halogen bulb for the spare bulb.

The use of a discharge lamp entails significant differences in operation compared with halogen lamps. The key parts are an electronic ballast housed in the cupola and a starter mechanism in the bulb holder.

2.2 Electronic ballast

The electronic ballast powers the bulb, ensuring the correct power supply voltage, automatically powering the spare bulb if the main bulb fails, counting the operating time of the main bulb and communicating with the control keypad.

2.3 Voltage source

The technical properties of discharge lamps are significantly different to those of halogen lamps. Halogen lamps are powered by a transformer at a constant voltage, whereas for discharge lamps the operating voltage varies depending on the operating conditions. The voltage supplied to the bulb is typically 20 V during the cold start-up phase, 80 V during normal operation, with peaks up to 20 kV during firing. In standby operation the ballast generates an output voltage of about 400 V.

For safety reasons, the entire power supply circuit is automatically cut off when the bulb change cover is open.

The voltage supplied to the discharge lamp in G8 Evolution units is controlled by the ballast.

This implies that the light output level provided by the G8 Evolution does not directly depend on the supply voltage – a major difference compared with conventional surgical lights.

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2.4 Specific light output level variation features

The light output level provided by the G8 Evolution is controlled solely by the ballast. Conventional external dimmers which reduce the power supply voltage will not work with this light unit as the ballast compensates the resulting voltage. Using an external dimmer may even be dangerous; if the supply voltage drops below the value which the ballast can compensate the discharge lamp will go out completely!

Accordingly the voltage supplied to the tube must always be > 24 V.

Discharge lamps have a slight reaction delay when the light output level is reduced using the keys on the control keypad. Once the light level has been adjusted, the chromaticity of the light differs somewhat from the full output level. Unlike halogen lamps, which contain a higher proportion of red light when their output level is reduced using a variator, the chromaticity of discharge lamps shifts towards green.

This difference is particularly noticeable if two G8 Evolution lights placed side by side are set to different light output levels.

It is therefore recommended that the light output level for G8 Evolution lights used in combination should be set to the same level.

2.5 Turning the light on and off (starter mechanism)

Discharge bulbs are lit by a trigger located on the bulb holder. They take several minutes to reach their intended brightness and colour.

Coloured flashes and strong colour deviations are completely normal when the lights are first turned on.

Unlike many other light units using discharge lamps, G8 Evolution surgical lights may be turned on even when hot and immediately reach their final brightness and colour.

A control keypad on the fork is used to turn the lamp on and off and vary its brightness. The keypad also comprises four indicator lights which show the service life of the main bulb. To avoid the light being turned off unexpectedly the Off button must be pressed and held for several seconds to turn off the lamp.

The service life of the bulb would be seriously reduced if turned off during the cold start-up phase. Accordingly, to protect the bulb, the light may only be controlled from the keypad once the bulb reaches its normal operating condition, after a hot or cold start-up (controls locked out for approx. 1 minute).

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
Suspension system:

- 1 Ceiling anchor plate with cover
- 2 Suspension tube and flange
- 3 Extension arm
- 4 Spring arm (standard or low ceiling)
- 5 Double fork (only for Comfort version)
- 6 Fork

Control keypad:

- 7 Bulb fail indicator
- 8 Lifecycle indicator (bulb service life)
- 9 Warning indicator (lit steadily to indicate a faulty main bulb or electronics unit)
- 10 ON / OFF button (press and hold for 2 seconds to turn the light off)
- 11 Light output level indicators (flash in turn approx. 60 seconds after start-up)
- 12 +/- keys for adjusting the light output level

Cupola:

- 13 Cupola
- 14 Round handle or circular strip (non-sterile)
- 15 Sterilisable handle
- 16 Reflector
- 17 Protective underside
- 18 Camera
- 19 Dummy camera (weight to compensate for optional camera)
- 20 Flat screen support
- 21 Rotation stop

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3 USE

3.1 Before starting

1. Insert the sterilisable handle (15). A click confirms that the handle is correctly in place.

-  Check the following points before turning on the cupola (13):
1. Protective underside (17) not cracked or crazed.
 2. Parts correctly attached to underside (17),
 3. Handle (15) correctly inserted.

If any of the above items are defective, the surgical light cannot be operated safely:

1. Turn off the power to all poles of the surgical light and ensure that it cannot be restarted.
2. Call your MAQUET service partner.

 Before each operation, check that the sterilisable handle (light colour on black handle) is in place. During operations this handle must only be used by the surgical team in order to guarantee aseptic conditions.

3.2 Turning on the light

1. Press the ON / OFF button (10).
 - The lamp is turned on.
 - **Cold start:** After a cold start, the indicators (11) flash until the operating temperature is reached (approx. 60 seconds).
 - The lifecycle indicators (8) are lit (see section 2.6).
 - The indicators (11) show the selected light output level (see section 2.5).
2. Orient the cupola as described in section 2.4.
3. Adjust the light output level as described in section 2.5.

 **Warning indicator (9):** If the warning indicator (9) is lit steadily, there is a fault on the main bulb or the electronic unit. The secondary halogen bulb is turned on. This provides a backup system to enable the operation to be completed. It is essential to restore normal operation as described in section 5 before operating another patient.

- The cupola may hum slightly until it reaches a stable state (5 to 10 minutes).

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
Suspension system:

- 1 Suspension tube
- 2 Extension arm
- 3 Spring arm
- 4 Spring arm (low ceiling)

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
3.3 Positioning the light

The cupola (5) on the upper extension arm (2) may be positioned directly beneath the suspension tube (1) as required for certain types of surgery.

Optimum mobility is achieved in an «M» position, with the extension arms (2) and spring arms (3) positioned as shown on the illustration.

Light energy can potentially dry tissue. Users must adapt lighting levels to the needs of the operation to be performed, in particular when several lights are used in combination.

- Use the sterilisable handle (7) or the non-sterile round handle or bar (6) to move the cupola (5) to the desired position.

The working distance is 70 to 140 cm from the surgical site.

Do not use the surgical light suspension arm to suspend or lift objects. Do not hang on the light.

The light must not be moved past the vertical position (cupola pointing upwards).

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

- 1 Bulb fail indicator
- 2 Lifecycle indicator (bulb service life)
- 3 Warning indicator (lit steadily to indicate a faulty main bulb or electronics unit)
- 4 ON / OFF button
(press and hold for 2 seconds to turn the light off)
- 5 Light output level indicators
(flash in turn approx. 60 seconds after start-up)
- 6 + / - buttons to reduce/increase light output level

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3.5 Setting the light output level

- The control buttons are not available when the indicators (5) are flashing.
- Due to the properties of discharge lamps the light does not react immediately when the light output level is reduced with the + / - buttons (6). At lower light levels, the colour of the light changes (higher colour temperature).
- It is preferable to select the same light output level on all cupolas. Colour reproduction is most faithful at the maximum output level.
- The light output may be adjusted to one of four levels using the + / - buttons (6); there are four indicators (5) corresponding to the four levels.
- The + / - buttons (6) are only active once the light has reached its service temperature.

Adjusting the diameter of the illuminated area:

- The diameter of the illuminated area (focus) is adjusted by turning the sterilisable handle.

3.6 Indicators on control keypad

Lifecycle indicator (2):

- The main bulb has an average service life of around 4500 hours (see section 2). This service life is shown by the indicators (2). If the rated service life is exceeded, the bulb will continue to operate but will lose power and may no longer start.

Bulb fail indicator (1):

- When the indicator above Lifecycle is lit, the bulbs must be replaced in order to maintain the light output specifications (section 6).

Warning indicator (3):

- If the main bulb or electronic unit fail, the warning indicator (3) is lit continuously and the spare bulb is activated (see section 5.1).

Light output level indicator (5):

- When the light is first turned on (cold start-up), the light output level indicators (5) flash until the maximum intensity is reached (approx. 60 seconds); thereafter the four indicators (5) show the set light output level.

3.7 Turning off the cupola

- Press and hold the On / Off button (4) for 2 seconds.
- The light features a system to prevent it being turned off accidentally. Accordingly the On / Off button (4) must be held down for 2 seconds before the light is turned off.
- To protect the discharge lamp, power down is delayed if the lamp has just been turned on (see section 2).

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4 CLEANING / DISINFECTION / STERILISATION

Users must contact their hospital's sanitary specialists. The recommended products and procedures must be applied. Should there be any doubt concerning the compatibility of active agents to be used, contact the local MAQUET customer service.

4.1 Cleaning and disinfecting the surgical light

 Check that the power is switched off and the light has cooled down before starting cleaning.

4.1.a General instructions concerning cleaning, disinfection and safety:

- Remove the sterilisable handles.
- Wipe the system with a cloth moistened with a surface cleaner. Follow the manufacturer's dilution and temperature recommendations.
- Rinse with a cloth moistened with water. Wipe dry.
- Wipe evenly with a cloth moistened with disinfectant. Follow the manufacturer's recommendations.
- Rinse with a cloth moistened with water to remove any residue (particularly when cleaners containing aldehydes, quaternary ammonium or surfactants are used).
- Wipe off with a dry cloth.
- Make sure that all liquid cleaning products used have been thoroughly wiped off.

4.1.b Examples of recommended products

GETINGE USA products: TEC QUAT 256.

ANIOS products: SURFA'SAFE; HEXANIOS G + R at 0.5% (AMMONIUM IV, POLYHEXANIDE); ANIOSYME P.L.A; SALVANIOS pH10; ANIOS DDSH.

Schülke & Mayr products: ANTIFECT PLUS.

4.1.c Examples of unsuitable products

 Solutions containing glutaraldehyde, phenol, iodine, bleach, alcohol or chloride ions must not be used.

 Do not use disinfection by fumigation methods.

4.2 Cleaning and sterilising the handles

4.2.1 Before cleaning

- Use a soft cloth immediately after use to wipe away soiling from the handle surface.
- Store handles in a place that keeps them moist to make further cleaning easier.
- Take care to store them in such a way that the inside does not get soiled.

4.2.2 Cleaning

- Soak the handles in a detergent solution.¹
 - Soak for 15 minutes to allow the solution to act, then clean by hand with a soft brush and a lint-free cloth.
- ¹ Never soak the handles in enzyme-based detergents as they may damage the material used; rinse thoroughly if these detergents are used.
- During cleaning, check regularly that the handles are fully clean and that no soiling remains on the inside or outside.
 - If any soiling remains, repeat cleaning or use an ultrasonic cleaning process.
 - Rinsing: Rinse thoroughly in clean water to completely eliminate the detergent solution.
 - Drying: Wipe with a clean lint-free cloth.

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4.2.3 Disinfection

Handles may be disinfected by machine (Clean MAQUET) and rinsed at a maximum temperature of 93°C.

Typical recommended cycles:

4.2.4 Sterilisation

After cleaning the handles must be steam sterilised as set out below:

- Check that each handle is clean before continuing the process.
- Wrap the handles with sterilisation wrapper material (double wrapper or equivalent).

Handles may also be placed in paper or plastic sterilisation bags for easier identification and reuse.

| Stage | Temperature | Time |
|----------------|-------------|---------|
| Pre-washing | 18 - 35° C | 60 sec |
| Washing | 46 - 50° C | 300 sec |
| Neutralisation | 41 - 43°C | 30 sec |
| Washing 2 | 24 - 28°C | 30 sec |
| Rinsing | 92 - 93°C | 600 sec |
| Drying | | 20 min |

- Place the handles on steriliser trays with the opening downwards.

| Country | Sterilisation cycle | Temperature [°C] | Time [min] | Drying [min] |
|-----------------|--------------------------|----------------------------------|------------|----------------|
| USA & Canada | Prevacuum ² | 132 - 135 | 10 | 16 |
| France | ATNC (Prion) (Prevacuum) | 134 | 18 | |
| Other countries | Prevacuum | Comply with national regulations | | |

- Package with biological and/or chemical indicators for monitoring the sterilisation process. Follow applicable regulations.
 - Run the sterilisation cycle according to the steriliser manufacturer's instructions.
-  To ensure correct sterilisation do not allow any soiling to penetrate inside the handle.
-  Handles are not guaranteed beyond 350 sterilisation cycles with the above sterilisation parameters.
- Dispose of the handles in the same way as other hazardous products in a hospital environment.

² This handle is made of a porous material.
³ Possible sterilisation bag suppliers:
Medical Action Industries
SBW Medical
Baxter International
⁴ For air removal and faster drying.

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5 MAINTENANCE

5.1 Preventive maintenance

To preserve your surgical light's original performance and reliability, annual maintenance and inspections should be performed as follows:

- by a MAQUET technician or MAQUET-approved distributor during the guarantee period,
- by a MAQUET technician or MAQUET-approved distributor or by the hospital's technical maintenance department after the guarantee period expires,

5.2 First level maintenance

5.2.a Daily inspection

- Check that the bulbs operate correctly.
- Check that the sterilisable handle clicks and locks in place correctly.
- Check that the arms move normally and do not sag.

5.2.b Monthly inspection

- Check that the light operates correctly from the backup power supply if a power cut were to occur.
- Check that the limit stops are in place.
- Test the fallback to the spare bulb (see section 5.3).

5.2.c Annual inspection (must be performed by an authorised technician)

- Check the bulb holder. Replace the bulbs.
- Check that the limit stops are in place on the flat screen support arms and cupolas.
- Clean the internal optical items (mirrors and lenses).
- Check that the light operates for one hour if the main power supply is cut (operation on backup supply).
- Replace the fan subassembly.

Safety items

Check the following points:

- Check that there are no cracks on the upper pivot of the DF Acrobat 3000 spring arm ref. **568101988** (if manufactured between 2001 and 2006).
- If there is the slightest sign of a crack, replace the spring arm.
- If no cracks are visible, install kit ref. 368104900.
- Attachment screws on suspension tube correctly tightened, seals in position.
- Arm(s) correctly assembled..
- Limit stop segments in place on cupola attachment points. Disassemble cupolas and lubricate sleeve.
- Attachment of all covers and caps on cupolas and arms.
- All visible screws correctly tightened.
- Camera or dummy unit firmly attached.

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Other checks

- Nominal light output level: see technical data.
- Earth continuity: max. 0.1 Ohm
- Suspension tube vertical.
- Balancing system adjusted correctly: brakes tight and spring arms adjusted.
- Sterilisable handle locking mechanism.
- Adjustment of rotation limit stops.

Appearance

Check the following points:

- General condition of arm covers.
- Condition of lenses.
- Condition of seals.
- Suspension system clean.

Dismantling certain elements may affect operation and safety.

Examples include:

- Servicing the electrical power supply,
- Servicing the suspension arm and balance system,
- Servicing the optical system of cupolas equipped with filters designed to eliminate radiation not visible to the patient. Surgical lights must never be used without these filters.

Contact the authorised MAQUET after-sales service department for this type of inspection.

5.3 Recycling

Discharge bulbs, like commercial fluorescent tubes, contain a proportion of mercury and must be disposed of as special waste.

Halogen bulbs can be disposed of as domestic waste.

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



1 2

3 4

- 1 Cover lock
- 2 Cover
- 3 Connectors
- 4 Bulb holder attachment screw
- 5 Bulb holder
- 6 Main bulb (discharge type)
- 7 Spare bulb (halogen type)
- 8 End of main bulb

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6 CHANGING BULBS

When changing a blown bulb:

- Switch off the power supply and leave the cupola to cool with the cover open for 45 minutes.
- Only use G BULB EVOLUTION bulbs.
- Handle bulbs carefully using a clean, dry cloth.
- Never touch bulbs with bare hands.

Grease on bulbs can shorten their life or even cause them to break.

6.1 Removing the main bulb and spare bulb from the cupola

Main bulb and spare bulb

- The cupola is equipped with one main bulb (6) (discharge type) and one spare bulb (7) (halogen type).
- If the main bulb (6) fails, the spare bulb (7) is automatically turned on. The bulb holder is moved to place the spare bulb (7) at the centre of the reflector.
- When the cupola is turned off, the main bulb automatically returns to the centre.

Removing the bulbs

1. Unscrew the cover lock (1) and open the cover (2).
2. Allow the cupola to cool for 45 minutes.
3. Disconnect the connectors (3).
4. Unscrew the two bulb holder mounting screws (4).
5. Carefully withdraw the bulb holder (5) from the cupola.
6. Grasp the end of the main bulb (8) and push the bulb in to release it from the spring (A).
7. Pull out one side of the bulb, then the other (B).
8. Remove the spare bulb (7) by pulling straight down.

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


- 1 Cover lock
- 2 Cover
- 3 Connectors
- 4 Bulb holder attachment screws
- 5 Bulb holder
- 6 Main bulb
- 7 Spare bulb
- 8 End of main bulb
- 9 Lifecycle indicator (bulb service life)
- 10 Warning indicator (lit steadily to indicate a faulty main bulb or electronics unit)
- 11 ON / OFF button (press and hold for 5 seconds to turn the light off)
- 12 Light output level indicators (flash in turn approx. 60 seconds after start-up)
- 13 + / - buttons to reduce/increase light output level

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6.2 Inserting the main bulb and spare bulb

To avoid any failure of the surgical light during use, always replace the spare bulb at the same time as the main bulb.

1. Insert the spare bulb (7) firmly into its socket.
2. Insert the main bulb (8) by pressing in one end, then releasing it to hold the other end. Check that it is correctly in place.
3. Carefully insert the bulb holder (5) into the cupola, ensuring that the male connectors (3) on the bulb holder are on the same side as the female connectors in the cupola (3).
4. Connect the two connectors (3).
5. Tighten the two bulb holder attachment screws (4).
6. Close the cover (2) and tighten the cover locking screw (1).
7. Check that the cover (2) is correctly in place.

The unit must be operated with the cover closed. When conducting maintenance, take appropriate precautions to avoid touching surfaces labelled as hot.

8. With the light off, press the + and - buttons (13) simultaneously and hold them for 5 seconds to reset the Lifecycle counter (back to 0 hours).
9. Press the ON / OFF button to turn on the light.
 - Only one indicator (Lifecycle) (9) on the control keypad should be lit. If this is not the case, repeat from operation 8.
 - After a cold start, the indicators (12) flash in turn until the operating temperature is reached (approx. 60 seconds).
 - If the warning indicator (10) is lit steadily, check that the main bulb (6) and spare bulb (7) are correctly in place.
 - If the problem persists please contact the technical department.

6.3 Testing failover to the spare bulb

This test is used to check that the failover function operates correctly.

1. Turn on the light by pressing the ON / OFF button.
2. Once the main bulb is stable (after 60 seconds) reduce the output to the minimum level by pressing the - button (13).
3. Hold down the - button (13) for 20 seconds until the spare bulb is selected (yellow light).
4. Check that the light beam from the spare bulb is satisfactory.
5. Turn off the surgical light by pressing the ON / OFF button (11) and holding for 2 seconds.
6. Turn the light back on and check that the main bulb comes on.

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7 GENERAL CHARACTERISTICS OF G8 EVOLUTION SURGICAL LIGHTS
(In accordance with standard IEC 60,601-2-41)

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* Factory set value

Note:

- The tolerated values are guaranteed on purchase of the product.
- The non-tolerated values are measured on a sample coming out of production.

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9 ELEKTROMAGNETINIO SUDERINAMUMO DEKLARACIJA

(In accordance with standard EN 60601-1-2, November 2001 edition)

| Table 201 - Guidance and manufacturer's declaration - electromagnetic emissions | | |
|--|----------------|--|
| G8 Evolution surgical lights are intended for use in the electromagnetic environment specified below. G8 Evolution customers or users should ensure that they are used in such an environment. | | |
| Immunity test | Compliance | Electromagnetic environment — directives |
| RF emissions CISPR 11 | Group 1 | G8 Evolution lights use RF energy only for their internal functions. Their RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. G8 Evolution lights are suitable for use in all establishments other than domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| RF emissions CISPR 11 | Class A | |
| Harmonic emissions IEC 61000-3-2 | Not applicable | |
| Voltage fluctuations and flicker emissions IEC 61000-3-3 | Not applicable | |

| Table 202 - Guidance and manufacturer's declaration - electromagnetic immunity | | | |
|---|---|---|--|
| G8 Evolution surgical lights are intended for use in the electromagnetic environment specified below. G8 Evolution customers or users should ensure that they are used in such an environment.. | | | |
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment — guidance |
| Electrostatic discharge (ESD) IEC 61000-4-2 | ± 6 kV contact ± 8 kV air | ± 6 kV contact ± 8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient / burst IEC 61000-4-4 | ± 2 kV for power supply lines ± 1 kV for input/output lines | ± 2 kV for power supply lines ± 1 kV for input/output lines | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | ± 1 kV differential mode ± 2 kV common mode | ± 1 kV differential mode ± 2 kV common mode | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11 | < 5% U_T (dip > 95% of U_T) for 0.5 cycles 40% U_T (dip = 60% of U_T) for 5 cycles 70% U_T (dip = 30% of U_T) for 25 cycles < 5% U_T (dip = 95% of U_T) for 5 cycles | < 5% U_T (dip > 95% of U_T) for 0.5 cycles 40% U_T (dip = 60% of U_T) for 5 cycles 70% U_T (dip = 30% of U_T) for 25 cycles < 5% U_T (dip = 95% of U_T) for 5 cycles | Mains power quality should be that of a typical commercial or hospital environment. If the G8 Evolution requires continuous operation during mains interruptions, it is recommended that it be powered from an uninterruptible power supply. |
| Power frequency (50/60 Hertz) magnetic fields IEC 61000-4-8 | 3 A | 3 A | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |
| Note: U_T is the voltage of the AC supply network before applying the test level. | | | |

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Table 204 - Guidance and manufacturer's declaration - electromagnetic immunity

G8 Evolution surgical lights are intended for use in the electromagnetic environment specified below. G8 Evolution customers or users should ensure that they are used in such an environment.

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic — guidance environment |
|-------------------------------|------------------------------|------------------|---|
| Conducted RF IEC 61000-4-6 | 3 V rms 150 kHz to 80 MHz | 3 V rms | <p>Portable and mobile RF communications equipment should be used no closer to any part of the BLUE 30/80, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[\frac{3.5}{\sqrt{f}} \right] \sqrt{P} = 1.17 \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \text{ from } 80 \text{ MHz to } 800 \text{ MHz} = 1.17 \sqrt{P}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \text{ from } 800 \text{ MHz to } 2.5 \text{ GHz} = 2.34 \sqrt{P}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic survey on site^a, should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |
| Radiated RF IEC 61000-4-3 | 3 V/m 80 MHz to 2.5 GHz | 3 V/m | |

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the G8 Evolution is used exceeds the applicable RF compliance level above, the G8 Evolution should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the G8 Evolution.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 206 - Recommended separation distance between portable and mobile RF communication devices and G8 Evolution lights.

G8 Evolution is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the G8 Evolution can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the G8 Evolution as recommended below, according to the maximum output power of the communications equipment.

| Maximum rated output power of transmitter W | Separation distance depending on transmitter frequency m | | |
|--|---|--|---|
| | 150 kHz to 80 MHz $d = [1.17] \sqrt{P}$ | 80 MHz to 800 MHz $d = [1.17] \sqrt{P}$ | 800 MHz to 2.5 GHz $d = [2.34] \sqrt{P}$ |
| 0.01 | 0.12 | 0.12 | 0.24 |
| 0.1 | 0.37 | 0.37 | 0.74 |
| 1 | 1.17 | 1.17 | 2.34 |
| 10 | 3.69 | 3.69 | 7.38 |
| 100 | 11.70 | 11.70 | 23.40 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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10 TROUBLESHOOTING

| Anomaly | Likely cause | Corrective action |
|---|--|---|
| Cupola drifts | Suspension tube not vertical | Check tube verticality and ceiling structure |
| Cupola too stiff or does not stay in position | Balancing incorrectly adjusted | Adjust the balancing |
| | Ceiling structure unstable | Call MAQUET technical department |
| | Brakes not tight enough | Adjust the brakes |
| Light flickers (normal during cold start-up) | Poor contact, problems with electronics or ageing of main bulb | Call MAQUET technical department |
| Insufficient light output level | Main bulb faulty, electronic unit faulty, low light output selected | Check bulbs and replace if necessary, check bulb position, call technical department |
| Light beam not round | Incorrect type of bulb(s) | Check bulbs and replace if necessary |
| | Incorrect bulb or bulb holder position | Check bulb position, call technical department |
| Bulb service life too short | Non-compliant bulbs used | Check that the bulbs recommended by MAQUET are used |
| Lifecycle indicator fully lit even after replacing bulb | Lifecycle indicator has not been reset | Reset the indicator |
| Greenish light | Main bulb too old, main bulb faulty | Check bulbs and replace if necessary, check bulb position |
| | Light output level too high | Select a lower light output level |
| Cupola does not turn off | Off button not held down long enough | Press the Off button and hold for at least 2 seconds |
| | Cupola has just been turned on | Wait until the light reaches its service temperature (60 seconds) before turning it off |
| Cupola does not turn on | Power supply cut | Check power supply and fuses |
| | No bulb, bulb incorrectly installed | Lampen überprüfen und ggf. ersetzen, Lampensitz prüfen |
| | Electronic unit faulty | Call technical department |
| Cupola cover crazed or cracked | Use of too aggressive cleaning solution, cleaning or disinfection when hot | Use a different cleaning solution |
| Service life of sterilisable handles too short | Sterilisation procedure too harsh | Check sterilisation procedure |
| Sterilisable handles damaged or cracked | End of sterilisable handle service life | Replace sterilisable handles |

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Our network is on hand to answer any of your questions.

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