**Axcel**

**AXL 5001, AXL 5501, AXL 5002**


**Éclairage Interventionnel**  
***Surgical Light System***  
**Corpuri de iluminat chirurgicale**

---

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*User's manual* **EN**

Manual de utilizare **RO**

## 1 - INTRODUCTION

Cher client,

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## 1 - INTRODUCTION

Dear customer,

*Thank you for choosing the AXCEL Operation lighting system.*

*This manual is for use of AXL 5001, AXL 5501 or AXL 5002 models.*

*MAQUET SAS has designed the range to ensure the most comfortable working conditions for the surgical team and optimum safety conditions for the patient.*

*AXCEL is the result of MAQUET SAS's ongoing contact with our customers and the previous experience of our engineering and design team. As a result, AXCEL ensures :*

- *Superb shadow dilution,*
- *Easy maneuverability*
- *Laminar flow design*
- *Incredible volume of light.*

*To get the best performances from the exceptional features of the AXCEL surgical light, please carefully read this manual.*

*Do not hesitate to contact MAQUET SAS with any questions you may have concerning AXCEL or our other products.*

*This document is the property of MAQUET SAS and may not be reproduced, in whole or in part, without our authorization.*

*This document was written with the assistance of our company's technical department in FRANCE. Your comments could be used to make corrections, and further details could be included on different installations you may come across to form an updated reference document for installing.*

*All correspondence should be addressed to:*

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**WARNINGS**

*Light energy can potentially dry tissue. The user must adapt the lighting level to the needs of the operation to be performed, particularly when several cupola are used in combination.*

*Light is a form of energy that, on account of certain wavelengths emitted, may not be suitable for certain pathologies.*

*Do not look directly at the light source due to its high intensity.*

*To fit bulb on bulb holder, push straight, up to limit stop.*

*To prevent any failures during operations, we recommend preventive replacement of the bulbs approximately every 600 hours to 800 hours, depending on the use.*

*After each sterilization and before each new use of the sterilizable handle :*  
*- Check that there are no cracks and that the lock operates correctly.*  
*- Fit the handle on the light and lock it by rotating ("Click").*

*Once the sterile handle is installed on the cupola, it should only be touched by the surgical team members.*

*Not suitable for use in the presence of a flammable anesthetics mixture with air or with oxygen or with nitrous oxide.*

*To avoid damaging the medical devices, do not use this light for purposes other than those stated in this document.*

*Do not use the suspension of the surgical light to carry or lift objects.*

*Do not hang on the surgical light.*

*Do not use cleaning solutions containing chloride ions.*

*Do not use disinfection by fumigation methods.*

*Disassembly of certain sub-assemblies can affect operation and safety, for example:*  
*- when servicing the electrical power supply,*  
*- when servicing the suspension arms and balancing system.*  
*- when servicing the optical system of the cupolas which are equipped with filters intended to eliminate radiations not visible by the patient. The unit should not be operated, under any circumstances, without these filters.*

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## WARNINGS

To avoid any risk of an electric shock, class I devices must be connected to a power supply system which is earthed.

Disconnect the batteries if the device is not scheduled to be used for more than three months.

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## Quality compliance

### Certification of MAQUET SAS quality system

LNE/G-MED certifies that the quality system created by MAQUET SAS for the design, manufacturing, marketing, installation and customer servicing of its surgical lights meets the requirements of the following international standards:

- ISO 9001:2000
- ISO 13485:2004.

### CE Marking

Compliance with the requirements of Directive 93/42/EEC dated of June 14th 1993, relating to medical devices has been assessed in accordance with Annex VII of this Directive. This **AXCEL™** Surgical lights range is a class I device in accordance with Annex IX of Directive 93/42/EEC.

### CE Marking affixed in 2007.

#### CE Marking products recognition :

AXCEL surgical light in single (ceiling mount AXL 5001 or wall mount AXL 5002) or double version (ceiling mount AXL 5501) for minor surgery, emergency, recovery, speciality, examination...

NOTE : With reference to installation rules complying with IEC 60 601-2-41 standard, in case of mains failure, AXCEL configurations can not be used in rooms where it is necessary to have emergency 24 V Direct current.

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## INSTRUCTIONS FOR SAFE USE

**Dear User**

*Please note that certain operations may only be carried out by staff with the appropriate qualifications :*

- ***The ceiling unit may only be operated by trained medical staff.***
- *This unit uses state of the art technology and is completely safe. Nevertheless dangers can arise, especially when it is operated by unqualified staff, when it is not well operated or not used as intended.*
- ***Cleaning of the appliance must be carried out by qualified sanitary specialists.***
- *Please read these operating instructions carefully before using the unit. You will then enjoy all the advantages it can offer and protect yourself and others from harm.*
- *Please follow these instructions to train staff in the operation, care and safe use of the unit.*
- *The unit is only intended for use as mentioned in the operating instructions and is only suitable for such use. Any other use can be dangerous to life and limb and/or cause damage to the product or to other property of the operator.*
- *For safety reasons, no alterations or adaptations may be made on your own initiative without the agreement of MAQUET.*
- *For your safety, contact your supplier should any problems arise that are not fully addressed by these instructions.*
- *The contents of the user manual may be changed at any time without notice.*

### Additional equipment:

- *The models named in this document may be equipped with other manufacturers' end appliances (e.g. monitors). For information on operating them, please refer to each manufacturer's user manuals.*
- *The ceiling service unit is designed for continuous operation.*

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**Transport and storage:**

- *Ambient temperature: -25°C to 70°C*
- *Relative humidity: 10 to 75%*
- *Atmospheric pressure: 500 - 1060 hPa*
- *Do not store outside.*
- *Do not subject to severe vibrations.*

**Operating conditions:**

- *Ambient temperature: 10°C to 40°C*
- *Relative humidity: 30 to 75%*
- *Atmospheric pressure: 700 - 1060 hPa*
- *The appliance is not suitable for use in explosive atmospheres.*
- *AXCEL must be installed and commissioned in accordance with the EMC information provided on page 29.*
- *Portable RF communication devices may affect the operation of this equipment.*

**Maintenance:**

- *Maquet's warranty is only valid and the safety and integrity of the product are only guaranteed if:*
  - *All inspection, maintenance and repair operations are performed by Maquet engineers or trained and authorised technical support technicians.*
  - *Only original accessories, consumables and spare parts are used.*

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## 2 - DESCRIPTION

*Because of our commitment to working with surgeons and nurses on a daily basis, we are able to identify and develop cost-effective solutions for the surgical light.*

**The innovative Axcel System key advantages are :**

**1) Improved shadow control for optimal visualization**  
*Axcel's original compact cross-shape cupola is an advantage since it can be positioned to direct the light beams around the surgeon's head, avoiding potential shadows with up to 100,000 lux illumination, Axcel is designed for a wide range of applications.*

**2) Uniform light distribution throughout the depth of field**  
*With Axcel, there is no need to refocus the light to illuminate deep cavities. Superior cavity penetration is based on four "external mirrors" that cross the individual light beams and one central lenses that provide 1 m (39.37 inches) depth of field.*

**3) Cool, comfortable lighting.**  
*Axcel is efficiently projecting useful (visible) light while eliminating the majority of the heat component. Light is projected onto the external dichroic coated mirrors where the infrared component is absorbed and is dissipated by conduction and convection.*

**4) True color rendition**  
*Axcel's color temperature was designed to reduce color distortion and reveals true tissue. Axcel incorporates a spectral filter system yielding a color temperature of 3,500 K and a high color rendition index (CRI) of 95.*

**5) Simple, durable design**  
*"Set it and forget it". The new Axcel surgical light was designed to be exceptionally durable and to require minimal maintenance.*

**6) Ease of positioning and movement**  
*Axcel is exceptionally light, combining mobility and stability. The small size of the Axcel cupola makes it easily maneuverable.*

**7) Cross-shaped design**  
*For laminar flow rooms, the aerodynamic shape of the Axcel cupola and its small area practically eliminate convection and turbulence effects, allowing the filtered air to circulate freely around the cupola.*

**8) Rounded, smooth surface for ease of cleaning**  
*The Axcel cupola is sealed and its surfaces are smooth and easy to clean.*

**9) Fast and easy bulb replacement**  
*The Axcel light bulb can be replaced quickly and easily without tools.*

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#### 4 - USE

Not suitable for use in the presence of a flammable anesthetics mixture with air or with nitrous oxide.

Do not use in an MRI environment.

Do not use the suspension of the surgical light to carry or lift objects.

Do not hang on the surgical light.

##### 4.1 - Switches on the Axcel cupola

Axcel is equipped with a multifunctional ergonomic handle to :

- Position the cupola
- Control the level of illumination and switch the lamp ON and OFF
- Protect against inadvertent contact with walls or other surgical equipment.

There is a membrane switch panel integrated into the light handle on the cupola yoke.

It includes three switches:

- The center switch toggles the lamp ON and OFF.
- The Dimmer control switches are located above and below the ON/OFF switch.

The plus sign will increase the illumination, while the minus sign will decrease the illumination level.

Light energy can potentially dry tissue. The user must adapt the lighting level to the needs of the operation to be performed, particularly when several cupola are used in combination.

Light is a form of energy that, on account of certain wavelengths emitted, may not be suitable for certain pathologies.

Do not look directly at the light source due to its high intensity.

In case of restoring of a power break, lighthoods switch on again automatically according to a level of median illumination whatever the original state.

##### 4.2 Moving and positioning

The surgical light must be prepositioned prior to any procedure to minimize subsequent handling.

By prepositioning the surgical light appropriately for each procedure, potential interactions with possible obstacles are limited (Serum holder, distribution arm, etc.).

Moving the cupola is achieved either through the sterilizable and removable central handle (1) or the external handle (2).

The dimming handle (3) is also designed to manipulate the unit.

Once the sterile handle is installed on the cupola, it should only be touched by the surgical team members.

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**DESCRIPTION**

*This sterilizable handle is suited for use with all Axcel surgical lighting lightheads.*

**USE OF STERILIZABLE HANDLE**

**INSTALLING AND REMOVING HANDLE**

Check handle's integrity (cracks, locking function, ...etc).

**a) Installing handle**

*Engage the sterilizable handle on its support to "click" it in to its locked position.*

*Turn the handle on its shaft until it clicks again.*

**b) Removing the handle**

*To remove the sterilizable handle, push and hold the button while pulling the handle downward to remove it from its support.*

**c) Option for once-only disposable handle, sterile type, DEVON®**

*To enable installation of sterile once-only disposable handles, DEVON® type, the handle support must be removed and replaced by a special adapter ref. DAX002.*

*For more details, refer to the installation procedure, supplied with DAX002 adapter.*

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## 5 - CLEANING AND MAINTENANCE

### CLEANING/DISINFECTION STERILIZATION

The user must contact the sanitary specialists at the medical center. The recommended products and procedures must be observed.

In the event of doubt concerning the compatibility of active agents to be used, contact the local MAQUET SAS customer service.

**5.1 Cleaning and disinfection of the surgical lighting system :**  
Before cleaning, ensure the power is turned off and the light is cooled down.

#### General instructions concerning cleaning, disinfection and safety :

- Remove the sterilisable handles.
- Wipe the system with a cloth moistened with a surface cleaner. Follow the manufacturer's dilution and temperature recommendations.
- Rinse the unit with a damp cloth and wipe.
- Wipe evenly with a cloth moistened with disinfectant. Follow the manufacturer's recommendations.
- Rinse with a cloth moistened with water to remove any residue (particularly when cleaners containing aldehydes, quaternary ammonium or surfactants are used).
- Wipe off with a dry cloth.
- Make sure that all liquid cleaning products used have been thoroughly wiped off.

#### a) Examples of recommended products

- **GETINGE USA Products**: TEC QUAT 256.
- **ANIOS Products**: SURFA'SAFE ; HEXANIOS G + Rat 0,5% ; ANIOSYME P.L.A. ; SALVANIOS pH10 ; ANIOS DDSH.
- **Schülke & Mayr Products**: ANTIFECT PLUS.

#### b) Examples of prohibited products

Solutions containing glutaraldehyde, phenol, iodine, bleach, alcohol or chloride ions must not be used.

Do not use disinfection by fumigation methods.

### 5.2 Sterilizable handles : cleaning and sterilization

#### 5.2.1 - Preparation for Cleaning

Use soft cloth to wipe away gross soils from light handle surface, immediately after its use.

Place handles in containment device for transporting and keep soils moist to ease further cleaning.

Avoiding situations that soil the inside surfaces of handles. This can help ease the required cleaning efforts.

Do not scratch the lens (video handles only).

#### 5.2.2 - Cleaning

Submerge handles while soaking and washing with a detergent<sup>1</sup> solution.

Let handles soak for 15 minutes and then manually wash soils off, using soft cloth and soft brushes.

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### 5.2.3 - Inspect Cleanliness of Surfaces

*Inspect handle periodically to assure all soils are removed from outer and also inner surfaces. If soils are seen on internal parts, then repeat cleaning or use an ultrasonic cleaning process, as necessary.*

#### Rinsing:

Thoroughly rinse handle using a brisk stream of clean tap water to remove cleaning solutions. Inspect handle for material properties (no cracks, chips, or crazing), cleanliness and assure it functions properly.

### 5.2.4 - Drying

Wipe dry with freshly cleaned soft cloth. Prior to reuse, the light handles are required to be steam sterilized.

### 5.2.5 - Disinfection

*Handles could be disinfected with a machine (Getinge washers for example) and rinsed with a temperature maxi. of 93°C.*

Stage	Exp. Temp.	Exp. Time
Prewash	18 - 35° C	60 sec
Wash	46 - 50° C	300 sec
Neutralization	41 - 43°C	30 sec
Wash 2	24 - 28°C	30 sec
Rinsing/Disinfection	92 - 93°C	600 sec
Drying		20 min

### 5.2.6 - Sterilization

Cleaned light handles are to be steam sterilized by using the following instructions and sterilization cycle:

Country	Sterilization cycle	Exp. Temp [°C]	Exp. Time [min]	Dry [min]
USA & Canada	Prevacuum <sup>2</sup>	132 - 135	10	16
France	ATNC (Prion) (Prevacuum)	134	18	
Autres pays	Prevacuum	Follow local / national regulations requirements		

- Re-inspect each handle to assure it is clean and free of defects before proceeding.
- Wrap handle with sterilization wrapper material (double wrapper or equivalent). Alternatively, handles may be placed in paper/plastic sterilization pouches<sup>3</sup> to ease identification after reprocessing and prior to reuse.
- Place on perforated tray with handle opening facing downwards<sup>4</sup>.
- Package with biological and/or chemical indicators for monitoring the sterilization process, per local regulations and facility policies.
- Run sterilization cycle, abiding to sterilizer manufacturer's instructions for load configurations.

**To ensure correct sterilization, do not allow any soiling to penetrate inside the handle. With the sterilization parameters indicated, the sterilizable handles are only guaranteed for no more than 50 utilisation cycles. Disposal of the sterilizable handles is subject to the standard rules governing disposal of hazardous products in a hospital environment.**

<sup>1</sup> A non-enzyme-based detergent is recommended (e.g., Getinge USA NeutraWash). Many enzyme based detergents can adversely affect the plastic material used in this handle and should not be used for extended periods of soaking and will need to be thoroughly rinsed off, if used.  
<sup>2</sup> This light handle is a porous instrument.  
<sup>3</sup> Possible suppliers of sterilization pouches: Medical Action Industries  
SBW Medical  
Baxter International  
<sup>4</sup> For air removal and drying purposes

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## 6 - REPLACING THE BULB

Make sure that the power supply is turned off.

1 - Turn the quarter-turn knob and open the upper cover.

- **CAUTION** : Before starting any maintenance operation, allow the cupola bulb holder to cool for 5 minutes.

2 - Remove the lamp holder taking care not to knock the hot globe of the bulb against hard surface.

- **CAUTION** : Never handle the bulb and the bulb holder with your bare hands. Always use a dry, non-greasy cloth to remove a faulty bulb.

Use a clean, dry cloth to grasp the bulb and pull it out of the lamp holder.

- Take the new bulb.
- Remove the bulb from its cover. Insert the bulb pins onto the bulb holder.
- Check the right position of the bulb.
- Reinstall the bulb holder, checking its position into the cupola.
- Close the cover and lock the assembly by the quarter-turn knob.

The unit must be operated with the cover closed. When performing maintenance, take necessary precautions with respect to hot surfaces which are indicated by appropriate markings.

To prevent any failures during operations, we recommend preventive replacement of the bulbs approximately every 600 hours to 800 hours, depending on the use.

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## 7 - ADJUSTMENTS

*The only adjustments to perform are to be done on the spring arm (A).  
No adjustment is necessary on main arm (B) nor on cupola (C).*

### **Adjustment of weight balance:**

*For the adjustment of the balance weight remove screw (1) at the locking sleeve (2) and move the sleeve back.*

### **Spring arm too strong:**

*Move arm until adjustment nut (3) becomes visible in slot (4).*

*Turn adjustment nut with rod (5) counter clock-wise, until spring arm is balanced and weight is carried properly in every position.*

### **Spring arm too weak:**

*Move arm until adjustment nut (3) becomes visible in slot (4).*

*Turn adjustment nut with rod (5) clock-wise, until spring arm is balanced and weight is carried properly in every position.*

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## 8 - MAINTENANCE

### 8.1 - Preventive maintenance

**Maintenance and inspection (to be performed by a duly trained and certified technician).**

*To preserve your surgical lights' original performance and reliability, annual maintenance and inspections should be performed as follows:*

- during the warranty period, by an MAQUET SAS technician or an MAQUET SAS-certified distributor,*
- outside the warranty period, by an MAQUET SAS technician or MAQUET SAS-certified distributor, or the hospital's engineering department..*

**MAQUET SAS will not accept liability resulting from handling operations which do not comply with the manufacturer's instructions.**

### 8.2 - First level maintenance

#### a) Daily inspection :

- *Check that bulbs operate correctly.*
- *Check that sterilizable handle clicks and locks in correctly.*
- *Check DIMMER ON/OFF and ± operation.*
- *Check arm's position.*

#### b) Monthly inspection :

- *Check Stop's presence.*

#### c) Annual inspection (to be done by an authorized technician):

*Bulbs have a maximum service life of 1 000 hours.  
The bulb holder's average service time is one year.*

## SAFETY POINTS

### Check the following :

- *Tightness of securing-screws on suspension/ceiling tube interface and rubber seal position.*
- *Mounting of Satellite arm(s).*
- *All lighthoods attachment stops.*
- *Disassembly of lighthoods and lubrication of sleeve.*
- *Attachment of all covers and caps.*
- *Replace the spring arm every 6 years.*

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**OTHER CHECKS:**

- *Nominal illumination: 100,000 lux*
- *Continuity of grounds 0.1 Ohm max.*
- *State of lenses and thermal filter.*
- *Condition of lamp holder and bulbs.*
- *Verticality of suspension tube*
- *Balancing adjustment.*
- *Locking mechanism of sterilizable handle.*

d) Cleaning

*Clean unit using damp cloth with water and soap.*  
*Do not use alcohol to clean underside.*

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**10 - STORAGE CONDITIONS**

Room temperature: -25 to 70°C  
 Relative humidity: 10-75%  
 Pressure: 500-1,060 hPa

- Store only in enclosed or covered rooms.
- Do not subject to strong vibrations.

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<b>Table 201 - Guidance and Manufacturer's Declaration - Electromagnetic Emissions</b>		
<i>AXCEL surgical lights are intended for use in the electromagnetic environment specified below. AXCEL customers or users should ensure that they are used in such an environment.</i>		
<b>Immunity test</b>	<b>Compliance</b>	<b>Electromagnetic environment — directives</b>
<i>RF emissions CISPR 11</i>	<i>Group 1</i>	<i>AXCEL lights use RF energy only for their internal functions. Their RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.</i>
<i>RF emissions CISPR 11</i>	<i>Class A</i>	<i>AXCEL lights are suitable for use in all establishments other than domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.</i>
<i>Harmonic emissions IEC 61000-3-2</i>	<i>Not applicable</i>	
<i>Voltage fluctuations and flicker emissions IEC 61000-3-3</i>	<i>Not applicable</i>	

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**Table 202 - Guidance and manufacturer's declaration - electromagnetic immunity**

AXCEL surgical lights are intended for use in the electromagnetic environment specified below. AXCEL customers or users should ensure that they are used in such an environment.

<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment — guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	< 5% $U_T$ (dip > 95% of $U_T$ ) for 0.5 cycles  40% $U_T$ (dip = 60% of $U_T$ ) for 5 cycles  70% $U_T$ (dip = 30% of $U_T$ ) for 25 cycles  < 5% $U_T$ (dip = 95% of $U_T$ ) for 5 cycles	< 5% $U_T$ (dip > 95% of $U_T$ ) for 0.5 cycles  40% $U_T$ (dip = 60% of $U_T$ ) for 5 cycles  70% $U_T$ (dip = 30% of $U_T$ ) for 25 cycles  < 5% $U_T$ (dip = 95% of $U_T$ ) for 5 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the AXCEL requires continuous operation during mains interruptions, it is recommended that it be powered from an uninterruptible power supply.
Power frequency (50/60 Hertz) magnetic fields IEC 61000-4-8	3 A	3 A	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note:  $U_T$  is the voltage of the AC supply network before applying the test level.

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**Table 204 - Guidance and manufacturer's declaration - electromagnetic immunity**

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Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment — guidance
Conducted RF IEC 61000-4-6  Radiated RF IEC 61000-4-3	3 V rms 150 kHz to 80 MHz  3 V/m 80 MHz to 2.5 GHz	3 V rms  3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the AXCEL, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = \left[ \frac{3.5}{V1} \right] \sqrt{P} = 1.17 \sqrt{P}$ $d = \left[ \frac{3.5}{E1} \right] \sqrt{P} \text{ from 80 MHz to 800 MHz} = 1.17 \sqrt{P}$ $d = \left[ \frac{7}{E1} \right] \sqrt{P} \text{ from 800 MHz to 2.5 GHz} = 2.34 \sqrt{P}$ <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, and <math>d</math> is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic survey on site<sup>a</sup>, should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the AXCEL is used exceeds the applicable RF compliance level above, the AXCEL should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the AXCEL.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m..

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**Table 206 - Recommended separation distance between portable and mobile RF communication devices and AXCEL lights.**

AXCEL is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the AXCEL can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the AXCEL as recommended below, according to the maximum output power of the communications equipment.

Maximum rated output power of transmitter <b>W</b>	Separation distance depending on transmitter frequency <b>m</b>		
	150 kHz to 80 MHz $d = [ 1.17 ] \sqrt{P}$	80 MHz to 800 MHz $d = [ 1.17 ] \sqrt{P}$	800 MHz to 2.5 GHz $d = [ 2.34 ] \sqrt{P}$
	0,01	0.12	0.12
0.1	0.37	0.37	0.74
1	1.17	1.17	2.34
10	3.69	3.69	7.38
100	11.70	11.70	23.40

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (*m*) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (*W*) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

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Caution: Disassembly of certain sub-assemblies can affect operation and safety, for example:

- when servicing the electrical power supply,
- when servicing the suspension arms and balancing system,
- when servicing the optical system of the copulas which are equipped with filters intended to eliminate radiation not visible by the patient. The unit should not be operated, under any circumstances, without these filters.

For this type of service, contact the MAQUET SAS-approved product support service.

## TROUBLESHOOTING

ANOMALY	PROBABLE CAUSE	CORRECTIVE ACTION
<input type="checkbox"/> <i>Cupola does not come on.</i>	1 - Bulb failure or missing bulb.	Switch off power supply and replace bulb(s).
	2 - Other cause.	Call MAQUET SAS Service engineering department.
	3 - Mains power failure	Check if another piece of equipment is operating on the same electrical power system.
<input type="checkbox"/> <i>All lighthead not lit</i>	- Each cupola has its own individual control.	Check that all the power supplies are properly switched on.
<input type="checkbox"/> <i>Low bulb service life</i>	- Improper bulb or overvoltage condition.	1) Ensure you are only using MAQUET SAS bulbs specified. Do not use bulbs used with surgical lights of previous generation (ANGENIEUX AX4-AX14) which are not compatible in power and size.  2) Verify voltage at bulb terminals.  $23,5 V_{EFF} \left( V_{AC+DC}^{RMS} \right)$
<input type="checkbox"/> <i>Sterilizable handle does not engage correctly on its support.</i>	1 - Average service life of sterilizable handle: 50 sterilization cycles at 134° ± 4°C for 20 minutes.	Check operation of handle lock mechanism (click) and overall handle.
	2 - Maximum service life exceeded / Handle is deformed.	Replace handle
<input type="checkbox"/> <i>Cupola drift</i>	- Suspension tube not levelled. - Unstable ceiling structure	- Check level and ceiling structure. - Call MAQUET SAS Service engineering department.
<input type="checkbox"/> <i>Cupola too loose or too difficult to handle.</i>	- Brake misadjustment. - Lack of lubrication.	- Proceed with Balance Adjustment (See Chap. 3.3). - Call MAQUET SAS Service engineering department.

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## INSPECTION SHEET

(copy to be sent to your MAQUET SAS Agency or Distributor)

MAQUET SAS project No. _____	Serial No. _____
Date of delivery or of dispatch FOB _____	Product Ref. _____
Date of installation _____	Operating theatre ref. _____
Date of commissioning _____	Operating room ref. _____

After Sales Department

Stamp

Customer identification

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

**CUSTOMER'S ACCEPTANCE :**

Customer (print name) : \_\_\_\_\_

Office : \_\_\_\_\_

Easy handiness :            yes    no

Correct stability :          yes    no

Installer (print name)	Signature	Date
Customer (print name)	Signature	Date

**Inspections to be carried out before placing in service**

**Inspection carried out**

- |   |  |
|---|--|
| 1 Ground continuity: from bulb cover hinge to cupola(s) to suspension tube ground screw         | <input type="checkbox"/>                             |
| 2 Operation of light bulbs  | <input type="checkbox"/>                             |
| 3 Maximum voltage at bulb-holder $23,5 V_{EFF} (V_{AC+DC}^{RMS})$                               | <input type="checkbox"/>                             |
| 4 Verticality of suspension tube  | <input type="checkbox"/>                             |
| 5 Balancing of intermediate arm (main cupola)<br>Balancing of suspension arm (satellite cupola) | <input type="checkbox"/><br><input type="checkbox"/> |
| 6 Click-in locking of sterilizable handles  | <input type="checkbox"/>                             |
| 7 General aspect and cleaning   | <input type="checkbox"/>                             |

(copy to be sent to MAQUET SAS Agency or Distributor)

## INSPECTION SHEET

(copy to be kept by the customer)

MAQUET SAS project No. _____	Serial No. _____
Date of delivery or of dispatch FOB _____	Product Ref. _____
Date of installation _____	Operating theatre ref. _____
Date of commissioning _____	Operating room ref. _____

After Sales Department

Stamp

Customer identification

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

**CUSTOMER'S ACCEPTANCE :**

Customer (print name) : \_\_\_\_\_

Office : \_\_\_\_\_

Easy handiness :            yes    no

Correct stability :            yes    no

Installer (print name)	Signature	Date
------------------------	-----------	------

Customer (print name)	Signature	Date
-----------------------	-----------	------

**Inspections to be carried out before placing in service**

**Inspection carried out**

- |   |                          |
|---|--------------------------|
| 1 Ground continuity: from bulb cover hinge to cupola(s) to suspension tube ground screw | <input type="checkbox"/> |
| 2 Operation of light bulbs  | <input type="checkbox"/> |
| 3 Maximum voltage at bulb-holder $23,5 V_{EFF} (V_{AC+DC}^{RMS})$                       | <input type="checkbox"/> |
| 4 Verticality of suspension tube  | <input type="checkbox"/> |
| 5 Balancing of intermediate arm (main cupola)   | <input type="checkbox"/> |
| Balancing of suspension arm (satellite cupola)  | <input type="checkbox"/> |
| 6 Click-in locking of sterilizable handles  | <input type="checkbox"/> |
| 7 General aspect and cleaning   | <input type="checkbox"/> |

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