

# Cross-Checks L

## CI 111

The Cross-Checks L Multicritical Process Variable indicator strips are designed for use in steam sterilizers operating at 121°C, 132°C and 134°C. When used as directed, the Cross-Checks L indicator strips give visible indication that sterilizing conditions were met.



### Key Features

- 121°C, 132°C and 134°C
- Requires all parameters for steam sterilization to be met before color change can occur

### Performance Characteristics

When all steam sterilization conditions are met, the check marks change color from white to brown/black

### Physical Characteristics

Dimensions: L 89mm x W 19mm

### Color Change

White to brown/black

### Stated Values

(as determined in a steam resistometer)

- 9.5 minutes at 121°C
- 3.5 minutes at 132°C
- 3 minutes at 134°C

### Storage

The storage conditions are printed on the product packaging. Store in resealable bag provided

### Shelf Life

36 months from date of manufacture. The expiry date is printed on the product packaging

### Standards Classification

ISO 11140-1, Type 4

### SteriTec Part No.

CI 111, Cross-Checks L

### Packaging

24 bags of 250