

Dry Heat Indicator Label

CI 113

Dry Heat Indicator Labels are designed for use in a dry heat sterilizer operating at 160°C. The label provides a clear identification of processed versus unprocessed items.

When used as directed, the Dry Heat Label gives visible confirmation that the item was exposed to a dry heat process.



Key Features

- For use in dry heat sterilizers operating at 160°C
- Color change in 10 minutes + / - 2 minutes

Performance Characteristics

When exposed to dry heat conditions, the indicator bar changes color.

Physical Characteristics

Dimensions: L 29mm x W 19mm

Color Change

Green to dark brown/black

Complete Color Change Time

10 minutes +/- 2 minutes at 160°C

Storage

The storage conditions are printed on the product packaging

Shelf Life

36 months from date of manufacture. The expiry date is printed on the product packaging

FDA 510(k)

K953122

Standards Classification

- FDA Process Indicator
- ISO 11140-1, Type 4

SteriTec Part No.

CI 113, Dry Heat Indicator Label

Packaging

10 boxes of 1,000

FTP-009 Rev. G