

Steri-Vapeur 18 Emulator

CI 117

The Steri-Vapeur 18 emulators are designed for use in steam sterilizers operating at 134°C for 4 minutes and 18 minutes. When used as directed, the Steri-Vapeur 18 gives visual confirmation that sterilization conditions were met at the location in the sterilizer chamber.



Performance Characteristics

The indicator has two different timing indicator inks printed on each indicator. The user can determine:

1. If the 4 MIN spot turns green, standard sterilization conditions have been met.
2. If both the 4 MIN and 18 MIN spots turn green, this indicates that an extended sterilization cycle has been performed.

Physical Characteristics

Dimensions: L 87mm x W 25mm

Color Change

Purple to green

Storage

The storage conditions are printed on the product packaging. Store in the resealable bag provided

Shelf Life

24 months from date of manufacture. The expiry date is printed on the packaging

Standards Classification

- 4 MIN - ISO 11140-1, Type 4
- 18 MIN - ISO 11140-1, Type 6

SteriTec Part No.

CI 117, Steri-Vapeur 18

Packaging

16 bags of 200