

Dry-Checks Dry Heat Indicator Strip

CI 136

The Dry-Checks Dry Heat Indicator Strips are designed for use in dry heat sterilizers at 160°C. The indicator provides clear identification of processed and unprocessed items. When used as directed, Dry-Checks give visible confirmation that the sterilized item was exposed to dry heat.



Key Features

- For use in dry heat sterilizers operating at 160°C
- Color change in 7 minutes

Performance Characteristics

When exposed to dry heat conditions, the indicator check mark changes color.

Physical Characteristics

Dimensions: L 97mm x W 19mm

Color Change

Green to black

Complete Color Change Time

7 minutes at 160°C

Storage

The storage conditions are printed on the product packaging. Store in the resealable bag provided

Shelf Life

36 months from date of manufacture. The expiry date is printed on the product packaging

Standards Classification

ISO 11140-1, Type 4

SteriTec Part No.

CI 136, Dry Checks
Dry Heat Indicators

Packaging

20 boxes of 250