

# Cross-Checks LS

## CI 139

The Cross-Checks LS Multicritical Process Variable Indicator strips are designed for use in pre-vacuum steam sterilizers operating at 121°C, 132°C and 134°C. When used as directed, the Cross-Checks LS indicator strips give visible indication that sterilizing conditions were met.



### Key Features

- 121°C, 132°C and 134°C
- Sterilant specific - must have steam to change
- Lead and toxic heavy metal free

### Performance Characteristics

When exposed to a steam sterilization process, the indicator bar changes color.

### Physical Characteristics

Dimensions: L 97mm x W 14mm

### Color Change

Purple to green

### Stated Values *(as determined in a steam resistometer)*

- 5 minutes at 134°C
- 5.5 minutes at 132°C
- 15 minutes at 121°C

### Storage

The storage conditions are printed on the product packaging. Store in the sealed bag provided

### Shelf Life

24 months from date of manufacture. The expiry date is printed on the product packaging

### Standards Classification

ISO 11140-1, Type 4

### SteriTec Part No.

CI 139, Cross-Checks LS

### Packaging

20 bags of 250