

Getinge Clean Heavy Soil Detergent

Strong, alkaline based cleaner

Getinge Clean Heavy Soil Detergent is for heavily soiled dried items such as gynecological and orthopedic instruments

Getinge Clean Heavy Soil Detergent is a strong high alkaline-based cleaner, suitable for alkaline-resistant materials including stainless steel instruments and glass. This product is not suitable for aluminum and anodized aluminum.

Getinge Clean Heavy Soil Detergent is suitable for all types of washer-disinfectors. It is effective for removing blood, fat, grease and oil.



Available in 5L, 10L, 25L and 205L containers

This document is intended to provide information to an international audience outside of the US.

For more detailed information, please contact your local Getinge representative at getinge.com/contact

www.getinge.com

GETINGE 