

# Getinge Clean Neutralizer Plus

Multi purpose acid neutralizer and cleaner

Getinge Clean Neutralizer Plus formulation includes a blend of orthophosphoric and nitric acid

Getinge Clean Neutralizer Plus removes rust, stains and hard water scale from stainless steel instruments, equipment and washer-disinfector chambers. It also restores stainless steel surface finish.

Neutralizer Plus can also be used as a process additive for pH neutralization of alkaline residues.

This formulation is suitable for use on stainless steel and glass.



Available in 5L containers

This document is intended to provide information to an international audience outside of the US.

For more detailed information, please contact your local Getinge representative at [getinge.com/contact](http://getinge.com/contact)