

# Getinge Clean Neutralizer

Mildly acidic, based on a citric acid

Getinge Clean Neutralizer is used to neutralize alkaline residues from the wash phase

A mild organic acid for use in washer-disinfectors for neutralization of alkaline detergents. This product removes hard water scale and other acid soluble minerals on medical devices.

Getinge Clean Neutralizer effectively reduces the pH of rinse water after alkaline processes. It is suitable for all acid resistant surfaces including stainless steel, ceramics, glass and plastics.



Available in 5L, 10L, 25L and 205L containers

This document is intended to provide information to an international audience outside of the US.

For more detailed information, please contact your local Getinge representative at [getinge.com/contact](http://getinge.com/contact)